

## Quad Canes, Large Base

HCPCS: E0105



FIGURE A

- Four point base, combined with offset handle provides additional stability and support
- Easy-to-use, one-button height adjustment with locking bolt, allows quick, safe fitting
- Available in Black (main image), Chrome (Figure A), Bronze (Figure B)
- No assembly required
- Handle height adjusts from 30" to 39"
- Handle available in foam rubber (main image) or vinyl (Figure C)
- Weight Capacity 300 lbs

### PRODUCTS

### Description

10300-4	Quad Canes, Large Base, CS
---------	----------------------------

### SPECIFICATIONS

### Description

Handle Height	29" - 37.5"
Warranty	Limited Lifetime
Weight Capacity	300 lbs