

V2™ Full-Face CPAP Mask

Redesigned for Improved Comfort and Seal



V2 Full Face
CPAP Mask

The V2 full face CPAP mask has been **redesigned to be ultimately more comfortable and to provide a stable seal** for a wider array of facial profiles. These improvements include incorporating a membrane inside the cushion for an **improved seal**. In addition, the cushion surface features a series of flanges around the base and is **constructed of a thinner, softer and more comfortable silicone material**. Two structural braces have been added to aid in cushion stability. Finally, **the CO₂ vents were moved** from the bridge of the nose to just above the anti-asphyxia valve, **pointing the airflow perpendicular to the body**.

Other improvements include **no forehead support, single-molded mask frame and cushion that are made of easy contouring silicone rubber**. The one piece design and soft material allows for greater flexibility through a broad range of 5 different mask sizes.

Performance Features

- Headgear uses **new quick-release strap clips** and a foam material for improved comfort
- **360-degree mask swivel** for expanded patient freedom
- Detachable swivel allows user to **disconnect from equipment without removing mask**
- **5 sizes of mask** are available (mask sizes tend to run large)
- **3 sizes of headgear** are available (small, medium, large)

Design Features

- Lighter weight than original design
- Single-molded mask frame and cushion with **integrated chin support**
- Anatomically contoured design with ribbed support to improve seal and comfort
- Integrated facial-sealing flange for improved fit
- **Anti-Asphyxia Valve** that opens at mask pressure of less than 3 cm H₂O
- Swivel port micro-vent holes that dissipate CO₂
- **Detachable swivel port** that allows the hose to be easily removed
- **Easily adjustable 4-strap headgear** made of durable and comfortable polyurethane foam
- Reusable mask is steam autoclaveable and dishwasher safe. Can also be cleaned via cold chemical and pasteurization processes
- Latex-free
- Multi-patient, multi-use products

V2 Ordering Information

(includes headgear)

PRODUCT DESCRIPTION	PART NUMBER
Large	113485
Medium	113486
Small	113487
Extra Small	113488
Petite	113489

Reimbursement Table

PRODUCT DESCRIPTION	CODE DESCRIPTION	HCPCS CODING	FREQUENCY	ALLOWABLE RANGE
Face Piece and Swivel Port	Full face mask used with PAP device, each	A7030	1 per 3 months	\$148.44 – \$174.64
Headgear	Headgear used with PAP device, each	A7035	1 per 6 months	\$31.27 – \$36.79
Face Piece	Face mask interface, replacement for full face mask, each	A7031	1 per 1 month	\$54.90 – \$64.59

The noted allowable rates are based off of the January 2012 HCPCS Fee Schedule. All information presented is for reference purposes only.

Accessories and Parts

PRODUCT DESCRIPTION	PART NUMBER	PRODUCT DESCRIPTION	PART NUMBER
Face Pieces		Face Piece with Swivel Port	
Large	669180	Large	113570
Medium	669181	Medium	113571
Small	669182	Small	113572
Extra Small	669183	Extra Small	113573
Petite	669184	Petite	113574
Headgear		Mask Sizing Gauge	
Large	201494		691143
Medium	201493		
Small	201492		



DeVilbiss Healthcare LLC | 100 DeVilbiss Drive | Somerset, PA 15501 USA
800.338.1988 | 814.443.4881 | www.DeVilbissHealthcare.com