

Now Available in
Medium Tall sizes!

JOBST® for Men



Effective leg
therapy in a
high-quality
dress sock

Medical Compression Stockings
15-20 mmHg*, 20-30 mmHg*, 30-40 mmHg*

#1 physician recommended

JOBST®

Comfort, Health and Style!

forMen features and benefits:

Stylish

- Fashionable ribbed design delivers gradient compression therapy, yet looks like a fine men's dress sock

Durable

- Reinforced heel resists shoe abrasion and strong yarns resist wear and tear

Comfortable

- All-day comfort band keeps socks up without binding or pinching
- Smooth toe seam and bubble shape add comfort

Cool and with odor control formula




- Microfibers wick away moisture six times faster than regular nylon for a dryer, cooler feeling
- Special formula fights odor by inhibiting growth of fungi and bacteria on the socks

Available colors:



Actual garment color may vary from printed swatch sample

Sizing chart:

Size	Ankle 	Calf 	Thigh 
S	7"-8 ¹ / ₄ " (18-21 cm)	11"-15" (28-38 cm)	15 ³ / ₄ "-24 ³ / ₈ " (40-62 cm)
M	8 ³ / ₈ "-9 ⁷ / ₈ " (21-25 cm)	11 ⁷ / ₈ "-16 ¹ / ₂ " (30-42 cm)	18 ¹ / ₈ "-27 ¹ / ₂ " (46-70 cm)
L	10"-11 ³ / ₈ " (25-29 cm)	12 ¹ / ₂ "-18 ¹ / ₈ " (32-46 cm)	21 ¹ / ₄ "-30 ³ / ₄ " (54-78 cm)
XL	11 ¹ / ₂ "-13" (29-33 cm)	13 ³ / ₈ "-19 ⁵ / ₈ " (34-50 cm)	23 ⁵ / ₈ "-32" (60-81 cm)
MT	8 ³ / ₈ "-9 ⁷ / ₈ " (21-25 cm)	11 ⁷ / ₈ "-16 ¹ / ₂ " (30-42 cm)	_____
LT	10"-11 ³ / ₈ " (25-29 cm)	12 ¹ / ₂ "-18 ¹ / ₈ " (32-46 cm)	_____
LFC	10"-11 ⁷ / ₈ " (25-30 cm)	18"-24" (46-61 cm)	_____
XLFC	12"-14" (30-36 cm)	18"-24" (46-61 cm)	_____

Care instructions:

Machine washable in lukewarm water. Do not bleach.
 Machine dry - low heat. Do not dry clean. Do not iron.
 Use a delicate fabric detergent. Rinse well. Do not wring.



Comfort, Health & Style!

Knee High ordering information:

Knee High (Closed Toe)



Size	Catalog Number/Color					Compression mmHg*
	Black	Khaki	White	Brown	Navy	
S	115000	115012	115008	115084	115004	15-20 mmHg
M	115001	115013	115009	115085	115005	15-20 mmHg
L	115002	115014	115010	115086	115006	15-20 mmHg
XL	115003	115015	115011	115087	115007	15-20 mmHg
LT	115068					15-20 mmHg
S	115088	115100		115096	115092	20-30 mmHg
M	115089	115101		115097	115093	20-30 mmHg
L	115090	115102		115098	115094	20-30 mmHg
XL	115091	115103		115099	115095	20-30 mmHg
MT	115255				115256	20-30 mmHg
LT	115260				115262	20-30 mmHg
LFC	115294					20-30 mmHg
XLFC	115295					20-30 mmHg
S	115108	115120		115116	115112	30-40 mmHg
M	115109	115121		115117	115113	30-40 mmHg
L	115110	115122		115118	115114	30-40 mmHg
XL	115111	115123		115119	115115	30-40 mmHg
MT	115257				115258	30-40 mmHg
LT	115265				115267	30-40 mmHg
LFC	115296					30-40 mmHg
XLFC	115297					30-40 mmHg

Knee High (Open Toe)



Size	Catalog Number/Color					Compression mmHg*
	Black	Khaki	White	Brown	Navy	
S	115432					20-30 mmHg
M	115433					20-30 mmHg
L	115434					20-30 mmHg
XL	115435					20-30 mmHg
LFC	115372					20-30 mmHg
XLFC	115373					20-30 mmHg
S	115452					30-40 mmHg
M	115453					30-40 mmHg
L	115454					30-40 mmHg
XL	115455					30-40 mmHg
LFC	115374					30-40 mmHg
XLFC	115375					30-40 mmHg

Fiber contents:

15-20 mmHg*: Closed toe: 76% Nylon and 24% Spandex/Elastan

20-30 mmHg*: Closed toe: 73% Nylon and 27% Spandex/Elastan; Open Toe: 71% Nylon and 29% Spandex/Elastan

30-40 mmHg*: Closed toe: 73% Nylon and 27% Spandex/Elastan; Open Toe: 71% Nylon and 29% Spandex/Elastan



Thigh High ordering information:

Thigh High (Closed Toe)



Size	Catalog Number/Color		Compression mmHg*
	Black	Khaki	
S	115516	115512	15-20 mmHg
M	115517	115513	15-20 mmHg
L	115518	115514	15-20 mmHg
XL	115519	115515	15-20 mmHg
S	115408	115400	20-30 mmHg
M	115409	115401	20-30 mmHg
L	115410	115402	20-30 mmHg
XL	115411	115403	20-30 mmHg
S	115412	115404	30-40 mmHg
M	115413	115405	30-40 mmHg
L	115414	115406	30-40 mmHg
XL	115415	115407	30-40 mmHg

Fiber contents:

15-20 mmHg*: 73% Nylon, 20% Spandex/Elastan, 7% Silicone

20-30 mmHg*: 71% Nylon, 25% Spandex/Elastan, 4% Silicone

30-40 mmHg*: 71% Nylon, 25% Spandex/Elastan, 4% Silicone

Special Care Instructions for the JOBST Thigh High Style:

The use of skin care products (such as skin creams, body lotion, or powder) as well as various skin types (such as dry or peeling skin) and hair growth may prevent the silicone band from keeping the thigh high in place.

To improve the performance of the silicone band under the above conditions, we recommend the following steps:

1. Turn the JOBST Thigh Highs inside out before washing.
2. Before the main wash, thoroughly hand wash the silicone band using a delicate fabric detergent. Concentrate on removing any lotion, powder, hair or other foreign material from the silicone.

Please refer to the product label and/or package insert for full instructions on the safe use of these products.

*The mean compression for the average ankle size.