


Med-Aire Plus 8" Alternating Pressure and Low Air Loss Mattress System



INSET A

- Active combination therapy system provides both alternating pressure and low air loss to optimize pressure redistribution, shear/friction reduction, and microclimate control
- Designed to prevent, treat and heal pressure ulcers in the home or long term care setting
- Digital control unit offers 4 cycle times and adjustable comfort weight settings for optimal immersion and patient satisfaction
- 3 diagnostic audio and visual alarms alert for issues with air pressure, the system, or power
- System includes 20 individual air cells to maximize immersion and envelopment with pillow feature for head stability & comfort
- Advanced Clinical and Care Giver Features: Auto-Firm, Seat Inflate, microprocessing feedback, lock-out and more!
- Cell-on-cell designed air cells prevent “bottoming out” and provide up to 24 hours of power outage protection
- 8" deep air cells are easily removed with quick disconnects for convenient servicing
- Removable 4 way stretch PU cover is fluid resistant, low shear and vapor permeable helping protect the skin from friction and moisture
- CPR valve allows for rapid deflation
- Tubing quickly and easily disconnects with EZ Lock/Quick Release design

PRODUCTS

PRODUCTS	Description
14029	Alternating Pressure/Low Air Loss Mattress System, E0277, 1/ea
14029-84	Alternating Pressure/Low Air Loss Mattress System, E0277, 1/ea

SPECIFICATIONS

SPECIFICATIONS	Description
Control Unit Airflow	8 LPM
Control Unit Cycle Time	10, 15, 20, 25 Minutes
Control Unit Type	Compressor
Control Unit Alarms	Visual / Audible
Mattress Cover Material	Fluid-Resistant Stretch
Weight Capacity	450 lbs.
Warranty	18 Month Limited

ACCESSORIES

ACCESSORIES	Description
14029MN	Mattress Only, 14029dp, 3 Tube System, 1/ea
14029XP	Pump Only for 14029/14029DP, 1/ea
14029C	Cover, 1/ea
14029-CELL	Standard Air Cell, 1/ea
14029-CELL 1	Ventilated Air Cell, 1/ea
14427	Electric Pump (Sold Separately) (Replaces Mason # 8151-EP), 1/ea
14029B	Mattress Carry Bag, 1/ea

PARTS

PARTS	Description
14029MN	Mattress Only, 14029dp, 3 Tube System, 1/ea

PARTS	Description
14029C	Cover, 1/ea
14029-base	Base, 1/ea
14029-CELL	Standard Air Cell, 1/ea
14029-CELL-1	Bladder with Laser Holes, 1/ea
14001-l	"L" Connector, 1/ea
14001-t	"T" Connector, 1/ea
14029-fftc	Quick release "T" connector, 1/ea
14029-fmtc	Bladder connector, 1/ea
14028v	One Way Valve For Low Air Loss Matt, 1/ea
14013	CPR plug, 1/ea
14028t	Tubing, 1/ea
14030xp-s	Pump socket, 1/ea
14007-3port	EZ lock Q/R connector, 1/ea
14007or-3	"O" Ring, 1/ea
14030hk	Hanging hook, 1/ea
14029fn	Intake filter, 1/ea
14029xp-fc	Intake filter cover, 1/ea
14530f	Fuse (1 amp 250V), 1/ea
15005PC	Power cord, 1/ea
14029B	Mattress Carry Bag, 1/ea
DS16C01207	Mattress Hose 3 port
14029-FFTCPLUG	Plug For Connector