

Med-Aire Plus 10" Bariatric Alternating Pressure and Low Air Loss Mattress Replacement System



INSET A

- Active combination bariatric therapy system provides both alternating pressure and low air loss to optimize pressure redistribution, shear/friction reduction, and microclimate control
- Designed to prevent, treat and heal all stages of pressure ulcers in the home or long term care setting
- Digital control unit offers 4 cycle times and adjustable comfort weight settings for optimal immersion and patient satisfaction
- 3 diagnostic audio and visual alarms alert for issues with air pressure, the system, or power
- System includes 20 individual air cells to maximize immersion and envelopment with pillow feature for head stability & comfort
- Advanced Clinical and Care Giver Features: Auto-Firm, Seat Inflate, auto-recovery, microprocessing feedback, lock-out and more!
- Cell-on-cell designed air cells prevent “bottoming out” and provide up to 24 hours of power outage protection
- 10" deep air cells are easily removed with quick disconnects for convenient servicing
- Removable quilted Nylon cover is fluid resistant, low shear and vapor permeable helping protect the skin from friction and moisture
- CPR valve allows for rapid deflation
- Tubing quickly and easily disconnects with EZ Lock/Quick Release design

PRODUCTS

| PRODUCTS | Description |
|----------|---|
| 14030 | Bariatric Alternating Pressure Mattress System, E0277, 1/ea |
| 14048 | Bariatric Alternating Pressure Mattress System, E0277, 1/ea |
| 14054 | Bariatric Alternating Pressure Mattress System (Replaces Mason # LS9000B-54), E0277, 1/ea |
| 14060 | Bariatric Alternating Pressure Mattress System (Replaces Mason # LS9000B-60), E0277, 1/ea |

SPECIFICATIONS

| SPECIFICATIONS | Description |
|-------------------------|--|
| Control Unit Airflow | 12 LPM |
| Control Unit Cycle Time | 10, 15, 20, 25 Minutes |
| Control Unit Type | Compressor |
| Control Unit Alarms | Visual / Audible |
| Mattress Dimensions | 42", 48", 54", 60" (W) X 80" (L) X 10" (H) |
| Mattress Cover Material | Fluid-Resistant, 4-way Stretch |
| Weight Capacity | 600 lbs. - 1000 lbs. |
| Warranty | 18 Month Limited |

ACCESSORIES

| ACCESSORIES | Description |
|-------------|--|
| 14427 | Electric Pump (Sold Separately) (Replaces Mason # 8151-EP), 1/ea |
| 14030XP | Replacement pump for 14030, 14048, 14054 & 14060, 1/ea |
| 14030MN | Replacement Mattress for 14030, 1/ea |
| 14048MN | Replacement Mattress for 14048, 1/ea |
| 14054MN | Replacement Mattress for 14054, 1/ea |

| ACCESSORIES | Description |
|-----------------|--------------------------------------|
| 14060MN | Replacement Mattress for 14060, 1/ea |
| 14030C | Replacement Cover, 1/ea |
| 14048C | Replacement Cover, 1/ea |
| 14054C | Replacement Cover, 1/ea |
| 14060C | Replacement Cover, 1/ea |
| 14048-CELL | Replacement Bladder, standard, 1/ea |
| 14030-CELL | Replacement Bladder, standard, 1/ea |
| 14051-CELL | Standard Cell, 1/ea |
| 14054-CELL | Replacement Bladder, standard, 1/ea |
| 14060-CELL | Replacement Bladder, standard, 1/ea |
| 14030-LASERCELL | Replacement Ventilated Bladder, 1/ea |
| 14048-LASERCELL | Replacement Ventilated Bladder, 1/ea |
| 14054-LASERCELL | Replacement Ventilated Bladder, 1/ea |
| 14060-LASERCELL | Replacement Ventilated Bladder, 1/ea |