

## Adjustable Tension General Use Wheelchair Back Cushion



INSET A



14301

- Perfect solution for patients with mild positioning needs
- Versatile back support system that improves patient comfort and posture
- 1" thick high density foam pad provides optimal pressure redistribution
- Secures directly over wheelchair canes for added strength, stability and safety
- 5 Adjustment straps featured for multi-purpose back support and customization to individual seating needs
- Handy back pocket is convenient for carrying keys, water bottle, medications, tissues and more
- Breathable cover is fluid resistant and easy to clean and maintain
- Cushion folds with the wheelchair for easy transport and storage
- Meets flammability standard TB117

### PRODUCTS

### Description

14300	Tension Adjustable Back Cushion, E2611, 1/ea
14301	Tension Adjustable Back Cushion, E2612, 1/ea

### SPECIFICATIONS

### Description

Cover	Fluid resistant, easy to clean, water proof Nylon
Weight Capacity	350 lbs
Cover material	Fluid resistant PU blend
Warranty	18 months