


Gel "E" Skin Protection 3" Gel/Foam Wheelchair Cushion



14891



14891

- Pressure redistributing cushion designed to assist in the prevention, treatment and management of pressure injuries
- Viscous gel bladder is encased in polyurethane foam shell and provides optimal pressure redistribution, support and comfort
- Cover top is a urethane-coated nylon that is low-shear, vapor permeable and water-resistant
- Cover base is waterproof vinyl for durability
- Carry handle included for convenient transport
- Meets flammability standard TB117

PRODUCTS

Description

PRODUCTS	Description
14891	3" Gel/Foam Cushion, E2603, 1/ea
14886	3" Gel/Foam Cushion, E2603, 1/ea
14878	3" Gel/Foam Cushion, E2603, 1/ea
14894	3" Gel/Foam Cushion, E2603, 1/ea
14902	3" Gel/Foam Cushion, E2604, 1/ea

SPECIFICATIONS

Description

SPECIFICATIONS	Description
Cover Material	Fluid Resistant Nylon Top Waterproof Vinyl Bottom
Weight Capacity	Up to 20" 275 lbs. 22" and up 500 lbs.
Warranty	18 Months Limited