

Premium Guard Gel Overlay



14893

- Best in class gel overlay optimizes pressure redistribution and is now available in up to 500 lbs weight capacity
- 15 gel bladders cover entire patient contact area for the prevention and treatment of pressure ulcers
- Full length “head-to-toe” gel bladders minimize “bottoming out” and “pitting”
- New convoluted foam topper vastly improves patient comfort and air circulation
- 4 elastic corner straps secure overlay to mattress
- Urethane coated nylon cover is water resistant and vapor permeable
- Vinyl base is stain resistant and waterproof
- Bladders contain an aqueous non-toxic gel and polystyrene beads to minimize gel migration

PRODUCTS

Description

14893	Gel Overlay, E0185, 1/ea
14901	Bariatric Gel Overlay, E0185, 1/ea

SPECIFICATIONS

Description

Cover Material	Fluid-Resistant Nylon Top Waterproof Vinyl Bottom
Warranty	18 Month Limited
Weight Capacity	34" W- 300 lbs 42" W- 500 lbs

ACCESSORIES

Description

14893C	34" (W) x 76" (L) x 3.5" (H), 1/ea
14901C	replacement cover, 1/ea