

Model 151 Atomizer

HCPCS:



- This 1.4 ounce capacity atomizer produces a fine mist that is ideal for nose and throat treatments, decubitus ulcer care, and wound cleansing
- A vacuum or venturi-type atomizer that can be operated by one of the following: compressed air, a No 1 hand bulb or the No 630 cut-off assembly
- Designed for use with aqueous solutions
- Precise control and application of topical anesthetics
- All types of oil or water solutions can be used that are compatible with rhodium and nickel plating
- Produces fine mist without force
- Small tube may be inserted well into nasal passages
- Removable spray tube for ease in sterilization
- Top consists of all metal construction
- Used in decubitus ulcer treatments, debridement of wounds, and wound cleansing
- Use by E.N.T doctors, research laboratories, veterinarians, burn centers, wound care facilities, emergency rooms

PRODUCTS

Description

151	Model 151 Atomizer, EA
-----	------------------------

ACCESSORIES

Description

150-OT-605	OUTLET TUBE ASSY 5 PK SEPE SALE PKG, PKG
------------	--

150-QB	1.4 OZ CLEAR GLASS, EA
--------	------------------------

251-S-605	WSHER 10-481 STK EQ 5 PK SEP SALE PKG, PKG
-----------	--

630-605	CUTOFF WITH TUBING 5-PK, PKG
---------	------------------------------