

## Chrome Oxygen Cylinder Cart

HCPCS:



- Dual lifting handles
- 5" Swivel casters, two with brakes
- Chrome plated finish, alloy carbon steel
- Module design saves on shipping

### PRODUCTS

### Description

18143	Chrome Oxygen Cylinder Cart, CS
18144	Chrome Oxygen Cylinder Cart, CS

### SPECIFICATIONS

### Description

Height	39.8"   39.8"
Length	22.6"   27.6"
Warranty	Limited Lifetime
Weight	23.8 lbs.   28.7 lbs.
Width	15.4"   20.5"