

Model 286-RD Atomizer

HCPCS:



- The Model 286-RD Atomizer is designed especially for applying a fine mist to nasal and oral cavities and for administering an oral vaccine
- Ultra-fine particle atomization of most medical solutions
- Straight flow tip is easily disassembled for cleaning
- Suitable for applying fine, medicated spray to nasal and oral cavities, or oral vaccines
- May be used with compressed air, the No 1 bulb or Model 630 Bleeder-Type cut-off
- Constructed of polypropylene material for spraying solutions incompatible with metal
- Popular with treatment of animals and veterinary clinics and is used for testing in laboratory research
- Used by E.N.T doctors, research laboratories, veterinarians, burn centers, wound care facilities, emergency rooms

PRODUCTS

Description

286-RD	Model 286-RD Atomizer, EA
--------	---------------------------

ACCESSORIES

Description

630-605	CUTOFF WITH TUBING 5-PK, PKG
1-RD	PLASTISOL BULB (5/32" OPENING), EA
251-C-605	SPRAY TIP 5 PK SEPARATE SALE PACKAGE, PKG
251-S-605	WSHER 10-481 STK EQ 5 PK SEP SALE PKG, PKG
286-F-605	OUTLET TUBE 5 PK SEPARATE SALE PACKAGE, PKG
286-RD-602	286 2 CAP AND 1 BULB PKG, PKG
286-QB	1.4 OZ AMBER GLASS, EA
1-RD-605	PLASTISOL BULB (5/32" OPENING) 5 PK SS, PKG