

Model 45 Nebulizer

HCPCS:



- Model 45 hand-held spray nebulizer delivers aerosol particles of uniformly consistent density deep into the lungs The break-resistant nebulizer easily disassembles for sterilization
- Used for treatment of asthmatic patients
- Designed for comfortable administration of fine mists of uniform density to reach deep into lung cavity
- Converts liquid (i e antibiotic, oil extracts, etc) into aerosol
- Used in respirator and mask test fit kits, and bulk shipped as item 4510
- Break-resistant nebulizer that easily disassembles for sterilization by using one of the following methods: autoclaving, boiling or cold sterilization
- Polycarbonate plastic nebulizer may be used with or without mouthpiece (No 25-31-605)
- Complete with mouthpiece, extra jet, and bulb
- Used by research laboratories

PRODUCTS

Description

PRODUCTS	Description
45	Model 45 Nebulizer, EA

SPECIFICATIONS

Description

SPECIFICATIONS	Description
Capacity (CTPT-3LPZ)	5cc
MMAD	5 Microns