

DV8 Steerable Aluminum Knee Walker

HCPCS: E0118



FIGURE A



FIGURE B




FIGURE C

- Knee walker can be steered for increased maneuverability
- Ideal for individuals recovering from foot surgery, breaks, sprains, amputation and ulcers of the foot
- Deluxe dual-braking system
- Sturdy dual frame adds stability and durability
- Leg pad with "channel" provides added stability and comfort
- Tiller can be easily folded for storage and transportation
- New foldable steering column allows for easy storage
- Tool-free height adjustment of seat and handle (Figure A - B)
- 8" casters are ideal for indoor/outdoor use
- Silver vein finish (Figure C) is easy-to-maintain

PRODUCTS

Description

790	DV8 Steerable Aluminum Knee Walker, EA
-----	--

SPECIFICATIONS

Description

Carton Shipping Weight	26 lbs.
Handle Height	33" - 37.5"
Knee Pad Height	17.5" - 21.5"
Unit Dimensions	17" (W) x 32" (D)
Warranty	Limited Lifetime, Lifetime Warranty on Brake Cable
Weight	19 lbs.
Weight Capacity	300 lbs

PARTS

Description

9503H79015-N	KneePadHeightKnob790New,1 ea, EA
9503H79016	Left Wheel Lock,790,1 ea, EA
9503H79017	Right Wheel Lock,790,1 ea, EA
9503H79018	Caster Complete Assy,790,1 ea, EA
9503H79019	Caster Axle,790,1 ea, EA
9503H79020	Axle Set Screws,790,4/pkg, PKG
9503H79021	Wheel Bearings,790,PhoenixHD3,1/pr, PR
9503H79001	Handle Bar,790,1 ea, EA
9503H79002	Hand Grips,790,1 Pair, PR
9503H79003	Left Brake Complete,790,1 ea, EA
9503H79004	Right Brake Complete,790,1 ea, EA
9503H79005	Handle Bar Column,790,1 ea, EA

PARTS**Description**

9503H79006	Handle Bar Height Adj Lock,790,1ea, EA
9503H79007	Steering Column,790,1ea, EA
9503H79008	Brake Cable Guide,790,1ea, EA
9503H79011	Steering Fork,790,1ea, EA
9503H79012	Column To Fork Locking Bolt,790,1ea, EA
9503H79013	Fork Bearings,790,1pair, PR
9503H79014	Knee Pad,790,1ea, EA
9503H79023	Brake Cable & Hardware Only for 790, 1ea, EA