

Economy Dual Channel TENS



- For chronic, acute and post operative pain
- Designed for rehabilitation and physical therapy
- Dual isolated channels
- Mode Selector: Constant
- Intensity control is adjustable from 0-80 mA, 500 load
- Pulse frequency adjustable from 2-150 Hz
- Pulse width is adjustable from 40-250 micro-seconds
- Wave form: asymmetrical, bi-phase, square pulse
- Plastic belt clip
- Comes complete with carrying case, 2 lead wires, 9-volt battery, 1 package of 4 pre-gelled electrodes and operating manual

PRODUCTS	Description
AGF-3E	TENS Unit Econ. w/Case, 1/bx

SPECIFICATIONS	Description
Warranty	3 Year Limited