

Aloha Nasal Pillow System

HCPCS: A7034, A7035



- Not only is the Aloha Nasal Pillow System extremely lightweight and quiet, but it offers an incredibly stable and comfortable seal
- Ball-and-swivel joint minimizes torque and pull to accommodate various sleeping positions
- Exclusive Arced-Track technology allows the angle of the pillow reservoir to be adjusted for a more customized and natural fit
- CPAP hose can be routed over the head or worn down the chest
- Wide-set headgear provides limited points of contact on the face and no interference with the eyes
- Angled exhalation ports minimize potential disruption or discomfort for the user or bed partner
- Minimalistic headgear made of soft, pliable, and breathable material

PRODUCTS

Description

ALO100	Aloha Nasal Pillow System, EA
--------	-------------------------------

ACCESSORIES

Description

ALO121	Aloha Replmt Nasal Pillows - Sm, EA
ALO125	Aloha Replmt Nasal Pillows - Lg, EA
ALO123	Aloha Replmt Nasal Pillows - Med, EA
ALO130	Aloha Nasal Pillows HdGr, EA
ALO131	Aloha Nasal Pillows Frm - Sm, EA
ALO133	Aloha Nasal Pillows Frm - Med, EA
ALO135	Aloha Nasal Pillows Frm - Lg, EA