



## Shower Chair

The **Armstrong Shower Chair** is designed for individuals who require seating support and stability when bathing or showering. The chair features a comfortable contoured textured seat that can support patient weights up to 400 pounds. The legs have one-inch, easy-snap button adjustments and non-skid feet for added sturdiness and safety. Tool-free assembly and limited lifetime warranty are standard.

### FEATURES

- Aluminum frame is lightweight, durable and corrosion-proof
- Blow-molded seat has a textured surface and is contoured with drainage holes to reduce slipping
- Quick tool-free assembly
- Two-hole handles in seat provide security when sitting or standing
- Legs are height adjustable with one-inch snap button adjustments
- TAA Compliant, made in Taiwan
- Available individually, 1 per case
- HCPCS Code: E0245\*

### SPECIFICATIONS

Leg Base Footprint:	17.25" W x 15" D
Seat Dimensions:	19.5"W x 11"D
Adjustable Seat Height:	14.75" – 18.75"
Weight Capacity:	400 lb
Carton Dimensions:	AMS-DLXSC: 25" x 15" x 5" AMS-DLXSC-4: 28" x 21" x 19"
Shipping Weight:	AMS-DLXSC: 6.6 lb AMS-DLXSC-4: 25.0 lb
Product Weight:	5.0 lb
Warranty:	Limited Lifetime on frame One year on all parts

Item #	Description	UPC	UOM
AMS-DLXSC	Deluxe Aluminum Shower Chair with Removable Back	815067072637	1/cs
AMS-DLXSC-4	Deluxe Aluminum Shower Chair with Removable Back	815067072934	4/cs

\* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.