



Steerable Aluminum Knee Walker



The Armstrong Steerable Knee Walker offers increased comfort and ease of mobility to patients who cannot put weight on their foot or ankle, but want to remain active. With fast and simple tool-free assembly, the knee walker is a great alternative to crutches or folding walkers that can have adverse effects such as back or under arm pain from extended use.

With its four large treaded wheels and comfortable padded knee platform, the Knee Walker is perfect for indoor and outdoor use. The offset post on the padded knee platform allows for right-to-left adjustment in seconds — no tools required. The Knee Walker is equipped with an easy folding mechanism that collapses using a thumb release lever and fits beneath the raised knee platform for ease of storage and transportation.



Superior treaded tires and drum-braking system



Adjustable hand brake and handle bar

FEATURES

- Ideal for individuals recovering from an ankle or foot injury
- Generous turning radius for increased maneuverability
- Large 8-inch treaded wheels with sealed bearings
- Adjustable hand brake ensures safety
- Height adjustable to meet a variety of patient profiles
- Durable metal, removable basket included
- Premium padded knee platform for optimal patient comfort
- Highly effective drum brake maximizes user safety and stopping ability
- 350 lb weight capacity
- Limited lifetime warranty on frame; 90 days on all other parts
- TAA-complaint, made in Taiwan
- Individually packaged, 1/cs
- HCPCS: E0118*

SPECIFICATIONS

Patient Height Range:	4'11" - 6'6"
Knee Platform Dimensions:	7" (W) x 14" (L)
Adjustable Handle Height:	32" - 41"
Knee Platform Height:	17" - 23"
Product Weight:	18.8 lb
Weight Capacity:	350 lb.

Item #	Description	UPC	UOM
AMS-KNW	Steerable Aluminum Knee Walker, 350 lb Weight Capacity	815067072613	1/cs

* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.