



## Elevated Toilet Seat

The **Armstrong Elevated Toilet Seat** is a must for anyone with bending or sitting difficulties. Constructed of heavy duty-blow molded plastic that is durable and easy to clean, this seat comes complete with an easy front locking mechanism and comfortable attachable arms to assist in patient stability. Available in multiple heights, with and with out arms and locking options to accommodate all needs. All toilet seats are constructed of heavy duty blow-molded plastic that is durable and easy to clean.

### FEATURES

- Adjustable front locking mechanism tightens easily with one hand and fits most toilets and prevents shifting of seat when transferring
- Elevates existing seat height by 4.5"
- Padded armrests for extra support
- Extra-large hygienic hole cutout
- Quick, secure, and tool-free installation
- Color: white
- TAA compliant, made in Taiwan
- Individually packaged, 1 per case
- HCPCS Code: E0244\*

### SPECIFICATIONS

Seat Width:	17.5"
Seat Depth:	16"
Width Between Arms:	19.75"
Seat Height:	4.5"
Seat Opening:	9" W x 10" D'
Weight Capacity:	350 lb
Carton Dimensions:	21" x 8" x 18"
Shipping Weight:	7.3 lb
Warranty:	Limited lifetime on frame One year on all parts

Item #	Description	UPC	UOM
AMS-LETS	Armstrong Elevated Toilet Seat	815067072583	1/cs

\* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.