



Drop-Arm Commode

The **Armstrong Drop-Arm Commode** features a durable steel frame with an easy-release pull pin mechanism allowing for safe lateral patient transfers to and from the commode. The anti-slip, soft padded armrests provide extra comfort and support. The frame can be configured to be used as a bedside commode, toilet safety frame or raised toilet seat. The tool-free removable back makes configuring the commode easy for any application.

FEATURES

- Durable steel frame provides strength and stability
- Equipped with easy-to-release arm mechanism allowing for safe lateral transfers to and from the commode
- Pull pin drop-arm allows easier side transfer
- Soft, slip-resistant, padded covered armrests
- Seat height readily adjusts in one-inch increments
- Bucket sets easily in frame
- Removable backrest
- Comes complete with 7.5 quart commode bucket with carry handle and splash guard
- Color: white frame, gray accessories
- TAA compliant, made in Taiwan
- Individually packaged, one per case
- HCPCS Code: E0163*

SPECIFICATIONS

Width Between Arms:	20"
Outside Width:	23"
Seat Dimensions:	16" W x 15.5" D
Seat Opening:	8" W x 12.5" D
Adjustable Seat Height:	20" – 24"
Adjustable Arm Height:	26" – 30"
Adjustable Back Height:	28.5" – 32.25"
Bucket Dimensions:	12.5" W x 12.5" D x 6.5" H
Weight Capacity:	350 lb
Carton Dimensions:	21.5" x 16.5" x 14"
Shipping Weight:	19.5 lb
Warranty:	Limited lifetime on frame One year on all parts

Item #	Description	UPC	UOM
AMS-PDAC	Armstrong Drop-Arm Commode	815067072590	1/cs

* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.