



## Deluxe Aluminum Transport Chair

The Armstrong Transport Chair is a lightweight and durable aluminum construction design with a glossy black finish. The handy back-release hinges allow the back to fold down compact for easy storage and transport. Comes complete with seat belt, tool-free adjustable swing-away composite footrests and heel loops. Supports patient weights up to 300 pounds.



### FEATURES

- Lightweight aluminum construction (only 25 lbs)
- Handy back-release hinge allow the back to fold down flat and compact for easy storage and transport
- Padded full-length permanent armrests
- Padded nylon upholstery is comfortable, lightweight and easy to clean
- Seat belt for added safety
- Nylon seat straps allow for ease of folding and lifting
- Push-to-lock rear wheel locks
- Composite 8" wheels with solid rubber tires are lightweight and maintenance free
- Tool-free, adjustable swing-away composite footrests with heel loops included
- TAA compliant, made in Taiwan
- Individually packaged, 1 per case
- HCPCS Code: E1038\*

### SPECIFICATIONS

Overall Width:	
AMS-TC17:	20" (open)   9" (closed)
AMS-TC19:	22" (open)   9" (closed)
Width between Hand Grips:	
AMS-TC17:	16.25" (outer) 14" (inner)
AMS-TC19:	18.25" (outer) 16" (inner)
Overall Height:	39.5"
Overall Depth:	36" with riggings 24" without riggings
Seat-to-Floor Height:	19"
Back Height:	19"
Seat Dimensions:	
AMS-TC17:	17" (W) x 16" (D)
AMS-TC19:	19" (W) x 16" (D)
Patient Weight Capacity:	300 lb
Product Weight:	
AMS-TC17:	23.4 lb
AMS-TC19:	23.8 lb
Carton Dimensions:	27" x 10" x 32"
Shipping Weight:	
AMS-TC17:	27.4 lb
AMS-TC19:	27.8 lb
Warranty:	Limited lifetime on frame One year on all parts

Item #	Description	UPC	Seat Width	Color	Wt. Cap.	UOM
AMS-TC17	Deluxe Aluminum Transport Chair	815067072644	17"	Black	300 lb	1/cs
AMS-TC17-C	Deluxe Aluminum Transport Chair with 3" Foam Cushion		17"	Black	300 lb	1/cs
AMS-TC19	Deluxe Aluminum Transport Chair	815067072651	19"	Black	300 lb	1/cs
AMS-TC19-C	Deluxe Aluminum Transport Chair with 3" Foam Cushion		19"	Black	300 lb	1/cs

\* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.