

SmartDose® Mini Auto-Adjusting Oxygen Conservers

HCPCS:



- Delivers a uniform oxygen pulse at increments: 1, 2, 3, 4, 5 plus CF present @ 2 LPM
- Continuously monitors breathing patterns and adjusts oxygen dose to match activity level
- Achieves maximum conservation efficiency by automatically adjusting the dose down when less oxygen is needed
- Delivers a quiet, gentle dose of oxygen using an efficient and compact design
- Allows delivered dose to increase up to 2 settings during patient exertion up to a maximum pulse dose setting of 7
- Can use 2 AA Alkaline or 2 AA rechargeable batteries for a duration of about 1 year for most patients *
- Operating pressure range: 500 – 3000 PSI
- Weighs just 15.5 ounces with batteries
- Auto shut-off feature conserves battery life
- Fits easily in virtually any carry bag

PRODUCTS

Description

CTOX-MN02	SmartDose® Mini Auto-Adjusting Oxygen Conservers, EA
-----------	--

SPECIFICATIONS

Description

Warranty	Limited 2 Year
----------	----------------

ACCESSORIES

Description

SEAL100	Seal Washer,Conservers/Regulators,10/pkg, PKG
OM-900-FEC	E-Clip for Conservers Handles, 10/bag, BAG
OM-900-FH	S-Handle for OM-900/800/500 Series, 1/e, EA