



## Cruiser X4 Wheelchair Lightweight, Dual-Axle Wheelchair

HCPCS: K0004



INSET A



INSET B



INSET C



INSET D



INSET E



INSET F



INSET G



INSET H

- Weighs 34 pounds or less when properly equipped in the basic configuration (Excluding front rigging)
- Easy to propel and transport
- Frame style eliminates seat guides and allows for custom backs and accessories
- Flip-back, removable arm is height adjustable from top of seat rail 10" to 14" and features a front mounted easy-to-use lever release mechanism (Figure A - C)
- Padded armrests provide added patient comfort
- Height adjustable segmented back
- Aluminum and Carbon steel blended frame with combination powder coated and silver vein finish is attractive, chip resistant and easy-to-maintain
- Durable, flame retardant nylon upholstery resists mildew and bacteria
- Patented adjustable depth seat 16" and 18" on chair depths (Figure D - F)
- 3 position quick release rear axles mounted on Composite Mag-style wheels with composite hand rims are lightweight and maintenance free (Figure G)
- Polyurethane rear tires and composite casters offer superior performance, smoother ride and are lightweight
- Dual over center cross links provide added rigidity (Figure H)
- Comes with swing-away footrests or elevating leg rests with tool-free adjustable length riggings
- Composite footplates with heel loops are durable and lightweight
- Precision sealed wheel bearings through-out ensure long lasting and reliability
- 8" front casters have 3 height adjustments and angle adjustments (Figure I)
- Triple rear axle positions provides easy transition of seat height from standard to hemi level
- Comes with push-to-lock wheel locks

### PRODUCTS

### Description

CX416ADDA-ELR	CruiserX4 16 Adj DeskArm ElvLg, 1/ea
CX416ADDA-SF	CruiserX4 16 Adj DeskArm SwgFt, 1/ea
CX416ADFA-ELR	CruiserX4 16 Adj FullArm ElvLg, 1/ea
CX416ADFA-SF	CruiserX4 16 Adj FullArm SwgFt, 1/ea
CX418ADDA-ELR	CruiserX4 18 Adj DeskArm ElvLg, 1/ea
CX418ADDA-SF	CruiserX4 18 Adj DeskArm SwgFt, 1/ea
CX418ADFA-ELR	CruiserX4 18 Adj FullArm ElvLg, 1/ea
CX418ADFA-SF	CruiserX4 18 Adj FullArm SwgFt, 1/ea
CX420ADDA-ELR	CruiserX4 20 Adj DeskArm ElvLg, 1/ea
CX420ADDA-SF	CruiserX4 20 Adj DeskArm SwgFt, 1/ea
CX420ADFA-ELR	CruiserX4 20 Adj FullArm ElvLg, 1/ea
CX420ADFA-SF	CruiserX4 20 Adj FullArm SwgFt, 1/ea

### SPECIFICATIONS

### Description

Width Open with 16"   18"   20"   Seat	23"   25"   27"
Width Closed	12"

SPECIFICATIONS	Description
Height	36"
Seat-to-Floor Height	17"   18"   19"
Back Height	16" - 17.5"
Seat Depth	16"/18" On Chair
Length with Riggings	43"
Product Weight without Riggings (16", 18", 20")	32 lbs.   34 lbs.   36 lbs.
Carton Shipping Weight (16", 18", 20")	46 lbs.   48 lbs.   50 lbs.
Weight Limit	300 lbs.
Warranty on Frame	Limited Lifetime
Warranty	Limited Lifetime

ACCESSORIES	Description
STDS860	Universal Forward Stabilizer, 1/pr
STDS850	Seat Belts, Auto-Clasp Type Closure, 1/ea
STDS851	Seat Belts, Hook and Loop Fastener Type Closure, 1/ea
STDS855	Seat Belts, Auto-Clasp Type Closure, Bariatric, 1/ea
STDS856	Seat Belts, Hook and Loop Fastener Type Closure, Bariatric, 1/ea
STDS1002	Standard 8" Caster, 1/ea
STDS834	Overhead Anti-Theft Device Single Pole, 1/ea
STDS806	Anti-Fold Bar, 1/ea
STDS1034	Wheelchair Cane/Crutch Holder, 1/ea
STDS831	Heel Loops, 1/bx
WASR	Limb Support, Right, 1/ea
WASL	Limb Support, Left, 1/ea
STDS601	Wheel Lock Extension, 1/ea
STDSSHK-K3	Super Hemi-Kit, 1/ea
STDSK234-PTL	Pull to Lock Brakes(left), 1/ea
STDSK234-PTR	Pull to Lock Brakes(right), 1/ea

PARTS	Description
STDS860	Universal Forward Stabilizer, 1/pr
STDS819	Anti-Tippers With Wheels, 1/pr
STDS850	Seat Belts, Auto-Clasp Type Closure, 1/ea
STDS851	Seat Belts, Hook and Loop Fastener Type Closure, 1/ea
STDS855	Seat Belts, Auto-Clasp Type Closure, Bariatric, 1/ea
STDS856	Seat Belts, Hook and Loop Fastener Type Closure, Bariatric, 1/ea
STDS806	Anti-Fold Bar, 1/ea
STDS1034	Wheelchair Cane/Crutch Holder, 1/ea
STDS831	Heel Loops, 1/bx
WASR	Limb Support, Right, 1/ea
WASL	Limb Support, Left, 1/ea
STDS835	Chart Carry Pocket, 1/ea
STDS601	Wheel Lock Extension, 1/ea
STDSSHK-K3	Super Hemi-Kit, 1/ea
STDS4SC4X-15	Back Cane, 1/ea
STDS312403	Hand Grip, 1/ea
STDS4SC4X-10	Back Upholstery 16", 1/ea
STDS4SC4X-11	Back Upholstery 18", 1/ea
STDS4SC4X-12	Back Upholstery 20", 1/ea
STDS3J2416-AD	Seat Upholstery 16", 1/ea
STDS3J2418-AD	Seat Upholstery 18", 1/ea
STDS3J2420-AD	Seat Upholstery 20", 1/ea
STDS3J2411	Upholstery Strap, 1/ea
STDS3J24X6-AD	Cross Frame 16", 1/ea
STDS3J24X8-AD	Cross Frame 18", 1/ea