

IntelliPAP® Standard CPAP System



INSET A



INSET B



INSET C



INSET D

- Whisper-quiet operation – 26 dBA!
- Exceptionally small footprint fits easily into the home environment and carry-on luggage
- Includes discreet carry bag that can be easily expanded to hold additional accessories (Figure A) which makes it great for travelers
- Optional integrated heated humidification system easily docks into the bottom of the unit without requiring a separate hose or power cord (Figure B-C)
- Auto on/off automatically starts the flow of air after breathing into the mask and stops the flow of air when mask is removed
- Delay pressure ramp eases the machine to the prescribed pressure – allowing the user to fall asleep at a lower pressure (0-45 minutes in 5-minute increments and 5 increments for pressure)
- SmartCode® compliance tracking encodes usage data that can be communicated to a clinician
- Optional SmartLink® Therapy Management System provides detailed patient therapy data, including oximetry (Figure D)

PRODUCTS

	Description
DV51D	IntelliPAP® Standard CPAP - Unit only, E0601, 1/ea
DV51D-HH	IntelliPAP® Standard and Heated Humidification System, E0601, E0562, 1/ea

SPECIFICATIONS

	Description
Dimensions	4.2" (H) x 6.5" (W) x 6.9" (D) - unit only;
Dimensions	6.4" (H) x 6.5" (W) x 8.4" (D) with humidifier
Weight	2.7 lbs - unit only; 4.45 lbs with humidifier
Pressure Range	3-20 cm H ₂ O
Automatic Altitude Adjustment	Up to 9,000 feet
Auto On/Off	Yes (Can be enabled or disabled)
SmartCode® Feature	Yes
SmartLink® Compatible	Yes
Unit warranty	3 years
Humidifier warranty	3 years
Chamber warranty	3 months

ACCESSORIES

	Description
7351D-616	6 ft. 22mm Tubing, Autoclave, A7037, 1/ea
7353D-601	Oxygen Adapter, A9999, 1/ea
DV51D-602	Air Inlet Filter Pkg (4-Pk) Gray, A7038, 1/ea
DV51D-603	Fine Particle Aire Filter Pkg (4-Pk) White, A7038, 1/ea
DV51D-604	Air Supply Port Plug, 1/ea
DV51D-605	Heater Connection Cover, 1/ea
DV51D-606	US Power Cord, A9999, 1/ea
DV51D-610	Carrying Case DV5X Series, A9999, 1/ea

ACCESSORIES

Description

DV51D-615	IntelliPAP Serial Cable (for direct connect), 1/ea
DV51D-619	IntelliPAP DC Power Cord, A9999, 1/ea
DV51D-620	IntelliPAP Outlet Cap, A9999, 1/ea
DV51D-626	Air Inlet Filter (2-Pk) Gray, A7038, 1/ea
DV51D-627	Fine Particle Air Filter (2-Pk) White, A7038, 1/ea
DV51D-631	Bacteria Filter/Elbow/Tubing Kit, 1/ea
DV51D-682	IntelliPAP Disinfection Kit, 1/ea
DV51D-683	DV5 Hospital Kit, 1/ea
DV51D-691	USB to Serial Adapter (for direct connect), 1/ea
DV51D-694	SmartLink USB Extension Cable, 1/ea
DV51D-696	IntelliPAP 12V Adapter with Clips, A9999, 1/ea
DV5C	IntelliPAP Humidifer Chamber, A7046, 1/ea
DV5C-614	Chamber Sealing Gasket, A9999, 1/ea
DV5HH	IntelliPAP Heated Humidifier, E0562, 1/ea
DV5HT	IntelliPAP Replacement Heated Tubing, A4604, 1/ea
DV5HTK	IntelliPAP Heated Tube Kit, A4604, 1/ea
DV5HTS	IntelliPAP Heated Tube On/Off Switch, 1/ea
DV5M	SmartLink Module, 1/ea
DV5M-FC-1	SmartLink Module with Data Card, 1/ea

PARTS

Description

DV51D-602	Air Inlet Filter Pkg (4-Pk) Gray, 1/ea
DV51D-603	Fine Particle Air Filter Pkg (4-Pk) White, 1/ea
DV51D-626	Air Inlet Filter (2-Pk) Gray, 1/ea
DV51D-627	Fine Particle Air Filter (2-Pk) White, 1/ea
DV51D-610	Carrying Case DV5X Series, 1/ea
DV51D-606	US Power Cord, 1/ea
DV51D-607	Power Cord Europe, 1/ea
DV51D-608	Power Cord UK, 1/ea
DV51D-609	Power Cord Australia, 1/ea
DV51D-619	IntelliPAP DC Power Cord, 1/ea
DV51D-620	IntelliPAP Outlet Cap, 1/ea
DV51D-696	IntelliPAP 12V Adapter with Clips, 1/ea
DV51D-604	Air Supply Port Plug, 1/ea
DV51D-629	6 ft. 22 mm Tubing, 1/ea
7353D-601	Oxygen Adapter, 1/ea
7351D-616	6 ft. 22mm Tubing, Autoclave, 1/ea
DV5HTK	IntelliPAP Heated Tube Kit, 1/ea
DV5HT	IntelliPAP Replacement Heated Tubing, 1/ea
DV5HTS	IntelliPAP Heated Tube On/Off Switch, 1/ea
DV5HH	IntelliPAP Heated Humidifier, 1/ea
DV5C	IntelliPAP Humidifer Chamber, 1/ea
DV5C-614	Chamber Sealing Gasket, 1/ea
DV51D-605	Heater Connection Cover, 1/ea
DV5SS2	SmartLink 2.0 Software, Data Card & Reader, 1/ea
DV5M-SS-2	SmartLink Module, 2.0 Software, Data Card & Reader, 1/ea
DV5FC-1	SmartLink Data Card, 1/ea
DV5FC-10	SmartLink Data Card (PK/10), 1/ea
DV5M	SmartLink Module, 1/ea
DV5CR	SmartLink Data Card Reader, 1/ea
DV51D-615	IntelliPAP Serial Cable (for direct connect), 1/ea
DV51D-690	ABOB Interface, 1/ea
DV51D-691	USB to Serial Adapter (for direct connect), 1/ea
DV51D-693	SmartLink USB Extension Cable modem, 1/ea