



DeVilbiss IntelliPAP® Bilevel S System

Once again, you talked. We listened.
Everyone gets what they need.

The IntelliPAP Bilevel S is the perfect choice for patients who require higher pressures to maintain airway patency. With a pressure range of 3-25 cm H₂O and the ability to adjust the IPAP and EPAP delta up to 22 cm H₂O, the IntelliPAP Bilevel S offers maximum adjustability for operation in CPAP and Bilevel modes.

The IntelliPAP Bilevel S system offers our patented Flow Rounding® technology to provide additional comfort during pressure transition. This helps to better mimic a patient's natural breathing pattern. IntelliPAP is built for years of dependable performance and patient satisfaction. It's just what you would expect from DeVilbiss.



DV55D Series

Performance

- **Quality** – Less than 1% failure rate and backed by an industry-leading 3-year warranty
- **Compliance** – Free remote tracking of compliance data through SmartCode® Technology
- **Comfort** – Optimal humidification and Flow Rounding Technology
- **Design** – Stacked design creates an exceptionally small footprint and whisper-quiet operation (26 dBA)*

DeVilbiss IntelliPAP Bilevel S System

Design

- Easy and intuitive setup
- Auto on/off** automatically starts flow of air with patient breath
- Eight-foot power cord provides greater positioning flexibility – no bulky power brick
- Reminders alert patient to replace or check parts
- Includes attractive travel bag
- Large, slip-resistant pads on bottom of unit to minimize movement on nightstand



Humidifier Chamber

- Overfill protection
- Corrosive resistant, stainless steel plate
- Unique baffling increases air water exchange for optimal humidity
- 400 ml capacity for 8 hours of use at the highest setting




Usage Data	
Last 90 Days	
% Days of at least 4 hours (not continuous) in a 24 hour period	84%
Day Count	90
Days at least 4 hours	76
A/EI	5.5
High Leak Flow Time	2%
NRI	1.75
EPI	2
SmartCode	H8ZJT-L1D9-7DRAB
Adherence Score (4 continuous hours)	67% (2000)

SmartCode

- Encrypted usage data codes that can be communicated by phone, e-mail or text to generate compliance and clinical reports

Selectable IPAP and EPAP Flow Rounding*** (0-5)



The Aloha™ Nasal Pillow System is the perfect fit with IntelliPAP



DeVilbiss Healthcare | 100 DeVilbiss Drive | Somerset, PA 15501 | USA
800.338.1988 | 814.443.4881 | www.DeVilbissHealthcare.com