



DeVilbiss IntelliPAP® AutoBilevel System

Comfortable therapy at higher pressures.

The IntelliPAP AutoBilevel system incorporates the comfort of AutoAdjust® and Bilevel technology to deliver comfortable therapy at higher prescription pressures. The IntelliPAP AutoBilevel offers our patented Flow Rounding® technology to provide additional comfort during pressure transition. This helps to better mimic a patient's natural breathing pattern.

The IntelliPAP AutoBilevel system can be used for conventional CPAP, AutoAdjust, standard Bilevel and AutoBilevel. The system allows clinicians to adjust the apnea and hypopnea detection to fit individual patient standards. The split-night feature, similar to a single-night lab study and titration, allows patient activity to be measured through the first part of the night without pressure changes.



DV57D Series

Performance

- **Quality** – Less than 1% failure rate and backed by an industry-leading 3-year warranty
- **Compliance** – Free remote tracking of compliance data through SmartCode® Technology
- **Comfort** – Optimal humidification and Flow Rounding Technology
- **Design** – Stacked design creates an exceptionally small footprint and whisper-quiet operation (26 dBA)*

DeVilbiss IntelliPAP AutoBilevel System


Design

- Easy and intuitive setup
- Auto on/off** automatically starts flow of air with patient breath
- Eight-foot power cord provides greater positioning flexibility – no bulky power brick
- Reminders alert patient to replace or check parts
- Includes attractive travel bag
- Large, slip-resistant pads on bottom of unit to minimize movement on nightstand



Humidifier Chamber

- Overfill protection
- Corrosive-resistant, stainless steel plate
- Unique baffling increases air water exchange for optimal humidity
- 400 ml capacity for 8 hours of use at the highest setting




SmartCode

- Encrypted usage data codes that can be communicated by phone, e-mail or text to generate compliance and clinical reports



The Aloha™ Nasal Pillow System is the perfect fit with IntelliPAP

Selectable IPAP and EPAP Flow Rounding*** (0-5)



DeVilbiss Healthcare | 100 DeVilbiss Drive | Somerset, PA 15501 | USA
800.338.1988 | 814.443.4881 | www.DeVilbissHealthcare.com