



WE HELP YOU BREATHE EASIER.

Inside this catalog, you'll find a wide selection of respiratory products for all your patients – from neonates to adults. Our portfolio spans the continuum from acute to long-term care.

Improve Clinical Outcomes

We are constantly developing products that can help improve outcomes in respiratory care.

Professional Insight.

Knowledge and insight from fellow respiratory professionals help us continuously improve our products.

Business and Value Oriented

Our wide selection includes a variety of products for every aspect of good respiratory care, so you can customize ordering and maximize savings.




Table Of Contents

● Oxygen	4-13
● Asthma Management & Medication Delivery	14-22
● Airway Clearance	23-33
● CPAP/BIPAP	34-41
● Filtration & Humidification	42-45
● Resuscitation	46-49
● Patient Monitoring & Accessories	50-54
● Oxygen Regulators & Carts	55-59

OXYGEN

Cannulas	5
Oxygen Tubing	6
Bubble Humidifiers	7
Accessories	7
Oxygen Masks	8
Aerosol Masks	9
Blenders	10
Medical Gas Hose Assemblies	11
Oxygen and Air Flowmeters	12
Neonatal, Perinatal and Specialty Flowmeters	13

SuperSoft Nasal Cannulas

Medical device-related pressure ulcers (MDRPU) have been receiving greater attention. MDRPUs tend to progress rapidly, as they often occur over areas with minimal fatty tissue. Findings indicate that patients that use a medical device are up to four times more likely to develop a pressure ulcer of any kind.¹

Even short term use of an oxygen cannula can lead to the development of pressure ulcers on the ears and cheeks. The Medline SuperSoft nasal cannula is made of an advanced material that is soft, flexible and non-tacky. Also features curved, non-flared nasal prongs. You will notice a difference right away:

- » Maintain skin integrity – Reduce pressure and friction on the ears and cheeks.
- » Improve patient compliance with added comfort – Patients are less apprehensive about wearing a cannula that feels soft and comfortable.
- » Reduce costs – No need for the additional time and expense of padding cannulas with foam ear cushions.
- » Not made with natural rubber latex.

Adult

Item No.	Description	Connector	Pkg
HCSU4514S	7' (2.1 m) Tubing	Universal	50/cs
HCSS4514S	7' (2.1 m) Tubing	Standard	50/cs
HCS4510S	No Tubing, (Cannula Harness only)	Male	50/cs
HCSU4504S	4' (1.2 m) Tubing	Universal	50/cs
HCSU4516S	14' (4.3 m) Tubing	Universal	50/cs
HCSU4515S	25' (7.6 m) Tubing	Universal	25/cs

Pediatric

Item No.	Description	Connector	Pkg
HCSU4518S	Pediatric, 7' (2.1 m) Tubing	Universal	50/cs
HCSU4517S	Infant, 7' (2.1 m) Tubing	Universal	50/cs
HCSU4521S	Preemie, 7' (2.1 m) Tubing	Universal	50/cs
HCSU4520S	Neonatal, 7' (2.1 m) Tubing	Universal	50/cs

Soft-Touch Cannulas

- » Curved, non-flared nasal prongs
- » Available in a wide variety of sizes and lengths
- » Crush-resistant tubing
- » Not made with natural rubber latex

Adult

Item No.	Description	Connector	Pkg
HCS4514	7' (2.1 m) Tubing	Standard	50/cs
HCS4510	Cannula Harness w/o Tubing	Male	50/cs
HCS4504B	4' (1.2 m) Tubing	Standard	50/cs
HCS4516B	14' (4.3 m) Tubing	Standard	50/cs
HCS4515B	25' (7.6 m) Tubing	Standard	25/cs
HCSU4514	7' (2.1 m) Tubing	Universal	50/cs
HCSU4504B	4' (1.2 m) Tubing	Universal	50/cs
HCSU4516B	14' (4.3 m) Tubing	Universal	50/cs
HCSU4515B	25' (7.6 m) Tubing	Universal	25/cs

Pediatric

Item No.	Description	Connector	Pkg
HCS4518	7' (2.1 m) Tubing	Standard	50/cs
HCSU4518	7' (2.1 m) Tubing	Universal	50/cs

Infant

Item No.	Description	Connector	Pkg
HCS4517	7' (2.1 m) Tubing	Standard	50/cs
HCSU4517	7' (2.1 m) Tubing	Universal	50/cs



OXYGEN



Crush-Resistant Tubing
This multi-channel interior feature helps to maintain a strong airflow



Standard Cannulas

- » Crush-resistant tubing with standard and universal oxygen connectors
- » Not made with natural rubber latex

Item No.	Description	Connector	Pkg
HCS4511B	Straight Prongs, 7' (2.1 m) Tubing	Standard	50/cs
HCSU4511B	Straight Prongs, 7' (2.1 m) Tubing	Universal	50/cs



HCS4511B

High-Flow SuperSoft Cannulas

- » Made of an exceptionally soft material for increased patient comfort
- » Adult liter flow from 6LPM to 15LPM
- » Pediatric liter flow up to 7LPM
- » For use with high-flow bubble humidifiers
- » Not made with natural rubber latex

Item No.	Description	Pkg
HCS4514HF	Adult, 7' (2.1 m) Tubing, Standard Connector	25/cs
HCS4516BHF	Adult, 14' (4.3 m) Tubing, Standard Connector	25/cs
HCS4515HF	Adult, 25' (7.6 m) Tubing, Standard Connector	10/cs
HCS4518HF	Pediatric, Harness w/2" Connector	10/cs
HCS4517HF	Infant, Harness w/2" Connector	10/cs



HCS4514HF

Oxygen Tubing

- » Unique, relaxed, crush resistant, low-memory material helps prevent kinking
- » Selected products in violet color to provide greater visual contrast to assist in patient safety
- » For standard and high flow applications
- » Not made with natural rubber latex

Item No.	Description	Connector	Pkg
HCS4507	7' (2.1 m) Oxygen Tubing	Standard	50/cs
HCS4507V	7' (2.1 m) Oxygen Tubing, Violet	Standard	50/cs
HCS4524	14' (4.3 m) Oxygen Tubing	Standard	50/cs
HCS4524V	14' (4.3 m) Oxygen Tubing, Violet	Standard	50/cs
HCS4525	25' (7.6 m) Oxygen Tubing	Standard	25/cs
HCS4525V	25' (7.6 m) Oxygen Tubing, Violet	Standard	25/cs
HCS455010	50' (15.2 m) Oxygen Tubing	Standard	10/cs
HCS4550VB	50' (15.2 m) Oxygen Tubing, Violet	Standard	10/cs
HCSU4507	7' (2.1 m) Oxygen Tubing	Universal	50/cs
HCSU4507V	7' (2.1 m) Oxygen Tubing, Violet	Universal	50/cs
HCSU4524	14' (4.3 m) Oxygen Tubing	Universal	50/cs
HCSU4524V	14' (4.3 m) Oxygen Tubing, Violet	Universal	50/cs
HCSU4525	25' (7.6 m) Oxygen Tubing	Universal	25/cs
HCSU4525V	25' (7.6 m) Oxygen Tubing, Violet	Universal	25/cs
HCSU455010	50' (15.2 m) Oxygen Tubing	Universal	10/cs
HCSU4550VB	50' (15.2 m) Oxygen Tubing, Violet	Universal	10/cs



HCS4507



HCS4525V

Medline Bubble Humidifiers

Choose from a wide variety of sizes and configurations.

Prefilled

Item No.	Description	Pkg
HCS00350	Prefilled Bubble Humidifier, 350 ml, with Adaptor	20/cs
HCS00550	Prefilled Bubble Humidifier, 550 ml, with Adaptor	12/cs
HCS00640	Prefilled Bubble Humidifier, 640 ml, with Adaptor	10/cs

Large Volume Nebulizer

HCS03128	O28 Nebulizer Adaptor, F102 range 28-98%	50/cs
HCS04000	1,070 ml Sterile Water Bottle	10/cs
HCS04128	Prefilled Aerosol Humidifier, 1070 ml, with O28 Nebulizer Adaptor	10/cs
M-X1065	1,000 ml Sterile Water Bottle	12/cs

Empty Bubble Humidifier

HCS64377	350 ml Empty Reusable 6 PSI Bubble Humidifier	50/cs
VTB6900	350 ml Empty Reusable Hi-Flow Bubble Humidifier	10/cs



Accessories

Medline's complete selection includes a wide variety of everyday connectors and adaptors.

Item No.	Description	Pkg
HCS0440	EZ Wrap Cushions for Cannula	50/cs
HCS0590	Swivel Connector for Oxygen Tubing	25/cs
HCS1642	Pressure Line Adapter, 5-7 mm ID	50/cs
HCS2503	Disposable Water Trap	25/cs
HCS4555	Green Christmas Tree Adaptor with Swivel	50/cs
HCS65100	Rigid Connector for Oxygen Tubing	200/cs
HCS02KIT1	Oxygen Setup Kit	10/cs
HSK962A	Circuit Connector, 22 OD/15 ID 22 OD	50/cs
HCS974	Aerosol Drainage Bag with "T" Adapter Hanger	50/cs
HCS974Y	Aerosol Drainage Bag with "Y" Adapter Hanger	50/cs
WST219	Aerosol Tee Adaptor	50/cs
B-F64019	Nose Clip	100/cs
AGISH749	Y Connector	5/ea
HCS676	15" Humidifier Connector, Tubing Adapter	50/cs
HCS1790	21" Humidifier Connector, Tubing Adapter	50/cs
HCS8986	Humidifier Connector with Rigid Black Elbow Connector	1/ea





Oxygen Masks

Choose from an industry-leading selection of masks for all your oxygen-delivery needs.

- » Soft, anatomical form for patient comfort
- » Crush-resistant tubing
- » A secure fit with elastic straps
- » Convenient choice of universal or standard connectors
- » Not made with natural rubber latex

Adult

Item No.	Description	Connector	Pkg
HCS4600B	Adult, Medium-Concentration Mask with 7' (2.1 m) Tubing	Standard	50/cs
HCSU4600B	Adult, Medium-Concentration Mask with 7' (2.1 m) Tubing	Universal	50/cs
HCS4610B	Adult, High-Concentration Mask with 7' (2.1 m) Tubing	Standard	50/cs
HCSU4610B	Adult, High-Concentration Mask with 7' (2.1 m) Tubing	Universal	50/cs
HCS4640B	Adult, Non-Rebreather Mask with 7' (2.1 m) Tubing	Standard	50/cs
HCSU4640B	Adult, Non-Rebreather Mask with 7' (2.1 m) Tubing	Universal	50/cs
HCS120E	Adult, Venturi Mask with 6 Jets and 7' (2.1 m) Tubing	Standard	50/cs
HCSU120E	Adult, Venturi Mask with 6 Jets and 7' (2.1 m) Tubing	Universal	50/cs
HCS121E	Adult, Venturi Mask with Single Dial and 7' (2.1 m) Tubing	Standard	50/cs
HCSU121E	Adult, Venturi Mask with Single Dial and 7' (2.1 m) Tubing	Universal	50/cs



Pediatric

Item No.	Description	Connector	Pkg
HCS4601B	Pediatric, Medium-Concentration Mask with 7' (2.1 m) Tubing	Standard	50/cs
HCSU4601B	Pediatric, Medium-Concentration Mask with 7' (2.1 m) Tubing	Universal	50/cs
HCS4642B	Pediatric, Non-Rebreather Mask with 7' (2.1 m) Tubing	Standard	50/cs
HCSU4642B	Pediatric, Non-Rebreather Mask with 7' (2.1 m) Tubing	Universal	50/cs



Universal Connector



Standard Connector



HCS4630B



HCS4630



HCS4632



HCS245E

Aerosol Masks

Extra patient comfort comes from an anatomical form and a soft, light and comfortable vinyl material.

- » A secure fit with elastic straps
- » Standard 22 mm connector
- » Not made with natural rubber latex

Adult


Item No.	Description	Pkg
HCS4630B	Adult Elongated Aerosol Mask	50/cs
HCS4630	Adult Elongated Aerosol Mask with Swivel Connector	50/cs
HCS4632	Adult Face Tent with 22 mm Connector	50/cs
HCS245E	Adult Trach Mask	50/cs

Pediatric

HCS4631B	Pediatric Elongated Aerosol Mask	50/cs
HCS4621B	Pediatric Trach Mask	50/cs

Corrugated Tubing

HCS4541	Corrugated Tubing, 6" (15.2 cm), Blue	50/cs
HSK600	Corrugated, Segmented Tubing, 100' (30.5 m), Opaque with cuttable cuff every 6" (15.2 cm)	1/cs
HSK626	Corrugated, Segmented Tubing, 6' (1.8 m), Opaque	50/cs



HSK600

Blenders

Our blenders require very little maintenance and last for years.

- » Compact size
- » Warranty: two years (manufacturer)



GTBH High-Low Flow Air/Oxygen Blender GTBL Low Flow Air/Oxygen Blender

DD - Left Port	EE - Right Port	F- Bottom Port	GG - Accessories
NN = None	NN = None	NN = None	NN = None
NI = NIST	NI = NIST	1 = NIST	01 = Blender Wall Mount
DI = DISS	DI = DISS	2 = DISS	02 = Blender Mount with Pole Clamp
01 = 0-1 LPM Flowmeter	01 = 0-1 LPM Flowmeter		03 = Blender Mount with Eclipse® or Fairfield Equipment Rail Bracket
02 = 0-200 cc Flowmeter	02 = 0-200 cc Flowmeter		04 = Blender Mount with Unimount Equipment Rail Bracket
03 = 0-3 LPM Flowmeter	03 = 0-3 LPM Flowmeter		WT = Water Trap
04 = 0-15 LPM Flowmeter	04 = 0-15 LPM Flowmeter		
05 = 0-12 LPM Flowmeter	05 = 0-12 LPM Flowmeter		
06 = 10-70 LPM Flowmeter	06 = 10-70 LPM Flowmeter		

EXAMPLE: GTHBDI04202 is a high flow blender, left port DISS, right port 1-15 LPM flowmeter with DISS bottom port and wall mount with pole clamp.

Multiple configurations available. Contact your Medline Representative for more information.

Medical Gas Hose Assemblies

Medline offers custom-configured medical gas hose assemblies in five easy steps:

- » Select hose type
- » Select gas service
- » Select the inlet fitting (gas service and type)
- » Select the outlet fitting (gas service and type)
- » Specify the hose length

GT = Hose Assembly

X = Hose Type (1/4" ID for pressure gases, 5/16" ID for vacuum)

- A = 1/4" ID Non-Conductive
- C = 5/16" ID Non-Conductive

AAA = Gas Service

- OXY = Oxygen
- VAC = Vacuum
- N2O = Nitrous Oxide
- INA = Instrument Air
- AIR = Medical Air
- NIT = Nitrogen
- CO2 = Carbon Dioxide
- WAG = Waste Anesthetic Gas Disposal

B = Inlet Fitting (Connects to the Source)


- A = 1/8" Male Pipe Thread
- C = 1/4" Male Pipe Thread
- E = DISS Female Hex Nut
- G = Chemetron® Adapter
- I = Puritan-Bennet Adapter
- K = Schrader Adapter, Swivel
- B = 1/8" Female Pipe Thread
- D = 1/4" Female Pipe Thread
- F = DISS Female Hand-tight
- H = Ohmeda® Adapter
- J = Schrader Adapter, Non-Swivel

C = Outlet Fitting (Connects to the Medical Device)

- A = 1/8" Male Pipe Thread
- C = 1/4" Male Pipe Thread
- E = DISS Male with Check
- G = DISS Female Hand-tight
- I = DISS Male, NO CHECK, with Grip
- K = Ohmeda Coupler
- M = Schrader Coupler
- B = 1/8" Female Pipe Thread
- D = 1/4" Female Pipe Thread
- F = DISS Female Hex Nut
- H = DISS Male, NO CHECK
- J = Chemetron Coupler
- L = Puritan-Bennet Coupler

DDD = Hose Length, in feet (e.g. 036, 048, 144, etc.)

EXAMPLE: GTAOXYFE060 1/4" ID non-conductive oxygen hose assembly with a DISS female hand-tight inlet and a DISS male with check outlet, 5' in length.

















Oxygen and Air Flowmeters

Antimicrobial*, environmentally friendly, and exclusively distributed nationwide by Medline.

- » Antimicrobial and environmentally friendly additive in the rotary knobs*
- » For both oxygen and medical air
- » Thoroughly tested for accuracy
- » Durable polycarbonate hood and metering tube
- » Available in nickel plated solid brass or aluminum
- » Stainless steel ball float
- » Warranty: ten years (manufacturer)

Flowmeter Ordering Matrix

FLOWMETER TYPE ADAPTER						
 Standard Inlet 1/8" NPT Female	GTFM187A15L00	GTFM187A1500P	GFMD187A15L00	GTFM197A15L00	GTFM197A1500P	GFMD197A15L00
 Ohmeda® Adapter	GTFM187A15LOH	GTFM187A150HP	GFMD187A15LOH	GTFM197A15LOH	GTFM197A150HP	GFMD197A15LOH
 Chemetron® Adapter	GTFM187A15LCH	GTFM187A15CHP	GFMD187A15LCH	GTFM197A15LCH	GTFM197A15CHP	GFMD197A15LCH
 DISS Hand-Tight Nut and Nipple Adapter	GTFM187A15LDH	GTFM187A15DHP	GFMD187A15LDH	GTFM197A15LDH	GTFM197A15DHP	GFMD197A15LDH
 DISS Nut and Nipple Adapter	GTFM187A15LDS	GTFM187A15DSP	GFMD187A15LDS	GTFM197A15LDS	GTFM197A15DSP	GFMD197A15LDS
 Puritan-Bennett Adapter	GTFM187A15LPB	GTFM187A15PBP	GFMD187A15LPB	GTFM197A15LPB	GTFM197A15PBP	GFMD197A15LPB
 Oxequip® Medstar Adapter	GTFM187A15LOE	GTFM187A150EP	GFMD187A15LOE	GTFM197A15LOE	GTFM197A150EP	GFMD197A15LOE
 Oxequip Twist Adapter	GTFM187A15LOQ	GTFM187A150QP	GFMD187A15LOQ	GTFM197A15LOQ	GTFM197A150QP	GFMD197A15LOQ

Options

- » Flowmeters listed above are aluminum. For brass, change the A to a B (example: GTFM187A15LOH to GTFM187B15LOH)
- » Dual flowmeter with 2 power take-offs, also available

Specifications

Weight (No Adapter): Aluminum Body: 4.24 oz/122 gm
 Brass Body: 6.64 oz/188 gm

Increment: 0-6 LPM in ½ LPM increment,
 6-15 LPM in 1 LPM Increment

Flush Flow Rate: 65 LPM (depending on supply flow)

Max Inlet Pressure: 100 psig

PTO Flow Rate: Over 100 LPM (depending on supply flow)








* This product does not protect users or others against disease-causing bacteria.

Neonatal, Perinatal and Specialty Flowmeters

The performance and accuracy you need for critical care applications for neonate through adult patients.

- » 200 cc and 1 LPM models are designed for perinatal and neonatal applications where extreme accuracy under very low flow is required
- » Available for oxygen and air
- » Acrylic block provides a strong flow tube
- » Maximum inlet pressure: 100 PSI
- » Outlet: DISS male (9/16"-18UNF)
- » Calibrated at 500 PSI
- » Warranty: five years (manufacturer)



ADAPTER \ FLOWMETER TYPE	200 CC Oxygen	1 LPM Oxygen	3 LPM Oxygen	16 LPM Oxygen	75 Oxygen LPM	16 LPM Air	75 LPM Air
 1/8" NPT Female	GTFM293200CXX	GTFM29301XX	GTFM29303XX	GTFM29316LXX	GTFM29375LXX	GTFM28316LXX	GTFM28375LXX
 Chemetron® Adapter	GTFM293200CCH	GTFM29301CH	GTFM29303CH	GTFM29316LCH	GTFM29375LCH	GTFM28316LCH	GTFM28375LCH
 DISS Hex	GTFM293200CDS	GTFM29301DS	GTFM29303DS	GTFM29316LDS	GTFM29375LDS	GTFM28316LDS	GTFM28375LDS
 DISS Hand-Tight	GTFM293200CDH	GTFM29301DH	GTFM29303DH	GTFM29316LDH	GTFM29375LDH	GTFM28316LDH	GTFM28375LDH
 Ohmeda® Adapter	GTFM293200COH	GTFM29301OH	GTFM29303OH	GTFM29316LOH	GTFM29375LOH	GTFM28316LOH	GTFM28375LOH
 Puritan-Bennett Adapter	GTFM293200CPB	GTFM29301PB	GTFM29303PB	GTFM29316LPB	GTFM29375LPB	GTFM28316LPB	GTFM28375LPB
 Oxequip® Medstar Adapter	GTFM293200COE	GTFM29301OE	GTFM29303OE	GTFM29316LOE	GTFM29375LOE	GTFM28316LOE	GTFM28375LOE

ASTHMA MANAGEMENT & MEDICATION DELIVERY

Aeromist Compact™ Nebulizer Compressor	15
Aeromist Buddies Nebulizer Compressors	16
Aeromist Compact Colors™ Nebulizer Compressors	18
Aeromist Trek Portable Nebulizer Compressor	19
Handheld Nebulizers	20
MDI Spacers	21
Valved Tee Adapters	21
Peak Flow Meters	22
Spirometers	22



Aeromist Compact™ Nebulizer Compressor


Our new, smallest nebulizer compressor is less bulky and will take up less of your valuable space.

- » Smallest foot print of our compressors
- » Low center of gravity
- » Four non-skid rubber feet
- » 7 ft power cord
- » Warranty: 5 years

Item No.	Description	Pkg
HCS70004	Aeromist Compact Nebulizer Compressor with Disposable Nebulizer kit	8/cs
HCS70004F	Replacement Filters	5/pk

Specifications

Method of nebulization: dry piston compressor
 Electrical Requirements: 120VAC/60Hz
 Maximum nebulizer capacity: 10 ml
 Sound noise level: ≤ 60 dBA



Aeromist Buddies™ Nebulizer Compressors

Kid-friendly nebulizers help put young patients at ease, while still providing safe and effective aerosol treatments. Fun-loving Monkey is sure to get them in the swing of things. Cute and calming Clownfish reels in patients with its inquisitive eyes, gently flowing fin and perky tail. The Doggy looks paws-itively fetching with its friendly eyes and red tongue. The Turtle calms kids' fears with his slow pace and silly face.

- » Kit includes a disposable nebulizer medicine cup, 7-ft. (2.13 m) of tubing, and mouthpiece, and pediatric mask
- » Durable hard plastic shell
- » Rugged compressor rated for continuous use
- » Compact and stable, with a small footprint, low center of gravity and four non-skid rubber feet
- » Warranty: 5 years

MonkeyNeb Nebulizer Compressor

Item No.	Description	Pkg
HCSMONKEY	MonkeyNeb Nebulizer Character Compressor	8/cs
HCS70004F	Replacement Filters	5/pk

FishNeb Nebulizer Compressor

- » Grab-N-Go backpack provides convenient storage

Item No.	Description	Pkg
HCSFISHNEB	FishNeb Nebulizer Character Compressor	8/cs
HCS70004F	Replacement Filters	5/pk

DoggyNeb Nebulizer Compressor

- » Grab-N-Go backpack provides convenient storage

Item No.	Description	Pkg
HCSDOGNEB2	DoggyNeb Nebulizer Character Compressor	8/cs
HCS70004F	Replacement Filters	5/pk

Turtle Nebulizer Compressor

Item No.	Description	Pkg
HCSTURTLE	Turtle Nebulizer Character Compressor	8/cs
HCS70004F	Replacement Filters	5/pk



NEW BUDDY



Specifications
Method of nebulization: dry piston compressor
Electrical Requirements: 120 VAC/ 60Hz
Maximum nebulizer capacity: 10 ml
Sound noise level: \leq 60 dBA



Aeromist Compact Colors™ Nebulizer Compressor

The Aeromist Colors nebulizer compressor provides patients with a more fun and stylish alternative than the clinical-looking models of the past.

- » Kits include a disposable nebulizer medicine cup, 7-ft. (2.13 m) of tubing and a mouthpiece
- » Color choices: electric blue, neon green, neon pink, or sophisticated charcoal
- » Durable hard plastic shell
- » Rugged compressor rated for continuous use
- » Compact and stable, with a small footprint, low center of gravity and four non-skid rubber feet
- » Warranty: 5 years

**BREAK AWAY
FROM THE BLAND!**

Aeromist Colors™
Nebulizer Compressors
Combo Pack HCS60004DP
includes 2 of each color.



Item No.	Description	Pkg
HCS70004G	Neon Green Compressor and Nebulizer Kit	8/cs
HCS70004BL	Electric Blue Compressor and Nebulizer Kit	8/cs
HCS70004P	Neon Pink Compressor and Nebulizer Kit	8/cs
HCS70004C	Charcoal Compressor and Nebulizer Kit	8/cs
HCS70004DP	Sample Pack, Two of each color	8/cs
HCS70004F	Replacement Filters	5/pk

Specifications

Method of nebulization: dry piston compressor
 Electrical Requirements: 120VAC/60Hz
 Maximum jet nebulizer capacity: 10 ml
 Sound noise level: ≤ 60 dBA

Aeromist Trek Portable Nebulizer Compressor

Better compliance with easy treatments, anytime and anywhere.

- » Suitable for travel with a sleek design, extra small size, light weight and battery operation
- » Durable hard plastic shell
- » Compact and stable, with a small footprint
- » Kit includes a disposable nebulizer medicine cup, 7-ft. (2.13 m) of tubing, mouthpiece, rechargeable battery, carry case, AC and DC adapters, adult and pediatric masks, and 5 replacement filters
- » Warranty: 3 years for the compressor, 6 months for the battery



Item No.	Description	Pkg
HCS2GO	Portable Battery-Operated Nebulizer Compressor Kit	each
HCS2GOBATT	Replacement Battery for Portable Nebulizer Compressors (Items HCS2GO and HCS2GONB)	each
HCS2GOF	Replacement filters	5/pk
HCS2GOAC	Replacement AC Adapter	each
HC22GDC	Replacement Car Charger	each

Handheld Nebulizer Kits

Medline offers a wide selection of reusable (up to 6 months use) and disposable (up to 2 weeks use) nebulizer kits for asthma medication delivery.

- » Ideal particle size of less than five microns, increasing the amount of medication delivered deep within lungs
- » Standard fittings accommodate most clinical applications
- » Not Made with Natural Rubber Latex
- » Medicine can be added without removing the cap
- » Kits include a medicine cup, 7-ft. (2.1 m) of kink-resistant extension tubing, and either a mask or a mouthpiece with a T adapter
- » Packaged 50/cs, except WST0382 is packaged 12/cs

DISPOSABLE

Item #	Upstream	Side Stream	Tee	Mouthpiece	Reservoir Tube	7-ft Tubing	Adult Mask	Pediatric Mask	Pacifier
HCS4482	●		●	●		●			
HCS4483	●		●	●	●	●			
HCS4485	●					●	●		
HCS4486	●					●		●	
HCS4492		●	●	●		●			
HCS4493		●	●	●	●	●			
HCS4495		●				●	●		
WST0382	●					●			●
WST0312	●					●		Character	
WST0214	●								

REUSABLE

Item #	Upstream	Side Stream	Tee	Mouthpiece	Reservoir Tube	7-ft Tubing	Adult Mask	Pediatric Mask	Pacifier
HCS4482R	●		●	●		●			
HCS4483R	●		●	●	●	●			

Did You Know?

The surface area of a person's lungs is roughly the same size as one side of a tennis court.



Medline MDI Spacers

Medline spacers are designed for patients who have difficulty with metered dose inhalers (MDIs). With spacers, more medication is deposited in the lungs and less is deposited in the mouth and throat.

- » A simple and economical way to improve effectiveness of inhaler
- » Easy to use and helps avoid some of the unpleasant taste of inhaled medication
- » Anti-static and compatible with any standard mask with a 22 mm fitting
- » Masks are made of a soft silicone and fit snugly to ensure optimal drug delivery
- » The collapsible model is portable and convenient

Item No.	Description	Pkg
HCSA100	MDI Spacer Only	10/bx
HCSA200	MDI Spacer Only, Collapsible	10/bx
HCSA112	MDI Spacer with Infant Mask	10/bx
HCSA113	MDI Spacer with Pediatric Mask	10/bx
HCSA114	MDI Spacer with Adult Mask	10/bx



HCSA200 is collapsible for portability.



Valved Tee Adapters

- » For inline medication delivery
- » 22 mm outer/inner diameter

Item No.	Description	Pkg
HCSA1008	Valved Tee Adapter, 22 mm OD/OD	30/cs
HCSA1009	Valved Tee Adapter, 22 mm OD/ID	30/cs



Peak Flow Meters

Peak flow meters help patients manage their asthma and gauge the effectiveness of their treatment program, which may help them predict their next attack.

- » Individually calibrated to ensure accuracy
- » Durable for years of reliable performance
- » The pocket model offers superb convenience and portability
- » Available in low range and full range versions – low range models are better for geriatric, pediatric or other patients with diminished flow rates



Peak Flow Meter Color Zones²

Green Zone

80- 100% of normal peak flow rate. Asthma is reasonably under control. Continue prescribed management program.

Yellow Zone

50- 80% of normal peak flow rate. Airways are narrowing and may require extra treatment. Alert physician.

Red Zone

Less than 50% of normal peak flow rate. Severe airway narrowing signal. Immediate decisions and actions need to be taken. Alert physician.

Item No.	Description	Pkg
HCSPF100	Medline Portable Peak Flow Meter, Full Range, 60-720 lpm	12/bx
HCSPF200	Standard Peak Flow Meter, Full Range, 60-850 lpm	12/bx
HCSPF210	Standard Peak Flow Meter, Low Range, 30-400 lpm	12/bx
HCSPF220	Disposable Pediatric Mouthpiece, One-Way Valve	20/bx
HCSPF262	Disposable One Way Valve Mouthpiece, Adult	100/cs

Spirometers

Help patients perform and monitor post-surgical breathing exercises.

- » One-way valve ensures patients inhale rather than exhale
- » Highly visible pistons and universal graphics indicate correct inspiratory flow rate
- » Two volume capacities – 2500ml and 5000ml

Item No.	Description	Pkg
HCS93250	Incentive Spirometer, 2500 ml	12/cs
HCS93500	Incentive Spirometer, 5000 ml	12/cs



AIRWAY CLEARANCE

Open Suction Catheters	24
Sleeved Catheters	24
Kits and Trays	25
Suction Catheter Kits w/ Sterile Solutions	26
Sterile Solutions	26
Closed Suction Catheters	27
Mucus/Specimen Traps	28
Suction Aspirator	28
Phoneris/Aeroflex System	29
TrachVox® Speaking HME	30
Tracheostomy Ties	30
Tracheostomy Care and Cleaning Trays	31
Suction Regulators	32

Open Suction Catheters

Medline's D419 series catheters offers a full range of straight packaging.

- » Available in a complete range of continuous adult and pediatric sizes
- » All Medline suction catheters feature easy-to-read depth markings
- » DeLee tip is flat, soft and smooth to minimize airway trauma for newborns, infants and children
- » Angled, soft and smooth whistle tip maximizes suctioning capability for adult patients, while minimizing trauma potential

Single Sterile-Straight Pack


Item No.	Description	Pkg
DYND41906	Suction Catheter, Single Sterile: 6 Fr	100/cs
DYND41908	Suction Catheter, Single Sterile: 8 Fr	100/cs
DYND41900	Suction Catheter, Single Sterile: 10 Fr	100/cs
DYND41901	Suction Catheter, Single Sterile: 12 Fr	100/cs
DYND41902	Suction Catheter, Single Sterile: 14 Fr	100/cs
DYND41903	Suction Catheter, Single Sterile: 16 Fr	100/cs
DYND41904	Suction Catheter, Single Sterile: 18 Fr	100/cs

Poly-Cath

DYND40990	Suction Catheter, Poly-Cath, 10 Fr	100/cs
DYND40992	Suction Catheter, Poly-Cath, 14 Fr	100/cs

Continuous Valve

DYND41960	Suction Catheter, Continuous Valve, 10 Fr	100/cs
DYND41962	Suction Catheter, Continuous Valve, 14 Fr	100/cs
DYND41964	Suction Catheter, Continuous Valve, 18 Fr	100/cs



DYND41900C



DYND41902

Sleeved Catheters

Our D407 series open-suction catheters are contained in a sleeve to help avoid cross-contamination.

Item No.	Description	Pkg
DYND40708F	Sleeved Suction Catheter: 8 Fr	50/cs
DYND40700F	Sleeved Suction Catheter: 10 Fr	50/cs
DYND40701F	Sleeved Suction Catheter: 12 Fr	50/cs
DYND40702F	Sleeved Suction Catheter: 14 Fr	50/cs



DYND40701F



Kits and Trays

Medline's D409 series kits and trays provide all the necessities for airway clearance.

- » Irrigation solutions can be poured into the tray or the pop-up cup
- » All kits and trays contain Medline Aloetouch® 3G gloves
- » Catheters have easy-to-read depth markings
- » DeLee tip is flat, soft and smooth to minimize airway trauma for newborns, infants and children
- » Angled, soft and smooth whistle tip maximizes suctioning capability for adult patients, while minimizing trauma potential

Suction Catheter Trays

Item No.	Description	Pkg
DYND40986	Suction Catheter Mini Tray with Gloves: 6 Fr	100/cs
DYND40988	Suction Catheter Mini Tray with Gloves: 8 Fr	100/cs
DYND40980	Suction Catheter Mini Tray with Gloves: 10 Fr	100/cs
DYND40981	Suction Catheter Mini Tray with Gloves: 12 Fr	100/cs
DYND40982	Suction Catheter Mini Tray with Gloves: 14 Fr	100/cs
DYND40983	Suction Catheter Mini Tray with Gloves: 16 Fr	100/cs
DYND40984	Suction Catheter Mini Tray with Gloves: 18 Fr	100/cs

Suction Catheter Kits

Item No.	Description	Pkg
DYND40976	Suction Catheter Pack with gloves and Pop-up cup, 6 Fr	50/cs
DYND40978	Suction Catheter Pack with gloves and Pop-up cup, 8 Fr	50/cs
DYND40970	Suction Catheter Pack with gloves and Pop-up cup, 10 Fr	50/cs
DYND40971	Suction Catheter Pack with gloves and Pop-up cup, 12 Fr	50/cs
DYND40972	Suction Catheter Pack with gloves and Pop-up cup, 14 Fr	50/cs
DYND40973	Suction Catheter Pack with gloves and Pop-up cup, 16 Fr	50/cs
DYND40974	Suction Catheter Pack with gloves and Pop-up cup, 18 Fr	50/cs

Suction Catheter Trays, Poly-cath

Item No.	Description	Pkg
DYND48982	Suction Catheter Mini Tray with Gloves, Poly-Cath: 14 Fr	100/cs

MEDLINE MANUFACTURES
OVER 123,000
 sterile procedure kits *every day.*

Suction Catheter Kits with Sterile Solutions

Our D414 series suction catheters come with everything needed for airway clearance, including a choice of irrigation fluids.

- » All kits and trays contain Medline Aloetouch 3G gloves
- » Catheters have easy-to-read depth markings
- » DeLee tip is flat, soft and smooth to minimize airway trauma for newborns, infants and children
- » Angled, soft and smooth whistle tip maximizes suctioning capability for adult patients, while minimizing trauma potential

Water

Item No.	Description	Pkg
DYND41476	Suction Catheter Tray with Sterile Water, 100 ml.; 6 Fr	36/cs
DYND41478	Suction Catheter Tray with Sterile Water, 100 ml.; 8 Fr	36/cs
DYND41470	Suction Catheter Tray with Sterile Water, 100 ml; 10 Fr	36/cs
DYND41471	Suction Catheter Tray with Sterile Water, 100 ml; 12 Fr	36/cs
DYND41472	Suction Catheter Tray with Sterile Water, 100 ml; 14 Fr	36/cs
DYND41474	Suction Catheter Tray with Sterile Water, 100 ml; 18 Fr	36/cs

Saline

Item No.	Description	Pkg
DYND41478SA	Suction Catheter Tray with Sterile Saline, 100 ml; 8 Fr	36/cs
DYND41470SA	Suction Catheter Tray with Sterile Saline, 100 ml; 10 Fr	36/cs
DYND41442	Suction Catheter Tray with Sterile Saline, 100 ml; 14 Fr	36/cs
DYND40442	Suction Catheter Tray with Saline and Drape: 14 Fr	36/cs



DYND41472



DYND40442



Sterile Solutions

These convenient, ready-to-use containers reduce the contamination risks and potential waste associated with bulk containers of sterile water and saline.

- » Ideal for device irrigation, suction, dressing moistening and wound debridement
- » Cup size provides the proper amount for suctioning

Water

Item No.	Description	Pkg
RDI30295	Sterile Water, 100 ml	48/cs
DYND40570	Sterile Water, 110 ml	48/cs
PCS1550	Sterile Water, 250 ml	24/cs
DYND500MLW	Sterile Water, 500 ml	18/cs

Saline

Item No.	Description	Pkg
RDI30296	Sterile Saline, 0.9%, 100 ml	48/cs
DYND40540	Sterile Saline, 0.9%, 110 ml	48/cs
PCS1650	Sterile Saline, 0.9%, 250 ml	24/cs
DYND500MLS	Sterile Saline, 0.9%, 500 ml	18/cs



RDI30295

RDI30296

DYND500MLW

DYND500MLS

Closed Suction Catheters

Medline's closed suction catheters provide enhanced patient care, while saving staff time.

- » Soft tip helps minimize damage to mucous membranes
- » Transparent connector enables easy observation
- » High-quality, translucent EVA sleeve helps reduce cross-contamination and facilitates catheter advancement
- » One-way needleless irrigation
- » Easy-to-read depth markings
- » Color-coded locking valves and day label
- » Compatible with all modes of ventilation




Pediatric

Elbow

Item No.	Description	Pkg
DYNCP5	Closed Suction Catheter with Elbow: 5 Fr	10/cs
DYNCP6	Closed Suction Catheter with Elbow: 6 Fr	10/cs
DYNCP7	Closed Suction Catheter with Elbow : 7 Fr	10/cs
DYNCP8	Closed Suction Catheter with Elbow: 8 Fr	10/cs

Y Connector

Item No.	Description	Pkg
DYNCPY5	Closed Suction Catheter with "Y": 5 Fr	10/cs
DYNCPY6	Closed Suction Catheter with "Y": 6 Fr	10/cs
DYNCPY7	Closed Suction Catheter with "Y": 7 Fr	10/cs
DYNCPY8	Closed Suction Catheter with "Y": 8 Fr	10/cs



Adult

Double Swivel

Item No.	Description	Pkg
DYNCS12	Closed-Suction Catheter, Double Swivel, ETT: 12 Fr	10/cs
DYNCS14	Closed-Suction Catheter, Double Swivel, ETT: 14 Fr	10/cs
DYNCS12T	Closed-Suction Catheter, Double Swivel, Trach: 12 Fr	10/cs
DYNCS14T	Closed-Suction Catheter, Double Swivel, Trach: 14 Fr	10/cs

T-Piece

Item No.	Description	Pkg
DYNCS12T	Catheter, "T" piece, Closed Suction, ETT: 12 Fr	10/cs
DYNCS14T	Catheter, "T" piece, Closed Suction, ETT: 14 Fr	10/cs
DYNCS12T	Catheter, "T" piece, Closed Suction, Trach: 12 Fr	10/cs
DYNCS14T	Catheter, "T" piece, Closed Suction: Trach: 14 Fr	10/cs

Mucus/Specimen Traps

Our mucus and specimen traps have leak-proof transport caps for more secure specimen collection.

- » Sterile plastic vials, graduated in 1cc increments
- » Includes threaded cap, tubing not made from natural rubber latex and patient ID label

Sterile DeLee Mucus Trap

Item No.	Description	Pkg
DYND44108	Trap with catheter and cap, capacity 20cc: 8 Fr	50/cs
DYND44110	Trap with catheter and cap, capacity 20cc: 10 Fr	50/cs

Sterile Mucus Specimen Trap

Item No.	Description	Pkg
DYND44140	Trap With Transport Cap, Capacity: 40 cc	50/cs
DYND44180	Trap With Transport Cap, Capacity: 80 cc	50/cs



Suction Aspirator and Disposable Canister Kits

Our easy-to-use aspirators are ultra-quiet.

- » Weight: 12.1-lb./5.49 kg
- » Maximum suction: -21 inHg/520 mmHg
- » Rugged construction is ideal for bedside carts
- » Disposable 850 cc or 1200 cc suction canister sets with float lids.
- » HCS7851 is a kit that also includes suction tubing (1 each of 6' and 18" length).
- » Replacement lids are not available


Item No.	Description	Pkg
HCS7000	Suction Aspirator	1/ea
HCS7850	Canister with Float Lid: 850 cc	12/cs
HCS7851	Canister with Float Lid and Tubing: 850 cc	12/cs
HCS7852	Canister with Filter Lid and Tubing: 850 cc	50/cs
HCS7853	Canister with Filter Lid: 850 cc	100/cs
HCS7120	Canister with Float Lid: 1,200 cc	12/cs



Tubing included with HCS7851 and HCS7852




Positive-Click jaws hold tight to help prevent disconnect. Releases easily with a twist



Longer, **easy-grip jaws** make the inner cannula easier to remove

Dual "ramps" make a positive-click connection with the ventilator connector jaws



Smaller corrugations, tapered inner walls and flush-fit with cannula let inline suction catheters pass with ease

360 degree swivel for better movement and less torsion on the inner cannula

Phoneris[®] Inner Cannula and Aeroflex[®] Ventilator Circuit Connector

The Phoneris/Aeroflex system uses an innovative design to comfortably and securely connect tracheostomy patients to ventilator circuits.

- » **Secure** – A tight connection helps avoid accidental disconnects
- » **Gentle** – Swivel enables comfortable movement and reduces torque on the stoma
- » **Easy** – Puts an end to improvised, unapproved securing devices

Item No.	Description	Pkg
HCSTIC4A	Phoneris Inner Cannula; Jackson 4	10/cs
HCSTIC6A	Phoneris Inner Cannula; Jackson 6	10/cs
HCSTIC8A	Phoneris Inner Cannula; Jackson 8	10/cs
HCSTIC10A	Phoneris Inner Cannula; Jackson 10	10/cs
HCSTIC4681A	Phoneris AeroFlex	50/cs

TrachVox® Speaking HME

TrachVox HME helps spontaneously breathing tracheostomy patients verbally communicate. When depressed, the HME changes the airway flow to make speaking possible. Additionally, it heats and humidifies inspired air and delivers supplemental oxygen through an integrated port.

- » Allows for force inspiration, expiration and suctioning
- » Low flow oxygen can be attached directly to the TrachVox oxygen port
- » Uses standard 15 mm connector

Item No.	Description	Pkg
EMS3131010	Trachvox Speaking valve with HME	30/cs



EMS3131010

Tracheostomy Ties

These soft, comfortable straps have adjustable hook-and-loop closures to keep tracheostomy tubes securely in place, ensuring uninterrupted airflow.

Premium Tracheostomy Ties with Technology from Pepper Medical

- » Soft, hypoallergenic, no stretch cotton ensures patient safety and comfort while reducing skin breakdown
- » Size range to fit neonate to adult patients

Item No.	Description	Pkg
HCS240	Tracheostomy Tube Holder, One Size Fits Most	10/cs
HCS241	Tracheostomy Tube Holder, Pediatric	10/cs
HCS242	Tracheostomy Tube Holder, Fits Neonate to Infant	10/cs

Standard Medline Tracheostomy Ties

- » Available in two sizes to fit infants, young children and adults

Item No.	Description	Pkg
ZPP301P	Tracheostomy Tube Holder, 1-Piece Pediatric	20/cs
ZPP501P	Tracheostomy Tube Holder, 2-Piece Pediatric	20/cs
ZPP501	Tracheostomy Tube Holder, 2-Piece Adult	20/cs
ZPP301	Tracheostomy Tube Holder, 1-Piece Adult	20/cs
ZPP401	Anti-Disconnect Tracheostomy Tube Holder	20/cs



Did You Know?

We can build a custom kit to your exact specifications. Medline is a leading manufacturer of custom procedure kits.

Tracheostomy Care and Cleaning Trays

Medline offers a wide variety of adult and pediatric tracheostomy care and cleaning trays in standard and customized configurations to meet diverse clinical needs.

- » Available in economy, basic and full-compartment kits
- » May be customized

Contents	DYND40582	DYND40622	DYND40610	DYND40620	DYND40589	DYND40615
Tray Compartments	3	1	3	1	3	3
Cotton Tip Applicator	2	2	5	5	2	5
Non-woven Gauze 4" x 4"	4	2	2	2	2	2
Twill Tape 34"	1	1	1	1	1	1
Aloetouch Gloves	2		2	2	2	2
Drape	1		1	1	1	1
Pipe Cleaners		2	2			
Pliable Cleaning Brush	1	1	1		1	1
6-ply Tracheostomy Dressing	1		1	1	1	1
Hydrogen Peroxide					2	2
Saline, 100ml					1	
Water, 100ml						
Suction Catheter, 14 Fr						
Pkg	20/cs	30/cs	20/cs	20/cs	20/cs	20/cs

Contents	DYND40580	DYND40650	DYND40511	DYND40594	CC3T3030A	CC4681A	CC3T4691A
Tray Compartments	3	3	3	3	1	3	3
Cotton Tip Applicator	2	2	4	2	2	3	2
Non-woven Gauze 4" x 4"	2	2	1	3	1	3	1
Twill Tape 34"	1	1	1		1	1	1
Aloetouch Gloves	2	2	2	2		2	2
Drape	1		1	1			1
Pipe Cleaners				3	3	3	3
Pliable Cleaning Brush	1		1	1	1	1	1
6-ply Tracheostomy Dressing	1	1		1	1	1	1
Hydrogen Peroxide		1		1			
Forceps						1	
Saline, 100ml				1			
Water, 100ml		1	1				
Suction Catheter, 14 Fr	1			1			
Pkg	20/cs	20/cs	20/cs	20/cs	30/cs	20/cs	20/cs



DYND40582



DYND40589

Suction Regulators

Exclusively distributed nationwide by Medline.

- » Antimicrobial and environmentally friendly additive in the rotary knobs and selector switches*
- » Extremely durable
- » Continuous-only and continuous/intermittent designs
- » Available with a variety of inlet and outlet fittings and suction trap assemblies
- » Warranty: ten years (manufacturer)

User-friendly mode selector

Large, user-friendly mode selector switch and a 1.25" (3.2 cm) diameter adjustment knob

Large diameter gauge and Easy-to-read numbers

Large 2.5" (6.4 cm) diameter gauge with glow-in-the-dark scale and arrow


Color-coded sticker

Color-coded to help facilities manage inventory and allocation across different departments




Antimicrobial additive

Rotary knobs and selector switches are made with an antimicrobial additive*



Continuous/Intermittent (881 Series)

- » Available with ranges of 0-160 mmHg and 0-300 mmHg
- » Settings include continuous, off and intermittent
- » Factory default settings for on/off cycles are 16 seconds (± 2) on, 8 seconds (± 2) off; on/off settings can be adjusted without opening the case (controls are not accessible without tools)



Continuous Only (882 Series)




























- » Available with ranges of 0-300 mmHg and 0-760 mmHg
- » Settings include regulated, off and full

*This product does not protect users or others against disease-causing bacteria.

Suction Regulator Ordering Matrix

Instructions: Select suction regulator, adapter then fitting (steps 1-3).
Combine codes to create model number, then write chosen sticker color. See example below.

STEPS:	1 SUCTION REGULATOR	2 ADAPTER	3 FITTING	MODEL NO.	COLOR				
Example:	GT881VR160	+	OH	+	T	=	GT881VR160OHT	in	YL
Please Write Your Selection:		+		+		=		in	

1 SUCTION REGULATORS	2 ADAPTERS		3 FITTINGS		STICKER COLOR OPTIONS
 <p>Continuous / Intermittent 0-160 GT881VR160</p>	<p>NO ADAPTER</p> <p>1/8" NPT Female</p> <p>OO</p>	 <p>Chemetron Adapter GT3103VACNT2 CH</p>	<p>NO FITTING</p> <p>1/8" NPT Female</p> <p>A</p>	 <p>Hose Barb GT880VR304 D</p>	<ul style="list-style-type: none">  Blank (White)  YL (Yellow)  BL (Blue)  RD (Red)  PL (Purple)  PC (Peach)  PK (Pink)  GY (Gray)  GN (Green)  BB (Baby Blue) <p>* Color must be specified in sales notes.</p>
 <p>Continuous / Intermittent 0-300 GT881VR300</p>	 <p>Ohio® Adapter GT3111UVACNT2 OH</p>	 <p>DISS Hand-Tight Adapter GT3222UVACNT2 DH</p>	 <p>Trap Fitting GT881VR14K B</p>	 <p>DISS Male GTDISSV18NPT C</p>	
 <p>Continuous Only 0-300 GT882VR300</p>	 <p>Puritan-Bennett Adapter GT3104VACNT2 PB</p>	 <p>DISS Nut-Nipple Adapter GT3102VACNT2 DS</p>	 <p>Complete Suction Trap with Trap Fitting GT880VT5 T</p>		
 <p>Continuous Only 0-760 GT882VR760</p>	 <p>Oxequip Twist Adapter GT3311VACNT2 OQ</p>	 <p>Oxequip Medstar Adapter GT3301VACNT2 DE</p>	 <p>Complete Suction Trap with DISS GT880VT6 S</p>	 <p>Complete Suction Trap with DISS Hand-Tight GT880VT7 H</p>	

ACCESSORIES	Suction Traps				
	 <p>Suction Trap with No Fitting on Top GT880VT0</p>	 <p>Suction Trap with Male Trap Fitting Only GT880VT1</p>	 <p>Suction Trap with DISS Female Nut and Nipple GT880VT2</p>	 <p>Suction Trap with DISS Female Hand-Tight GT880VT3</p>	 <p>Suction Trap with Tubing Nipple GT880VT4</p>
 <p>Vacuum Adapter DISS to Hose Barb GTDISSMVHCHK</p>	 <p>Extension Adapter 2-1/2" Long GTEXTADPT18</p>	 <p>Security Cable GTFMS2</p>			

CPAP/BIPAP

CPAP System	35
AIRGel™ CPAP Masks	36
Silicone CPAP Masks	38
CPAP Accessories	39
NIV Masks	40



CPAP200

Continuous Positive Airway Pressure (CPAP) by Medline

In obstructive sleep apnea treatments, Medline's CPAP unit delivers precise, reliable therapy to your sleep apnea patients.

Description	Value
Power Supply	12V DC
Current Consuming	Max 1.5A
Air Pressure Range	3 cmH ₂ O-20 cmH ₂ O
Ramp Time	10-60 min
Air Pressure Calibrate	± 3 cmH ₂ O
Working Flow	3 ± 0.3m ³ /h @ 9 cm H ₂ O
Noise	< 33 dB @ 9 cmH ₂ O (1m)
Working Humidity	10% RH to 90% RH
Working Temperature	5°C to 36°C (41°F to 104°F)
Storage Humidity	< 95% RH
Storage Temperature	-20°C to 60°C (-4°F to 140°F)
Dimension (LxWxH)	144 x 120 x 91 mm
Weight	460g
Life of Use	2 years
Safety Criterion	CE, FDA



CPAP210

CPAP220

Item No.	Description
CPAP200	CPAP w/ SD Card and Heated Humidifier
CPAP210	CPAP w/ SD Card
CPAP200F	Filter Replacement (CPAP200), 5 pk
CPAP200PC	Replacement Cords (CPAP200)
CPAP220	Humidifier (CPAP210)
CPAP200SD	Replacement SD Card

AIRgel™ CPAP Masks

Amazing comfort and customized fit technology. The patented AIRgel replaceable cushioning minimizes discomfort and skin irritation while effectively sealing airflow. It is easily detached and reattached for easy cleaning, saving time and money.

- » Proprietary AIRgel cushioning minimizes discomfort and skin irritation while effectively sealing airflow
- » The flexible and formable mask provides a custom fit—flex, fit and forget!

Aura Nasal Mask

- » The shapeable shell of each mask can be easily molded to fit individual needs
- » Flexible spacebar is designed to provide superior stability without touching the skin
- » Active headgear connector moves with the patient to help maintain an effective seal

Aura-Nasal

Item No.	Description	Mask Size	Pkg
CPAP50176	Nasal Mask With Adjustable Gel Cushion	Small	ea
CPAP50177	Nasal Mask With Adjustable Gel Cushion	Medium	ea
CPAP50178	Nasal Mask With Adjustable Gel Cushion	Large	ea

Mojo® Full-Face Gel Mask

- » Full-face design in four sizes
- » Adjustable tension wheel minimizes discomfort while ergonomic cushions adapt to unique facial features for a comfortable fit
- » Magnetic headgear clips hold firmly and disconnects with one hand easy on-and-off

Item No.	Description	Headgear Style	Mask Size	Pkg
CPAP50835	Vented	One Size	Small	ea
CPAP50836	Vented	One Size	Medium	ea
CPAP50837	Vented	One Size	Large	ea
CPAP50838	Vented	One Size	XL	ea
CPAP55045	Headgear only	One Size	One Size	ea



iQ Blue® Nasal Gel Mask

- » One-size nasal mask design fits a wide variety of patients
- » Clearest field of vision possible with no forehead pad required – patients can keep their glasses on while reading or watching TV
- » Easy on-and-off thanks to quick release, breathable headgear

iQ Blue Nasal Gel

Item No.	Mask Style	Headgear Style	Mask Size	Pkg
CPAP50655	Vented	Stable Fit	One Size	ea
CPAP55044	Headgear Only	Stable Fit	One Size	ea



CPAP50655

MiniMe® 2 Pediatric Gel Mask

- » Nasal mask in two sizes, designed for children ages 2 to 12
- » Touchless spacebar eliminates the need for any forehead contact
- » Soft, breathable headgear for all-night comfort
- » Available in vented and non vented designs
- » Adjustable cushion with memory wire for a true custom fit.

Item No.	Description	Headgear Style	Mask Size	Pkg
CPAP60217	Vented w/Leak Port	Small	Small/Medium	ea
CPAP60218	Vented w/Leak Port	Large	Medium/Large	ea
CPAP60257	Non-Vented	Small	Small/Medium	ea
CPAP60258	Non-Vented	Large	Medium/Large	ea



CPAP60217

Custom Fit Technology


The shapable shell mask can be easily molded to fit individual needs.



Silicone CPAP Masks

A consistent seal without excessive tightness, delivering the prescribed therapy with enough comfort to help improve compliance.

- » OmniClip and silicone forehead pad for a custom fit at the bridge of the nose for patients of all head shapes and sizes
- » SilentVent elbow flushes out exhaled air quietly, without an additional inline exhalation port
- » Headgear attaches with convenient clips, includes a handy wash bag to keep it clean and fresh
- » Cushions, forehead pads, elbows, cones and clips are individually replaceable
- » Available in four sizes, with and without headgear




Item No.	Size	Description	Pkg
CPAP9400	Small	Full Face Mask with Head Strap » SilentVent Anti-Asphyxia Elbow, 22mm Male » OmniClip with Silicone Forehead Pad	1/ea
CPAP9401	Medium		
CPAP9402	Large		
CPAP9403	X-Large		
CPAP9404	Small	Full Face Mask, No Head Strap » SilentVent Anti-Asphyxia Elbow, 22mm Male » OmniClip with Silicone Forehead Pad	1/ea
CPAP9405	Medium		
CPAP9406	Large		
CPAP9407	X-Large		
CPAP9408	Small	Nasal Mask with Head Strap » SilentVent Elbow, 22mm Male » Silicone Forehead Pad » OmniClip with Silicone Forehead Pad	1/ea
CPAP9409	Medium		
CPAP9410	Large		
CPAP9411	X-Large		
CPAP9412	Small	Nasal Mask, No Head Gear Strap » SilentVent Elbow, 22mm Male » Silicone Forehead Pad » OmniClip with Silicone Forehead Pad	1/ea
CPAP9413	Medium		
CPAP9414	Large		
CPAP9415	X-Large		

CPAP Accessories

Nasal Pillows

Item No.	Description	Pkg
SUHCM012	Breeze Nasal Pillows, Includes: » One Small, Medium and Large Pillow » Headgear	ea
SUHCM011LP	Nasal Pillows for Breeze, Large	ea
SUHCM011MP	Nasal Pillows for Breeze, Medium	ea
SUHCM011SP	Nasal Pillows for Breeze, Small	ea



Citrus II® Cleaning Supplies

Item No.	Description	Pkg
SUH1001	Citrus II® Spray Cleaner, 1.5 oz.	24/cs
SUH1002	Citrus II® Spray Cleaner, 8 oz.	12/cs
SUH1003	Citrus II® Mask Wipes	62 Wipes/ Dispenser, 12/cs
SUH1004	Citrus II® Cleaner Concentrate, 32 oz.	ea
SUH1005	Citrus II® CPAP Wipes	12 bx/cs
SUHCAP1003S	CPAP Cleaning Wipes, Unscented	62 Wipes/ Dispenser, 12/cs



Tubing

Item No.	Description	Pkg
SUH0610	6' Tubing with Flex Connector	10/bx
SUHTUB010	10' Tubing with Flex Connectors	ea
SUHTUB008	8' Tubing with Flex Connectors	ea



Chin Straps

Item No.	Description	Pkg
SUHCS001	Standard Chin Strap	ea
SUHCS003	Adjustable Chin Strap	ea
SUHCS004	Deluxe Chin Strap	ea
SUHCS006	Premium Chin Strap	ea

Pillows

Item No.	Description	Pkg
COACA3042CP	Standard Memory Foam CPAP Pillow	6/cs
COACA3043CP	Zippered Memory Foam CPAP Pillow	4/cs
COACA3044CP	Fiber-Filled CPAP Pillow	6/cs



Non Invasive (NIV) Ventilation Masks

Enhance patient comfort with our full face and nasal masks.

AIR°gel
Maximum comfort, minimum weight. Our patented AIR°gel minimizes pressure points and skin irritation.



Touchless Spacebar™
Unique design eliminates any contact with the patient's forehead.

Size Indicator
available in two sizes.

Quick-Release Headgear
Lightweight, breathable headgear slides off easily with one hand.

Advanced Cushion Technology (A.C.T.)
Our ergonomic cushion technology delivers an effective seal for the most comfortable fit.

Active Headgear Connector
Moves as you move for a reliable seal.

360° Rotating Elbow
For optimal tube positioning and easier mobility

Vented and Non-Vented Versions Available
For CPAP, Bi-level, and ventilator therapy.




CPAP50595

Innova Disposable Nasal Mask

The Innova® Hospital disposable mask feature patented AIR°gel with Advanced Cushion Technology and Active Headgear Connector™. Together, they provide an optimal fit, maximum comfort, and freedom to move without breaking the seal. Move Freely with Our Disposable Hospital Nasal Mask

- » Vented and non-vented versions are available for easy hook up to CPAP, Bi-level or ventilator therapy
- » Easy-to-use with a quick release headgear connector
- » Touchless Spacebar™ eliminates any contact with the forehead, yet retains stability.

Item No.	Description	Size	Pkg
CPAP50180	Disposable Vented Nasal Mask with Headgear, Adult	S/M	1/ea
CPAP50181	Disposable Vented Nasal Mask with Headgear, Adult	M/L	1/ea
CPAP50594	Disposable Non-vented Nasal Mask with Headgear, Adult	S/M	1/ea
CPAP50595	Disposable Non-vented Nasal Mask with Headgear, Adult	M/L	1/ea



Veraseal 2 Disposable Mask

Bring AIRgel technology to the acute care environment for a giant step forward in patient comfort. Three color-coded elbow designs make choosing the right mask for CPAP/Bi-level, single or dual limb circuit NIV systems simple. Fits quickly and easily while enhancing patient comfort.

- » Designed with minimal intentional leaks and elbows that work with CPAP/Bi-level and single or dual limb circuits
- » Flexible spacebar is designed to provide superior stability without touching the skin
- » Active headgear connector moves with the patient to help maintain an effective seal

CPAP/BIPAP (Vented)

Item No.	Description	Size	Pkg
CPAP50827BX	Non Invasive Disposable Hospital Mask with Headgear, Vented	Small	10/cs
CPAP50828BX	Non Invasive Disposable Hospital Mask with Headgear, Vented	Medium	10/cs
CPAP50829BX	Non Invasive Disposable Hospital Mask with Headgear, Vented	Large	10/cs
CPAP50830BX	Non Invasive Disposable Hospital Mask with Headgear, Vented	X-Large	10/cs

BIPAP (Non-Vented)

Item No.	Description	Size	Pkg
CPAP50933BX	Non Invasive Disposable Hospital Mask with Headgear, Non-Vented	Small	10/cs
CPAP50934BX	Non Invasive Disposable Hospital Mask with Headgear, Non-Vented	Medium	10/cs
CPAP50935BX	Non Invasive Disposable Hospital Mask with Headgear, Non-Vented	Large	10/cs
CPAP50936BX	Non Invasive Disposable Hospital Mask with Headgear, Non-Vented	X-Large	10/cs

BIPAP (Non-Vented)

Item No.	Description	Size	Pkg
CPAP50937BX	Non Invasive Disposable Hospital Mask with Headgear, Non-Vented, w/ AAV	Small	10/cs
CPAP50938BX	Non Invasive Disposable Hospital Mask with Headgear, Non-Vented, w/ AAV	Medium	10/cs
CPAP50939BX	Non Invasive Disposable Hospital Mask with Headgear, Non-Vented, w/ AAV	Large	10/cs
CPAP50940BX	Non Invasive Disposable Hospital Mask with Headgear, Non-Vented, w/ AAV	X-Large	10/cs

FILTRATION & HUMIDIFICATION

Heat and Moisture Exchangers	43
Filter HMEs	44
Filters	45
Pulmonary Function Filters	45

Heat and Moisture Exchangers

- » Lightweight
- » Choose from a wide range of tidal volumes, from 90 ml to 1,500 ml



	DYNJAAHME10	DYNJAAHME12	DYNJAAHME6	DYNJAPHME14
Orientation	Straight	Straight	Straight	Angled
Vt (ml)	120-1500	120-1500	150-1500	120-1500
Vd (ml)	34	34	55	42
Weight (g)	20	20	27	30
CO2 Port	YES	YES	NO	YES
Resistance @ 15LPM (cm H ₂ O)	0.27	0.27	-	0.23
Resistance @ 30LPM (cm H ₂ O)	0.81	0.81	0.3	-
Resistance @ 60LPM (cm H ₂ O)	2.6	2.6	0.69	-
Moisture Output (Mg H ₂ O/L)	31	31	31	31
BFE(%)	Foam	Foam	Foam	Foam
Pkg	50/cs	50/cs	50/cs	50/cs

	DYNJAPHME23	DYNJAPHME20	DYNJAPHME18	DYNJAAHME30
Orientation	Angle	Straight	Straight	Barrel
Vt (ml)	>45	90-1500	>10	>25
Vd (ml)	13	25	3.2	11
Weight (g)	9	18	2	5
CO2 Port	YES	YES	NO	N/A
Resistance @ 15LPM (cm H ₂ O)	0.29	0.12	4.77	-
Resistance @ 30LPM (cm H ₂ O)	-	-	-	0.51
Resistance @ 60LPM (cm H ₂ O)	-	-	-	1.16
Moisture Output (Mg H ₂ O/L)	31	35	26	26
BFE(%)	Foam	Foam	Foam	Paper
Pkg	50/cs	50/cs	50/cs	50/cs

Filter HMEs

Medline's wide range of choices ensures compatibility with your equipment.

- » Choose from both angled and straight designs
- » Electrostatic and HEPA configurations available
- » Translucent housing makes it easy to see possible obstructions
- » HME media foam composition



	DYNJAAHME9	DYNJAAHME11	DYNJAAHME5	DYNJAAHME1	DYNJAAHME13
Orientation	Straight	Straight	Straight	Straight	Angled
Vt (ml)	120-1500	120-1500	150-1500	150-1500	120-1500
Vd (ml)	33	33	55	55	41
Weight (g)	20	20	27	27	32
CO2 Port	YES	YES	NO	NO	YES
Resistance @ 15LPM (cm H ₂ O)	-	-	-	-	1.2
Resistance @ 30LPM (cm H ₂ O)	1.2	1.2	0.98	0.98	-
Resistance @ 60LPM (cm H ₂ O)	2.75	2.75	2.28	2.28	-
Moisture Output (Mg H ₂ O/L)	30	30	31	31	29
Filter Media	Electro-Static	Electro-Static	Electro-Static	Electro-Static	Electro-Static
BFE(%)	99.999	99.999	99.999	99.999	99.999
VFE(%)	99.999	99.999	99.999	99.999	99.999
Pkg	50/cs	50/cs	50/cs	50/cs	50/cs

	DYNJAPHME21	DYNJAPHME17	DYNJAPHME19	DYNJAAHMEHEPA
Orientation	Angled	Straight	Angled	Straight
Vt (ml)	>45	90-1500	90-1500	200-1500
Vd (ml)	13	26	32	55
Weight (g)	9	19	18	40
CO2 Port	YES	YES	YES	YES
Resistance @ 15LPM (cm H ₂ O)	2.77	1.2	1.2	-
Resistance @ 30LPM (cm H ₂ O)	-	-	-	1.4
Resistance @ 60LPM (cm H ₂ O)	-	-	-	3.2
Moisture Output (Mg H ₂ O/L)	31	35	36	31
Filter Media	Electro-Static	Electro-Static	Electro-Static	HEPA
BFE(%)	99.999	99.999	99.999	99.9999
VFE(%)	99.999	99.999	99.999	99.9999
Pkg	50/cs	50/cs	50/cs	50/cs

Filters

Exceptional bacterial and viral filtration efficiency – up to 99.9999 percent!

- » Translucent housing makes it easy to see vomitus or other possible obstructions in the system

	DYNJAAHEPA	DYNJAABV1
Orientation	Straight	Straight
Vt (ml)	200-1500	90-1500
Vd (ml)	52	30
Weight (g)	40	23
CO2 Port	YES	NO
Resistance @ 15LPM (cm H ₂ O)	-	-
Resistance @ 30LPM (cm H ₂ O)	1.48	0.77
Resistance @ 60LPM (cm H ₂ O)	3.33	1.78
Filter	HEPA	Electro-Static
BFE(%)	99.9999	99.999
VFE(%)	99.9999	99.999
Pkg	50/cs	50/cs



DYNJAABV1



DYNJAAHEPA

Pulmonary Function Filters

Medline offers a variety of pulmonary function filters that are compatible with all spirometry machines. Please contact your Medline rep or dial 1-800-Medline(633-5463), for ordering information.

- » Independent testing shows >99.99% microbiological filtration effectiveness
- » Low resistance performance ensures efficacy of results
- » Low functional volume
- » Flexible, comfortable, disposable mouth piece helps with patient compliance and improves effectiveness
- » Comfortable single-patient nose clip prevents patient-to-patient touch contamination



RESUSCITATION

Manual Resuscitators	47
Manual Resuscitator Ordering Guide	48
Hyperinflation Systems	49



Manual Resuscitators

Innovations provide clinicians with improved ease of use. Available in infant, pediatric and adult sizes.

A Better Bag

- » Delivers up to 99% FIO₂
- » Maneuverable – The body of the bag swivels 360 degrees for multiple clinician hand positions, helping to reduce clinician fatigue
- » Translucent – The translucent, lightweight bag body and components allow for visualization of vomitus or occlusions in the system
- » Textured – Provides an enhanced, non-slip grip
- » Repositionable – The patient port swivels 360 degrees, allowing the clinician freedom to move around the patient without placing excessive torque or stress on the tracheal tube or mask
- » Manometer allows for a clear visual reading of the centimeters of water pressure being delivered to the patient

An Improved Mask

- » Firmly Sealing– Premium material provides a soft cushion and a firm seal



Fittings and Accessories

Item No.	Description	Pkg
CPRMP19	19 mm Adjustable PEEP Valve	20/cs
CPRMP30	30 mm Adjustable PEEP Valve	20/cs
HCS64184	Mouth-to-Mask Resuscitator with Valve	12/cs
CPRMFILTER	Expiratory Filter	50/cs
CERFLOCAP	CO ₂ Detector with Flow Indicator	24/cs

Product Specifications

	Adult	Pediatric	Infant
Bag Volume	1,900 ml	1,000 ml	300 ml
Stroke Volume	950-960 ml	650-670 ml	140-170 ml
Dead Space with Mask	10 ml	10 ml	10 ml
Inspiratory Resistance	<5 cm H ₂ O	<5 cm H ₂ O	<5 cm H ₂ O
Exhalation Resistance	<5 cm H ₂ O 50 LPM Flow	<5 cm H ₂ O 50 LPM Flow	<5 cm H ₂ O 50 LPM Flow
Attainable Delivery Pressure	>80 cm H ₂ O	>80 cm H ₂ O	>80 cm H ₂ O
Ventilatory Frequency	>60 BPM (No Load)	>60 BPM (No Load)	>90 BPM (No Load)
Body Mass Use Range	>40 kg (88-lbs.)	<40 kg (88-lbs.)	<10 kg (22-lbs.)
Minimal Tidal Volume	600 ml	500 ml	150 ml



Manual Resuscitator Ordering Guide

Item No.	Adult	Pediatric*	Infant*	Reservoir Bag	Reservoir Tube	Filter (F)	PEEP Valve (P)	CO2 Indicator (C)	CO2 FLOCAP	Integrated Manometer	Pkg
CPRM1116	•			•							6/cs
CPRM1116M	•			•						•	6/cs
CPRM1116K	•			•					•		6/cs
CPRM1116P	•			•			•				6/cs
CPRM1116F	•			•		•					6/cs
CPRM1116C	•			•				•			6/cs
CPRM1116FC	•			•		•		•			6/cs
CPRM1116PC	•			•			•	•			6/cs
CPRM1116FPC	•			•		•	•	•			6/cs
CPRM1116FP	•			•		•	•				6/cs
CPRM1126F	•				•	•					6/cs
CPRM1126	•				•						6/cs
CPRM1126M	•				•					•	6/cs
CPRM1126K	•				•				•		6/cs
CPRM1126P	•				•		•				6/cs
CPRM1126C	•				•			•			6/cs
CPRM1126FC	•				•	•		•			6/cs
CPRM1126PC	•				•		•	•			6/cs
CPRM1126FPC	•				•	•	•	•			6/cs
CPRM1126FP	•				•	•	•				6/cs
CPRM2216		•		•							6/cs
CPRM2216M		•		•						•	6/cs
CPRM2216K		•		•					•		6/cs
CPRM2216F		•		•		•					6/cs
CPRM2216P		•		•			•				6/cs
CPRM2216C		•		•				•			6/cs
CPRM2216FC		•		•		•		•			6/cs
CPRM2216PC		•		•			•	•			6/cs
CPRM2216FPC		•		•		•	•	•			6/cs
CPRM2216FP		•		•		•	•				6/cs
CPRM2226		•		•							6/cs
CPRM2226M		•		•						•	6/cs
CPRM2226K		•		•					•		6/cs
CPRM2226F		•			•	•					6/cs
CPRM2226P		•			•		•				6/cs
CPRM2226C		•			•			•			6/cs
CPRM2226FC		•			•	•		•			6/cs
CPRM2226PC		•			•		•	•			6/cs
CPRM2226FPC		•			•	•	•	•			6/cs
CPRM2226FP		•			•	•	•				6/cs
CPRM3312			•	•							12/cs
CPRM3312M			•	•						•	12/cs
CPRM3312K			•	•					•		12/cs
CPRM3312F			•	•		•					12/cs
CPRM3312P			•	•			•				12/cs
CPRM3312C			•	•				•			12/cs
CPRM3312FC			•	•		•		•			12/cs
CPRM3312PC			•	•			•	•			12/cs
CPRM3312FPC			•	•		•	•	•			12/cs
CPRM3312FP			•	•		•	•				12/cs
CPRM3322			•		•						12/cs
CPRM3322M			•		•					•	12/cs
CPRM3322K			•		•				•		12/cs
CPRM3322F			•		•	•					12/cs
CPRM3322P			•		•		•				12/cs
CPRM3322C			•		•			•			12/cs
CPRM3322FC			•		•	•		•			12/cs
CPRM3322PC			•		•		•	•			12/cs
CPRM3322FPC			•		•	•	•	•			12/cs
CPRM3322FP			•		•	•	•				12/cs


*Includes 40 cm pressure relief valve

Hyperinflation Systems

Help prevent or treat lung collapse with a hyperinflation system that performs as well as they do.

- » Bag design uses a “dipped” process to enhance tactile sensitivity and improve control
- » Uniform inflation and the highest level of performance
- » 100% compliant with ISO 5362 standards
- » Variety of sizes to more accurately fit every patient

Swivel Port ...
Optimal mobility in confined spaces.




Pleated Breathing Bag
Four sizes to meet the needs of all your patients.

Pressure-Release Dial
Clinician adjusts dial to achieve and maintain targeted base pressure.

Manometer
Color-coded dial to measure delivered pressure.

Pressure-Release Valve
At 40 or 25 cm/H2O the valve releases pressure to help prevent excessively high inspiratory pressures.



Oxygen Tubing ...
Choose from standard or universal connectors.



Mask ...
Available in all sizes, both inflatable and pre-inflated.

Hyperinflation Systems

7' (2.13 m) of O₂ tubing included*

Item No.	Breathing Bag	O ₂ Connector	Mask	Pkg.
DYNJAAHY005	0.5 L	Standard	N/A	20/cs
DYNJAAHY600	0.5 L w/Tail-End Bleed Valve	Universal	N/A	50/cs
DYNJAAHY30	0.5 L	Standard	1 Neonatal	20/cs
DYNJAAHY31	0.5 L	Standard	2 Infant	20/cs
DYNJAAHY29	0.5 L	Universal	2 Infant (Pre-Inflated)	20/cs
DYNJAAHY01	1 L	Standard	N/A	20/cs
DYNJAAHY601	1 L w/Tail-End Bleed Valve	Universal	N/A	50/cs
DYNJAAHY01M2	1 L	Universal	2 Infant	20/cs
DYNJAAHY02	2 L	Standard	N/A	20/cs
DYNJAAHY02M3	2 L	Universal	3 Child	20/cs
DYNJAAHY03	3 L	Standard	N/A	20/cs
DYNJAAHY03M5	3 L	Universal	5 Adult	20/cs



DYNJAAHY601

Hyperinflation Systems with Manometers

7' (2.13 m) of O₂ tubing included*

Item No.	Breathing Bag	O ₂ Connector	Mask	Pkg.
DYNJAAHY5005	0.5 L	Standard	N/A	50/cs
DYNJAAHY5053P	0.5 L	Universal	N/A	50/cs
DYNJAAHY32	0.5 L	Standard	1 Neonatal	50/cs
DYNJAAHY33	0.5 L	Standard	2 Infant	50/cs
DYNJAAHY35	0.5 L	Universal	2 Infant	50/cs
DYNJAAHY34	0.5 L	Universal	2 Infant (Pre-Inflated)	50/cs
DYNJAAHY501	1 L	Universal	N/A	50/cs
DYNJAAHY501M2	1 L	Universal	2 Infant	50/cs
DYNJAAHY502	2 L	Universal	N/A	50/cs
DYNJAAHY502M3	2 L	Universal	3 Child	50/cs
DYNJAAHY503	3 L	Universal	N/A	50/cs
DYNJAAHY503M5	3 L	Universal	5 Adult	50/cs



DYNJAAHY5005

*Custom options are available

PATIENT MONITORING & ACCESSORIES

Fingertip Pulse Oximeters	51
Handheld Oximeters	52
Accessories	53
Pulse Oximeter Probes	54

Finger Pulse Oximeters

Get quick and precise SpO₂ and pulse readings in seconds with our high-quality fingertip pulse oximeters. Small and lightweight, these clip sensors automatically take readings when a finger is inserted. They fit a wide range of finger sizes, have a long battery life and are solidly constructed for worry-free use. All models include AAA batteries, lanyard and warranty.



pulSTAT™ Finger Pulse Oximeter

Patented multi-directional display offers 6 viewing modes.

- » SpO₂, pulse rate, pulse bar, plethysmograph
- » Low power consumption provides 30 hours of use
- » Bright, digital OLED with adjustable readout sizes
- » Soft padding for patient comfort
- » Hinged battery cover

HCSM70D, 1/ea



Easy-Grip Finger Pulse Oximeter

Hinges adjust to comfortably accommodate larger fingers.

- » SpO₂, pulse rate, pulse bar, plethysmograph
- » Backlit LED display with 4 viewing modes
- » Slip-resistant grip
- » Hinged battery cover

HCSM70, 1/ea



Soft-Touch Finger Pulse Oximeter

Soft padding enhances patient comfort.

- » SpO₂, pulse rate, pulse bar, plethysmograph
- » Backlit LED display with 4 viewing modes

HCSM70C, 1/ea



Basic Finger Pulse Oximeter

Digital display lights up for easy reading in all lighting conditions.

- » SpO₂, pulse rate, pulse bar
- » Bright LCD screen

HCSM70B, 1/ea



High-Impact Finger Pulse Oximeter

Rugged casing withstands accidental drops.

- » SpO₂, pulse rate and pulse bar readings
- » Multi-directional display with full-color OLED screen
- » Durable silicone boot construction resists cracking and chipping
- » Hinged battery cover

HCSFSATSPORT, 1/ea



Pediatric Finger Pulse Oximeter

Child-friendly frog design makes kids smile.

- » SpO₂, pulse rate, pulse bar, plethysmograph
- » Patented multi-directional display has 6 viewing modes
- » Bright, digital OLED with adjustable readout sizes
- » Soft padding for patient comfort
- » Hinged battery cover

HCSM70P, 1/ea

BCIWW1026SYS
(Printer Sold Separately)



BCIWW1030EN

BCIWW1025SYS
(Docking Station Sold Separately)



BCIWW1000EN

Handheld Oximeters

Dependable and easy-to-use pulse oximeters with alarms for spot checks and continuous monitoring.

SPECTRO2® 30 Digital Pulse Oximetry System

- » Visual and audible alarms, including a remote alarm through the nurse call system
- » Low perfusion
- » Motion tolerance
- » 4 AA batteries

SPECTRO2® 10 Pulse Oximetry System

- » Large, bright LED display
- » Simple, one-handed operation
- » Durable construction
- » Comfort clip finger probe
- » 4 AA batteries

Medline MD300

- » SpO₂ and pulse rate
- » Adjustable display: 2 readout sizes and 10 brightness levels
- » Includes 1 adult finger probe
- » Automatic power off
- » Includes protective silicone cover

Medline MD400

- » SpO₂, pulse rate, pulse bar, waveform
- » Adjustable audible and visual alarms
- » Two display modes
- » Comfortable finger clip probe
- » Data transfer to computer via USB cable
- » Low battery indicator



HCSMD400



HCSMD300

Item No.	Description	Pkg
BCIWW1030EN	SPECTRO ₂ 30 Continuous Pulse Oximetry System with 20" (50.8 cm) Comfort Clip Finger Probe	ea
BCIWW1000EN	SPECTRO ₂ 10 Spot Check Pulse Oximetry System with 20" (50.8 cm) Spot-Check Probe	ea
HCSMD300	Handheld Spot Check Oximeter with 30" (70 cm) Adult Finger Probe	ea
HCSMD400	Handheld Continuous Economy Pulse Oximeter with Adult Finger Probe	ea

BCIWW1018B



BCIWW1018R



BCIWW1018G



BCIWW1018Y



Accessories

Item No.	Description	Pkg
BCIWW1025SYS	SPECTRO ₂ Docking Station » Includes AC Power Adapter (BCIWW1095) » For BCIWW1000 and BCIWW1030 Oximeters	ea
BCIWW1090	Rechargeable Battery » For BCIWW1000 and BCIWW1030 Oximeters	ea
BCIWW1018B	SPECTRO ₂ Protective Glove, Blue	ea
BCIWW1018G	SPECTRO ₂ Protective Glove, Green	ea
BCIWW1018Y	SPECTRO ₂ Protective Glove, Yellow	ea
BCIWW1018R	SPECTRO ₂ Protective Glove, Red	ea
BCIWW1089	SPECTRO ₂ USB Cable	ea
BCIWW1095	SPECTRO ₂ Mounting Bracket with Clamp » For Models BCIWW1000and BCIWW1030 Oximeters	ea
BCIWW1098C	SPECTRO ₂ Mounting Bracket with Clamp	ea
BCIWW1027	SPECTRO ₂ Printer Paper » For Model BCIWW1026SYS	ea
BCIWW1026SYS	SPECTRO ₂ Detachable Printer » Includes 5-pk Printer Paper » Requires Docking Station (BCIWW1025SYS)	ea



BCIWW1026SYS



BCIWW1095

7 - Pin Nellcor



9 - Pin BCI



Pulse Oximeter Probes

Medline has probes for a wide variety of pulse oximeters, including handheld and tabletop models, such as the popular 7-pin Nellcor® and 9-pin BCI®.

- » Probe choices for unique requirements, materials and applications
- » Disposable probes in neonatal through adult sizes
- » Reusable probes in pediatric and adult sizes
- » Micro-Foam for easy sensor application and repositioning, Medplast for maximum adhesive for long-term use in one location, or Transpore for transparency with repositioning

Finger Clip
PCL0603120



Adult

Item No.	Category	Type	Compatible With	Pkg
PCL0603104	Reusable	Finger Clip	BCI	ea
PCL0202301	Disposable	Micro-Foam Style	BCI	24/cs
PCL0603120	Reusable	Finger Clip	Nellcor	ea
PCL1003120	Reusable	Soft Boot	Nellcor	ea
PCL0342121	Disposable	Medplast Style	Nellcor	24/cs
BCI3078	Reusable	Ear	BCI	ea



Medplast
PCL0342121

Pediatric

Item No.	Description	Type	Compatible With	Pkg
PCL0623104	Reusable	Finger Clip	BCI	ea
PCL1023104	Reusable	Soft Boot	BCI	ea
PCL0202302	Disposable	Micro-Foam Style	BCI	24/cs
PCL0322302	Disposable	Transpore Style	BCI	24/cs
PCL0623120	Reusable	Finger Clip	Nellcor	ea
PCL1023120	Reusable	Soft Boot	Nellcor	ea
PCL0202122	Disposable	Micro-Foam Style	Nellcor	24/cs
PCL0322122	Disposable	Transpore® Style	Nellcor	24/cs



Micro-Foam
PCL0202301

Infant

Item No.	Description	Type	Compatible With	Pkg
PCL0202303	Disposable	Micro-Foam Style	BCI	24/cs
PCL0202123	Disposable	Micro-Foam Style	Nellcor	24/cs



Transpore
PCL0322101

Neonatal

Item No.	Description	Type	Compatible With	Pkg
PCL0202304	Disposable	Micro-Foam Style	BCI	24/cs
PCL0342124	Disposable	Medplast Style	Nellcor	24/cs
PCL0202124	Disposable	Micro-Foam Style	Nellcor	24/cs
PCL0703120	Disposable	Multi Site Sensor	Nellcor	ea



Soft Boot
PCL1003120

Universal

Item No.	Description	Type	Compatible With	Pkg
PCL0703104	Reusable	Multi Site Sensor	BCI	ea

OXYGEN REGULATORS & CARTS

Conservers	56
Regulators	57
Cylinder Transport Bags	58
Cylinder Storage Racks and Carts	58
Cylinder Transport Carts	59


EconO₂mizer II

Single-Lumen Pneumatic Conservers

Less noise for better compliance

Among the quietest pneumatic conserving devices available.

- » Very quiet with no hissing during pulse delivery
- » Aids in compliance with ambulatory patients
- » 6 conserving flow settings (1, 2, 2.5, 3, 4, & 5)
- » 3 continuous settings (2, 3, & 4)
- » Features all-brass construction in high-pressure zones
- » Weight: 13-oz.
- » Size: 4.8" x 2.4" x 1.5" (12.2 cm x 6.1 cm x 3.8 cm)
- » 5.7:1 conserving ratio




HCS2CON1L

Item No.	Description	Pkg
HCS2CON1L	Single-Lumen Pneumatic Conservers-870 connection	ea
HCSTHANDLE	Replacement T Handle, Black	ea
HCS166970	Replacement Tubing Barb Port Outlet to both HCSO2CON1 and HCSO2CONS	ea

EconO₂mizer Dual-Lumen Conserving Device

A breakthrough in conserving-device technology!

- » All-brass diaphragm and 4:1 conserving ratio for fewer service calls
- » Light and compact
- » 0-6 LPM
- » Dual mode for conserving and continuous flows
- » 11 settings for each mode (0,1, 1.5, 2, 2.5, 3, 3.5, 4, 4.5, 5, 5.5 & 6) (no additional regulator required)
- » Solid brass core



Dual Mode Switch

HCSO2CONS

Item No.	Description	Pkg
HCSO2CONS	Conservers, 0-6 Settings, 870 Connection	20/cs
HCSTHANDLE	Replacement T Handle, Black	20/cs



CPRB87251B2C

HCS8715M

HCS5415M

HCS8715MA

Medline Regulators

For all medical and emergency applications.

- » Fitted with brass sleeves around key parts
- » Laser-etched knob and body

Economy

Item No.	Connection	Flow Range	Flow Increments	Pkg
HCS8708M	870 CGA Connection	0-8 LPM	0.12, 0.5, 1, 1.5, 2, 2.5, 3, 4, 5, 6, 8	ea
HCS8715M	870 CGA Connection	0-15 LPM	0.12, 0.5, 1, 1.5, 2, 3, 4, 6, 8, 10, 15	ea
HCS8725M	870 CGA Connection	0-25 LPM	0.5, 1, 2, 3, 4, 5, 6, 8, 10, 15, 25	ea
HCS5408M	540 CGA Connection	0-8 LPM	0.12, 0.5, 1, 1.5, 2, 2.5, 3, 4, 5, 6, 8	ea
HCS5415M	540 CGA Connection	0-15 LPM	0.12, 0.5, 1, 1.5, 2, 3, 4, 6, 8, 10, 15	ea

Super Mini

HCS8715MA	870 connection	0-15 LPM	0.12, 0.5, 1, 1.5, 2, 3, 4, 6, 8, 10, 15	ea
-----------	----------------	----------	--	----

Accessories

Item No.	Description	Pkg
HCS166840	Brass Seal Washer for O ₂ regulators	50/bx
HCS0700	Cylinder Wrench	50/bx

EMS Regulators

Easily transport ventilator patients.

- » Maximum inlet pressure: 3,000 PSI
- » Chrome plated brass body
- » Solid brass core

Item No.	50 PSI Power Take-Off	Hose Barb	Liters Per Minute	Cylinder Style	Pin Index	Pkg
CPRB34501C	1			CGA346	Air	ea
CPRB87251B	2	●	0-5-1-2-3-4-5-6-8-10-15-26	CGA870	Oxygen	ea
CPRB87251B2C	2	●	0-5-1-2-3-4-5-6-8-10-15-26	CGA870	Oxygen	ea
CPRB54251B2C	2	●		CGA540	Oxygen	ea
CPRB54402C	2			CGA540	Oxygen	ea
CPRB54501C	1			CGA540	Oxygen	ea



HCSRNR2001W

Handle adjusts
from 36" (91.4 cm)
to 42" (106.7 cm)



D-Tank Height
HCS53006

Rack'n Roll Modular Cylinder Cart

For easy moving and storing of aluminum medical cylinders.

- » Completely modular
- » Heavy-duty rubber tires (removable for stationary cart)
- » Telescoping with adjustable holder height, enabling each tank tray to accommodate different sizes
- » Comfortable handgrip
- » Six-cylinder carts can link to make larger systems (12, 18, 24 cylinders)

Item No.	Description	Pkg
HCSRNR2001W	Rack'N Roll Modular Cart System	ea

D & E Cylinder Rolling Cart

- » Cart accommodates D and E gas cylinders
- » Helps patients to easily transport their own oxygen

Item No.	Description	Pkg
HCS53006	D-E Cylinder Rolling Carts	6

Oxygen Cylinder Bags

Medline is your one source for all your related components.

Item No.	Description	Pkg
HCSM6BAG6	Carry Bags, Shoulder Style, for M6 Cylinders	6/cs
HCSCBAG6	Carry Bags, Shoulder Style, for C Cylinders	6/cs
HCSDBAG6	Carry Bags, Shoulder Style, for D Cylinders	6/cs
HCSM6BKPK	Backpack, Cylinder, F/M6	ea



HCSM6BAG6



HCSCBAG6



HCSDBAG6



HCSM6BKPK

Cylinder Transport

- » Streamlined and compact designs
- » Store cylinders in minimal space

Item No.	Description	Tank Size	Pkg
RRI1500102	Cart for 2 Cylinders with Adjustable Handle	D, E	ea
RRI1500130	Cart for 4 Cylinders	D, E	ea
RRI1500172	Cart for 12 Cylinders, QuietRide® Casters	D, E	ea
RRI1500175	Cart for 12 Cylinders, Large Wheels	D, E	ea
RRI1500240	Stand for 1 Cylinder	D, E	ea
RRI1500261	Rack for 6 Cylinders	D, E	ea
RRI1500265	Composite for 8 Cylinders	D, E	ea
RRI1500275	Composite for 12 Cylinders	M6	ea
RRI1500271	Rack for 12 Cylinders	D, E	ea









Medline Industries, Inc. One Medline Place, Mundelein, IL 60060

Medline United States
 1-800-MEDLINE (633-5463)
 medline.com | info@medline.com

Medline Canada
 1-800-396-6996
 medline.ca | canada@medline.com

Medline México
 01-800-831-0898
 medlinemexico.com | mexico@medline.com

FOLLOW US    

REFERENCES 1. Black JM, Cuddigan JE, Walko MA, Diddier LA, Lander MJ, Kelpo MR. Medical device related pressure ulcers in hospitalized patients. International Wound Journal. 2010; 7(5):358-365. 2. American Lung Association website. Measuring Your Peak Flow Rate. Available at: <http://www.lung.org/lung-health-and-diseases/lung-disease-lookup/asthma/living-with-asthma/managing-asthma/measuring-your-peak-flow-rate.html>. Accessed September 27, 2016.

Some products may not be available for sale in Mexico or Canada. We reserve the right to correct any errors that may occur within this brochure. © 2016 Medline Industries, Inc. Aeromist Compact, Aeromist Buddies, and Aeromist Colors are trademarks, and Aloetouch, DoggyNeb, and Medline are registered trademarks of Medline Industries, Inc. Aeroflex and Phoneris are registered trademarks of Lazarus Medical, LLC. Airgel, iQ Blue, MiniMe and Mojo are registered trademarks of Sleepnet Corporation. Oxequip is a registered trademark of Allied Healthcare Products, Inc. Citrus II is a registered trademark of Beaumont Products, Inc. Eclipse is a registered trademark of SeQual Technologies Inc. Gentec is a registered trademark of Genstar Technologies Co., Inc. Nellcor is a registered trademark of Covidien LP Corporation. Ohio is a registered trademark of Ohio Medical Corporation. Ohmeda is a registered trademark of GE Healthcare. QuietRide is a registered trademark of The Mail Store, Inc. Trachvox is a registered trademark of Engineered Medical Systems, Inc. Transpore is a registered trademark of 3M Corporation. MKT1655429 / LIT085 / 2.5M / JP / 15