

# GV200-4 | GV210-4 WALKERS



Item #	Product Description	Product Height	UOM
GV200-4	Adult	32" - 39"	4/cs
GV210-4	Adult with 5" Wheels	32" - 39"	4/cs

## FEATURES

- Easy push-button mechanisms may be operated by fingers, palms or side of hand
- Each side operates independently to allow easy movement through narrow spaces and greater stability while standing
- Sturdy 1" diameter, aluminum construction ensures maximum strength while remaining lightweight
- Vinyl contoured hand grip
- GV210-4 comes with 5" wheels

## Specifications

Width GV200-4	17" inside hand grip   19" inside back legs
Width GV210-4	17" inside hand grip   20" inside back legs
Overall Width GV200-4   GV210-4	22"   24"
Depth (Opened at Base) GV200-4   GV210-4	15.5"   17.5"
Weight GV200-4   GV210-4	5.5 lbs.   6 lbs.
Carton Shipping Weight GV200-4   GV210-4	26 lbs.   30 lbs.
Weight Limit	350 lbs.
Warranty	Limited Lifetime

# GV220-2WW WALKER



Item #	Product Description	Product Height	UOM
GV220-2WW	Adult	32" - 39"	2/cs

## FEATURES

- For Bariatric use
- Steel legs and side braces
- Precision design provides additional strength while adding minimal weight
- Wider and deeper frame design accommodates individuals up to 450 lbs
- Easy push-button mechanisms may be operated by fingers, palms or side of hand
- Each side operates independently to allow easy movement through narrow spaces and greater stability
- Comes with heavy duty wheels

## Specifications

Width	20.5" inside hand grip   23" inside back legs
Overall Width	27"
Depth (Opened at Base)	23"
Weight without Wheels	5.7 lbs.
Weight with Wheels	9.5 lbs.
Carton Shipping Weight	24 lbs.
Weight Limit	450 lbs.
Warranty	Limited Lifetime

# BARIATRIC WALKER TRAY

## Bariatric Walker Tray with Cup Holders



Item #	Product Description	UOM
10125HD	Bariatric Walker Tray with Cup Holders, Gray	1/ea

### FEATURES

- Fits most 2-button bariatric walkers
- Allows personal items to be carried from room to room
- Contains two cup holders
- Easily folds down when not in use
- Easy-to-install
- Size: 20.0 x 11.6 x 2.0

### Specifications

Warranty

Limited Lifetime