



drive

DeVilbiss
HEALTHCARE



2018 Bath Safety Product Catalog

Table of Contents

| | |
|---------------------|---|
| Commodes | 1 |
| Shower Chairs | 2 |
| Transfer Benches | 3 |
| Raised Toilet Seats | 5 |
| Bath Lift | 6 |
| Grab Bars | 7 |
| Shower Spray | 8 |





Folding Steel Commode with Deep Seat

Heavy Duty Drop Arm Commode

| Item # | Product Description | UOM |
|-----------|--------------------------------------|------|
| GV33048-1 | Folding Steel Commode with Deep Seat | 1/ea |

| Item # | Product Description | UOM |
|-----------|-------------------------------|------|
| GV33035-1 | Commode, Heavy Duty, Drop Arm | 1/ea |

- Tool-free Assembly
- Easily opens and folds
- Folds flat for convenient storage and transportation
- Durable plastic snap-on seat and lid
- Assembles without pushing buttons
- Grey powder-coated steel-welded construction increases strength and durability
- Plastic arms for added comfort
- Comes complete with 8 qt commode bucket with carry handle, cover and splash shield

- Tool-free Assembly
- Easy-to-release arm mechanism allows for safe lateral patient transfers to and from commode
- Durable, heavy-duty, gray powder-coated steel tubing is sturdy and easy-to-maintain
- Large, durable snap-on seat
- Padded armrests provide comfort and stability
- Pail easily slides in from the front of the commode
- Comes complete with 12 qt commode bucket and cover that removes from front

| specifications | |
|--------------------|------------------------|
| Seat Height | 15.5" - 21.75" |
| Width Between Arms | 18" |
| Overall Width | 22.25" |
| Seat Dimensions | 13.25" (W) x 16.5" (D) |
| Weight Capacity | 350 lbs |

| specifications | |
|--------------------|------------------------|
| Seat Height | 17.5" - 22" |
| Width Between Arms | 25.75" |
| Overall Width | 26.25" |
| Seat Dimensions | 23.25" (W) x 18.5" (D) |
| Weight Capacity | 900 lbs |



Shower Chair with Back

Tool Free

Bariatric Shower Chair with Back

Heavy Duty, Tool Free

| Item # | Product Description | UOM |
|------------|-----------------------------------|------|
| GV3302KD-1 | Shower Chair with Back, Tool Free | 1/ea |

- Easy, tool-free assembly of back, seat and legs
- Aluminum frame is lightweight, durable and corrosion proof
- Depress push pins for height adjustment, the audible "click" will be heard when secure
- Legs are height adjustable in 1" increments, and angled with suction-style tips to provide additional stability
- Easy to clean

specifications

| | |
|--------------------|-----------------|
| Seat Height | 15"-19.25" |
| Width Between Arms | NA |
| Overall Width | 20.25" |
| Seat Dimensions | 11.25" x 20.25" |
| Weight Capacity | 400 lbs |

| Item # | Product Description | UOM |
|-------------|--|------|
| GV33021KD-1 | Shower Chair with Back Heavy Duty, Tool Free | 1/ea |

- Easy tool-free assembly of back and seat
- Frame has width/depth cross braces attached with aircraft-type rivets
- Durable, blow-molded plastic bench
- Aluminum frame is lightweight, durable and corrosion proof
- Drainage holes reduce slipping
- Euro-style clip for leg height adjustment
- Height adjusts in 1" increments
- Easy to clean

specifications

| | |
|--------------------|-----------------|
| Seat Height | 14.5"-19.25" |
| Width Between Arms | NA |
| Overall Width | 20.25" |
| Seat Dimensions | 20.25" x 12.25" |
| Weight Capacity | 500 lbs |



Shower Chair with Padded Arms with Back

Transfer Bench Tool-free Assembly Back, Legs and Arms

| Item # | Product Description | UOM |
|-------------|--|------|
| GV33445KD-1 | Shower Chair with Back, Padded Arms, Tool Free | 1/ea |

- Easy tool-free assembly of back, legs and arms
- Removable padded armrests provide comfort and added security
- Blow-molded bench and back provide comfort and strength
- Drainage holes in bench reduce slipping
- Aluminum frame is lightweight, durable and corrosion proof
- Leg height adjusts in 1" increments
- Easy-to-clean

specifications

| | |
|--------------------|----------------|
| Seat Height | 16.25"- 21.25" |
| Width Between Arms | 17.25" |
| Overall Width | 20.75" |
| Seat Dimensions | 16.25" x 13.5" |
| Weight Capacity | 300 lbs |

| Item # | Product Description | UOM |
|-------------|--|------|
| GV33011KD-1 | Transfer Bench with Back, Arm, Tool Free | 1/ea |

- Easy tool-free assembly of back, legs and arm
- Durable blow-molded plastic bench and backrest
- Height adjusts in 1/2" increments with unique "Dual Column" extension legs
- Reversible to accommodate any bathroom
- Easy to clean

specifications

| | |
|--------------------|-------------|
| Seat Height | 17.5"-21.5" |
| Width Between Arms | NA |
| Overall Width | 35" |
| Seat Dimensions | 26.5" x 16" |
| Weight Capacity | 400 lbs |



Padded Transfer Bench

Tool-free Assembly Back, Legs and Arms

Bariatric Transfer Bench

with Cross-Frame Brace

| Item # | Product Description | UOM |
|------------|--|------|
| GV3305KD-1 | Transfer Bench, Padded with Back, Arm, Tool Free | 1/ea |

- Tool-free assembly of back, legs and arm
- Comfortable cushioned seat and backrest
- Height adjusts in 1/2" increments with unique "Dual Column" extension legs
- Easy-to-clean

specifications



| | |
|--------------------|-------------|
| Seat Height | 18"-22.375" |
| Width Between Arms | N A |
| Overall Width | 30" |
| Seat Dimensions | 23.5" x 16" |
| Weight Capacity | 400 lbs |

| Item # | Product Description | UOM |
|-----------|--|------|
| GV33009-1 | Transfer Bench, Heavy Duty with Back, Arm, Tool Free | 1/ea |

- Durable, blow-molded plastic bench and backrest
- Aluminum frame is lightweight, sturdy and corrosion resistant
- Extra middle legs for added strength
- Seat height adjusts in 1" increments
- No tools required for installation
- Easy-to-clean

specifications

| | |
|--------------------|-------------|
| Seat Height | 16.5"-20.5" |
| Width Between Arms | NA |
| Overall Width | 30.5" |
| Seat Dimensions | 26.5" x 16" |
| Weight Capacity | 600 lbs |



| | |
|--|---|
| <h3>Elevated Toilet Seat</h3> <p>with Arms or Removable Arms</p> | <h3>Elevated Toilet Seat</h3> <p>without Arms</p> |
|--|---|

| Item # | Product Description | UOM |
|-----------|--|------|
| GV33027-1 | Elevated Toilet Seat with Arms or Removable Arms | 1/ea |

- Designed for individuals who have difficulty sitting down or getting up from the toilet
- Adds 5" to height of toilet seat
- Tool-free arms can be removed or added
- Arm sleeves and receivers are made of metal
- Heavy-duty molded plastic construction
- Locking device with larger, heavy-duty "worm screw" and locking plate for a safe, secure fit
- Lightweight and portable
- Easily attaches to toilet bowl, tool-free
- Easy-to-clean
- Fits Standard & most Elongated toilets

| specifications | |
|--------------------|----------------|
| Seat Height | 5.0" |
| Width Between Arms | 17.25" |
| Overall Width | 20.5" |
| Seat Dimensions | 20.5" X 16.25" |
| Weight Capacity | 300 lbs |

| Item # | Product Description | UOM |
|-----------|-----------------------------------|------|
| GV33026-1 | Elevated Toilet Seat without Arms | 1/ea |

- Designed for individuals who have difficulty sitting down or standing up from the toilet
- Heavy-duty molded plastic construction provides strength and durability
- Locking device with larger, heavy-duty "worm screw" and locking plate for a safe, secure fit
- Lightweight and portable
- Fits Standard & most Elongated toilets
- No tools required for installation
- Easy-to-clean

| specifications | |
|--------------------|----------------|
| Seat Height | 4.75" |
| Width Between Arms | NA |
| Overall Width | 20.5" |
| Seat Dimensions | 20.5" x 16.25" |
| Weight Capacity | 300 lbs |



Auto Bath Lift

| Item # | Product Description | UOM |
|-----------|---------------------|------|
| 477200252 | Auto Bath Lift | 1/ea |

- The lightest bathlift on the market at only 20.5 lbs
- The Bellavita lowers to 2.3" (Figure A)
- Seat height reaches 18.8" for deeper tubs
- The backrest reclines to 50° at its lowest position
- Completely tool-free set up makes it easy-to-install and operate When not in use, back easily folds down for space-saving storage (Figure B)
- Watertight hand control floats
- The easy-to-use quick charge lithium-ion battery is accommodated in the hand control (Figure C)
- Two-piece construction makes the lift easy to assemble/ disassemble for transport

- Padded seat and backrest provide added comfort
- Washable hygienic covers are sold separately
- Hygienic cutaway and slimline profile for easy handling
- 4 quick-release suction cups for anti-slip bath base
- Headrest available as an option
- Simplify transfers, positioning and rotating when using your bath lift with a Vitaturn Swivel and Transfer Aid
- Suction cups adhere to smooth surfaces only, and not textured surfaces

| specifications | |
|---------------------|--|
| Overall Dimensions | 32.6" (L) (with backrest in a reclined position) x 27.5" (W) (seat incl. side flaps) |
| Seat Dimensions | 13.7" (W) x 16.5" (D) x 22.4" (L) |
| Sitting Height | 2.3" - 18.8" |
| Backrest Dimensions | 25.1" (H) x 12.5" - 14.1" (W) |
| Backrest Reclines | 50° |

| | |
|-----------------------------------|---------------------|
| Baseplate Dimensions | 12.2" (W) x 22" (L) |
| Load Capacity Max. | 300 lbs. |
| Total Weight (excl. hand control) | 20.5 lbs. |
| Shipping Weight | 25.3 lbs |
| Warranty | Limited 3 Year |



White Powder-Coated Grab Bars

Chrome Knurled Grab Bar

| Item # | Product Description | UOM |
|----------|---------------------|------|
| 12016-3 | 16" Grab Bar | 3/bx |
| 12024-3 | 24" Grab Bar | 3/bx |
| RTL12012 | 12" Grab Bar | 3/bx |

- Durable, white powder-coated steel is attractive and easy to clean
- Mounts vertically or horizontally
- 1 ½" distance from wall
- Excess water should be wiped off grab bar after each use to prevent rusting
- Attractive retail packaging

| Item # | Product Description | UOM |
|----------|---------------------|------|
| 12112-3 | 12" Grab Bar | 3/bx |
| 12116-3 | 16" Grab Bar | 3/bx |
| RTL12118 | 18" Grab Bar | 3/bx |

- Knurled grab bar diameter provides a no-slip grip
- Mounts vertically or horizontally
- 1 ½" distance from wall
- Attractive retail packaging

specifications

Warranty Limited Lifetime

specifications

Warranty Limited Lifetime



Deluxe Handheld Shower Spray with Diverter Valve

| Item # | Product Description | UOM |
|--------|---------------------|------|
| 12037 | White Shower Spray | 1/bx |

- On/Off switch built into handle for easy access
- Extra-long 80" white reinforced nylon hose provides extra convenience
- Comes with diverter valve (Figure A) that allows either the regular shower head to operate or be used as a handheld shower spray (Figure B)
- Useful when shower chair, bench or stool is being used

specifications

Warranty

Limited Lifetime

Drive DeVilbiss Healthcare
99 Seaview Boulevard
Port Washington, NY 11050
t: 877.224.0946 | f: 516.998.4601
www.drivemedical.com

