

Nitro Aluminum Rollator, Tall Height, 10" Casters

HCPCS: E0143 E0156-1,E0143+E0156



FIGURE A



FIGURE B

- Tall height ideal for users 6'0" and above
- Attractive, Euro-style design
- Brake cable inside frame for added safety
- Handle height easily adjusts with unique push button (Figure A)
- Back support height easily adjusts with tool-free thumb screw (Figure B)
- Caster fork design enhances turning radius
- Large 10" front casters allow optimal steering and rolling comfort
- Easily folds with one hand to ultra-compact size for storage (Figure C)
- Plastic clip holds the Nitro closed when frame is folded (Figure D)
- Cross-brace design allows for side-to-side folding and added stability
- Seat is durable and comfortable (Figure E)
- Handsome and removable zippered storage bag with unique attachment to keep bag securely in place when Nitro is open or folded (Figure F)
- Lightweight, aluminum frame

PRODUCTS

	Description
RTL10266-T	Nitro Aluminum Rollator, Tall Height, 10" Casters, EA
RTL10266BK-T	Nitro Aluminum Rollator, Tall Height, 10" Casters, EA
RTL10266WT-T	Nitro Aluminum Rollator, Tall Height, 10" Casters, EA

SPECIFICATIONS

	Description
Carton Shipping Weight	30 lbs.
Handle Height	36" - 41"
Seat Dimensions	18" (W) x 23.5" (H) x 8.5" (D)
Seat Height	23.6"
Unit Dimensions	23" (W) x 27.75" (L)
Warranty	Limited Lifetime 5 Year Limited
Weight	17.5 lbs.
Weight Capacity	300 lbs

ACCESSORIES

	Description
10266-CH	Nitro Cup Holder, 1/ea, EA

PARTS

	Description
1026627-TALL	CrossBraceAssemblyforNitroTall, 1/ea, EA
1026634	CrossFrameEndCapsToteBagW/XFrameAtt,1/ea, EA
1026601-RNR	HandBrakeAssembly,Right,10266/1,ea, EA

PARTS	Description
1026601-GR	HandBrakeAssmby,Right,10266BK/10266BK-T, EA
1026633-R	ToteBag,CrossFrameAttachment,10266,1/ea, EA
1026633-BK	ToteBag,CrossFrameAttachment,10266BK,1ea, EA
1026608-TALL	SeatUpholW/AttachingHardwareForNitroTall, EA
1026609-FSW-TALL	SeatLiftingStrapW/AttachHrdwre,10266WT-, EA
1026609-FSR-TALL	SeatLiftingStrapW/AttachHrdwre,10266-T, EA
1026609-FSG-TALL	SeatLiftingStrapW/AttachHrdwre,10266BK-, EA
1026601-WR	HandBrakeAssmby,Right,10266WT, 1/ea, EA
1026632	Brake Cable Only, Nitro, 1/ea, EA
1026601-GR-T	HandBrakeAssmby,Right,NitroTall,Black,1, EA
1026601-G-T	HandBrakeAssmby,Left,NitroTall,Black,1e, EA
1026601-WR-T	HandBrakeAssmby,Right,NitroTall,White,1, EA
1026601-RNR-T	HandBrakeAssmby,Right,NitroTall,Red,1e, EA
1026601-RN-T	HandBrakeAssembly,Left,NitroTall,Red,1e, EA
1026632-T	Brake Cable & Sleeve, Nitro Tall, 1/ea, EA
1026633-WT	ToteBag,CrossFrameAttachment,10266WT,1ea, EA
1026646	CrossFrameBolt,Nitro/NitroHemiHD/NitroH, EA