



F-22 ROLLATOR



Item #	Color	UOM
GV974RD-T	Red 	1/cs
GV974BK-T	Black 	1/cs

FEATURES

- Fully assembled out of the box.
- Attractive, Euro-style design.
- Brake cable inside frame for added safety.
- Handle height easily adjusts with unique push button tool free.
- Back support height easily adjusts to meet patients needs tool-free.
- Caster fork design enhances turning radius.
- Large 10" front casters allow optimal steering and rolling comfort.
- Curb hopper for terrain transitions.
- Easily folds with one hand to ultra-compact size for storage.
- Cross-brace design allows for side-to-side folding and added stability.
- Seat is durable and comfortable.
- Handsome and removable zippered storage bag with unique attachment to keep bag securely in place when F-22 is open or folded.
- Lightweight, aluminum frame.
- Manufactured in Taiwan, TAA Compliant.

Specifications

Handle Height	36" - 41"
Unit Dimensions	23" (W) x 27.75" (L)
Seat Dimensions	18" (W) x 23.5" (H) x 8.5" (D)
Seat Height	23.6"
Weight	17.5 lbs.
Carton Shipping Weight	30 lbs.
Weight Limit	300 lbs.
Warranty (Frame)	Limited Lifetime
Warranty (Brake and Handle Housing)	5 Year Limited

F-22 HD ROLLATOR



Item #	Color	UOM
GV974HD-BR	Black ●	1/cs

FEATURES

- Fully assembled out of the box.
- Attractive, Euro-style design.
- Brake cable inside frame for added safety.
- Handle height easily adjusts with unique push button tool free.
- Back support height easily adjusts to meet patients needs tool-free.
- Caster fork design enhances turning radius.
- Large 10" front casters allow optimal steering and rolling comfort.
- Easily folds with one hand to ultra-compact size for storage.
- Curb hopper for terrain transitions.
- Cross-brace design allows for side-to-side folding and added stability.
- Seat is durable and comfortable.
- Handsome and removable zippered storage bag with unique attachment to keep bag securely in place when F-22 HD is open or folded.
- Lightweight, aluminum frame.
- Manufactured in Taiwan, TAA Compliant.
- Based on customer feedback, we have changed the frame color to black.
- This is a running change and there is no change of the order code GV974HD-BR.
- The only change is the frame color, no other features or specifications have changed.

Specifications

Handle Height	36" - 41"
Unit Dimensions	28" (L) x 28.25" (W)
Seat Dimensions	23" (W) x 10" (D)
Seat Height	23.5"
Weight	20.5 lbs.
Weight Limit	450 lbs.
Warranty (Frame)	Limited Lifetime
Warranty (Brake and Handle Housing)	5 Year Limited