

M3 WHEELCHAIR

Flexibility and Convenience Combined



The sturdy, lightweight M3 wheelchair is an ideal choice for flexibility and convenience. With a durable carbon steel frame and adjustable components, the M3 helps patients function at their best and can be customized to their changing needs.



Available with
Treaded Tires

FLEXIBILITY & CONVENIENCE

BUILT-IN SEAT RAIL EXTENSIONS
AND EXTENDABLE UPHOLSTERY
EASILY ADJUST SEAT DEPTH
FROM 16" TO 18"



DUAL AXLE POSITIONS
PROVIDES EASY TRANSITION
OF SEAT HEIGHT TO THE HEMI
POSITION

ANTI-TIPPERS INCLUDED



FOR MORE INFORMATION



Visit Us at
www.DriveMedical.com



Call Us Today!
877.224.0946



Contact Your Sales
Representative

FEATURES:

- 300lbs. Weight Capacity
- Available in 16", 18" & 20" seat widths
- Nylon upholstery is durable, lightweight, and easy to clean
- Carbon steel frame with attractive finish
- Swing-away footrests or elevating legrests with tool-free adjustable length riggings
- Removable flip-back armrests available in desk or full length padded arms
- Composite, mag-style wheels are lightweight and maintenance free
- Precision sealed wheel bearings in front and rear ensure long lasting performance and reliability
- Easy to propel and transport
- Push-to-lock wheel locks
- Country of Origin- Taiwan TAA Compliant
- Limited Lifetime Warranty

ITEM #	DESCRIPTION	UOM
GVM316DDA-SF	16" Light Weight Wheelchair, Detachable Desk Arm, Swing Away Foot Rests	1/ea
GVM316DFA-SF	16" Light Weight Wheelchair, Detachable Full Arm, Swing Away Foot Rests	1/ea
GVM318DDA-SF	18" Light Weight Wheelchair, Detachable Desk Arm, Swing Away Foot Rests	1/ea
GVM318DFA-SF	18" Light Weight Wheelchair, Detachable Full Arm, Swing Away Foot Rests	1/ea
GVM320DDA-SF	20" Light Weight Wheelchair, Detachable Desk Arm, Swing Away Foot Rests	1/ea
GVM320DFA-SF	20" Light Weight Wheelchair, Detachable Full Arm, Swing Away Foot Rests	1/ea
GVM316DDA-ELR	16" Light Weight Wheelchair, Detachable Desk Arm, Elevating Legrests	1/ea
GVM316DFA-ELR	16" Light Weight Wheelchair, Detachable Full Arm, Elevating Legrests	1/ea
GVM318DDA-ELR	18" Light Weight Wheelchair, Detachable Desk Arm, Elevating Legrests	1/ea
GVM318DFA-ELR	18" Light Weight Wheelchair, Detachable Full Arm, Elevating Legrests	1/ea
GVM320DDA-ELR	20" Light Weight Wheelchair, Detachable Desk Arm, Elevating Legrests	1/ea
GVM320DFA-ELR	20" Light Weight Wheelchair, Detachable Full Arm, Elevating Legrests	1/ea
GVM316DDA-SF-T	16" Light Weight Wheelchair, Detachable Desk Arm, Swing Away Foot Rests, Treaded Tires	1/ea
GVM318DDA-SF-T	18" Light Weight Wheelchair, Detachable Desk Arm, Swing Away Foot Rests, Treaded Tires	1/ea
GVM320DDA-SF-T	20" Light Weight Wheelchair, Detachable Desk Arm, Swing Away Foot Rests, Treaded Tires	1/ea
GVSTDS1034	Cane Holder	1/ea
GVSTDS601	6" Wheel Lock Extension	1 pair
GVSTDS803	Oxygen Tank Holder E tank	1/ea
GVSTDS819	Anti-tippers	1 pair
GVLK3JELR	Elevating Legrests	1 pair