



Combine True Low Air Loss and Pulsation Therapy with One Mattress System!

drive Mason LS9000 10" True Low Air Loss Mattress System with Pulsation

Item #	Product Description	Product Dimensions	Product HCPCS	UOM
LS9000	True Low Air Loss Mattress System	80" (L) x 36" (W) x 10" (H)	E0277	1/ea

- The Masonair™ LS9000 system offers an extraordinary therapeutic mattress system for the prevention and treatment of pressure ulcers.
- The Powerloft Pro 900 pump produces 1180 liters of air flow per minute offering a cool and dry sleeping surface.
- The 10 Adjustable Comfort Settings allows pressure to be customized by patient's weight.
- The Auto-Firm feature allows for a 1 minute fill rate.
- The exclusive Pulsate mode reduces the air flow every 30 seconds to 50% of the comfort setting. This feature will intermittently reduce pressure and increase blood flow.
- The internal pressure control sensor continuously monitors surface pressure for optimal pressure redistribution.
- For use in the Homecare, Long Term Care, and Acute Care Markets.
- 2" convoluted foam base in case of power outage.
- CPR feature allows for rapid deflation of mattress system
- Static function can suspend the pulsate mode.
- Visible and audible alarms with Alarm Mute feature (3 minute auto-return feature)
- CPR valve allows for rapid deflation.

specifications

Pump Power	110VAC / 60 Hz
Pump Airflow	1180 LPM
Pump Cycle Time	30 Second Pulsation
Pump Type	Blower
Pump Alarms	Visual/Audible
Mattress Dimensions	80" x 36" x 10"
Mattress Cover Material	Fluid-Resistant Nylon
Weight Capacity	1000 lbs.
Warranty	1 Year Limited