

## Premier One General Use Foam Cushion

HCPCS: E260



RETAIL BOX



M8061-PB

- The Premier One General Use Foam cushion is designed to provide comfort and stability with a simple, economical design
- Aerated foam helps keep the user cool and dry
- The low-shear water-resistant stretch cover helps maintain sensitive skin integrity
- Convenient 12" tie straps secure the cushion to the wheelchair
- The waterproof vinyl bottom protects the wheelchair seat
- Meets flammability standard TB117

PRODUCTS	Description
M8061-PB	14X14X3 PREMIER ONE FOAM CUSHION, E2601, 1/ea, 4/cs
M8062-PB	14X16X3 PREMIER ONE FOAM CUSHION, E2601, 1/ea, 4/cs
M8063	Premier One General Use Foam Cushion, E2601, 1/ea, 4/cs
M8064	Premier One General Use Foam Cushion, E2601, 1/ea, 4/cs
M8065	Premier One General Use Foam Cushion, E2601, 1/ea, 4/cs
M8070	Premier One General Use Foam Cushion, E2601, 1/ea, 4/cs
M8077	Premier One General Use Foam Cushion, E2602, 1/ea, 4/cs
M8078	Premier One General Use Foam Cushion, E2602, 1/ea, 4/cs
M8079	Premier One General Use Foam Cushion, E2602, 1/ea, 4/cs
M8080	PremierOneBarFoamCushion24x20x31, E2602, 1/ea, 4/cs

SPECIFICATIONS	Description
Cover Material	Fluid-Resistant Stretch Top   Waterproof Vinyl Bottom
Weight Capacity	250 lbs.   400 lbs.
Warranty	1 Year Limited