



YOUR
HOME.
**YOUR
WAY.**



Home Medical Equipment.

 **MOBILITY**

 **SAFETY**

 **COMFORT**



YOUR FIRST CHOICE IN HOME MEDICAL EQUIPMENT



Please stop by or call to ask about other products not featured in this catalog.

Town Name Location 1

123 Main Street
Anytown, IL
00000
T: 800-123-4567
F: 800-987-6543
www.sitename.com

Town Name Location 2

123 Main Street
Anytown, IL
00000
T: 800-123-4567
F: 800-987-6543
www.sitename.com

Town Name Location 3

123 Main Street
Anytown, IL
00000
T: 800-123-4567
F: 800-987-6543
www.sitename.com

Table of Contents

- 7 Wheelchairs
- 12 Walking Aids
- 19 Bath & Toilet
- 26 Beds
- 27 Aids for Daily Living
- 30 Bed Accessories

Our Staff

Duis aute irure dolor in reprehenderit in voluptate velit esse cillum dolore eu fugiat nulla pariatur. Excepteur sint laborum.

Our Knowledge

Duis aute irure dolor in reprehenderit in voluptate velit esse cillum dolore eu fugiat nulla pariatur. Excepteur sint laborum.

Our Location

Duis aute irure dolor in reprehenderit in voluptate velit esse cillum dolore eu fugiat nulla pariatur. Excepteur sint laborum.



ENJOY YOUR HOME SAFELY AND COMFORTABLY.

With a few adjustments and the right equipment you can continue to live your life in the home you love. Medline offers a wide range of options that can help you stay safe, stay mobile and be comfortable—in every room of your home.

Bed Assist Bar pg. 30


Get out of bed more comfortably

Bath Bench pgs. 20–22

Avoid falling in the bath

Bedside Commode pgs. 21 & 25

Make bathroom trips easier



Explore our home diagram to see how Medline can help you live your life your way.

 **MOBILITY**

 **SAFETY**

 **COMFORT**

Transfer Bench pgs. 21-22

Get in and out of the tub with ease

Ultra Reacher pg. 29


Reach household items without straining

Rollator pgs. 13-16

Walk on your own with the ability to sit and rest

Transport Chair pgs. 7-8

Have a lightweight wheelchair option



KEEP MOVING.

Maintaining your mobility means maintaining your independence—sometimes you just need a little help. Medline offers an extensive selection of walkers, rolling walkers, canes and more on many of your favorite retail websites.



HYBRID 2

- » Combination wheelchair and transport chair
- » Two chairs in one
- » Rear wheels remove with just the push of a button
- » Weighs only 33 lbs. (15 kg), 24 lbs. (11 kg) when in transport mode
- » Converts easily in seconds
- » Comfortable nylon upholstery
- » Includes anti-tippers
- » 300-lb. (136 kg) weight capacity



In transport mode

MDS806250H2

Item Number	Arm Style	Leg Style	Seat W x D: Inch (cm)
MDS806250NH2	Desk Length Removable	Swing-Away Detachable Footrest	16 x 16 (41 x 41)
MDS806250H2	Desk Length Removable	Swing-Away Detachable Footrest	18 x 16 (46 x 41)
MDS806300H2	Desk Length Removable	Swing-Away Detachable Elevating Legrest	18 x 16 (46 x 41)

FREEDOM 2

- » Lightest transport chair on the market!
- » Back folds down for easy storage and transport
- » Back side features bottle holder, key ring, cell phone pocket, identification holder and pocket
- » Comfortable nylon upholstery
- » Convenient side bag for storage
- » Restaurant-style perm. armrests
- » Convenient cup holder
- » Seat belt for safety
- » Tool-free adjustable, detachable footrests
- » Patented stow-n-go clips to safely lock footrests in place
- » Weighs less than 15 lbs. (7 kg)
- » 300-lb. (136 kg) weight capacity



MDS808200F2R

Item Number	Color	Seat W x D: Inch (cm)
MDS808200F2R	Red	19 x 16 (48 x 41)
MDS808200F2B	Blue	19 x 16 (48 x 41)
MDS808200F2P	Pink	19 x 16 (48 x 41)
MDS808200F2S	Metallic Silver	19 x 16 (48 x 41)

TRANSLATOR

- » Converts in seconds, just flip the back to instantly convert from transport chair mode to rollator/walker mode
- » Easy-to-adjust push-button footrests lock to sides of translator when not in use
- » Breathable nylon upholstery
- » Height-adjustable push grips
- » Comfortable hand brake
- » Restaurant-style permanent armrests
- » Convenient side carrying case
- » Convenient cup holder
- » Roomy, under-seat basket
- » Strong carrying handles
- » 250-lb. (113 kg) weight capacity



MDS808200TR

Item Number	Color	Leg Style	Seat W x D: Inch (cm)
MDS808200TRR	Red	Swing-Away Detachable Footrest	19 x 16 (48 x 41)
MDS808200TR	Blue	Swing-Away Detachable Footrest	19 x 16 (48 x 41)

EXCEL® K1 BASIC

- » Elevating legrest models feature a notched, stainless steel ratchet bar to lock the legrests securely in place
- » Durable tig-welded frame with silver vein finish
- » Comfortable nylon upholstery
- » Smooth-rolling solid, flat-free tires
- » Hemi-height adjustable axles raise/lower seat 2" (5 cm)
- » Weighs 38 lbs. (17 kg)
- » 300-lb. (136 kg) weight capacity



MDS806450EE

Item Number	Arm Style	Leg Style	Seat W x D: Inches (cm)
MDS806250NEE	Desk Length Removable	Detachable Footrest	16 x 16 (41 x 41)
MDS806300NEE	Desk Length Removable	Detachable Elevating Legrest	16 x 16 (41 x 41)
MDS806150EE	Full Length Permanent	Detachable Footrest	18 x 16 (46 x 41)
MDS806200EE	Full Length Permanent	Detachable Elevating Footrest	18 x 16 (46 x 41)
MDS806250EE	Desk Length Removable	Detachable Footrest	18 x 16 (46 x 41)
MDS806300EE	Desk Length Removable	Detachable Elevating Legrest	18 x 16 (46 x 41)
MDS806400EE	Desk Length Removable	Detachable Footrest	20 x 16 (51 x 41)
MDS806450EE	Desk Length Removable	Detachable Elevating Legrest	20 x 16 (51 x 41)

EXCEL® K2B

- » Strong vinyl upholstery
- » Durable tig-welded frame with silver vein finish
- » Padded, upholstered armrests and calf pads are standard
- » Hemi-height adjustable axles raise/lower seat 2" (5 cm)
- » Convenient chart pocket on the back
- » Weighs 38 lbs. (17 kg)
- » 300-lb. (136 kg) weight capacity



MDS806250EV

Item Number	Arm Style	Leg Style	Seat W x D: Inches (cm)
MDS806250NEV	Desk Length Removable	Swing Away	16 x 16 (41 x 41)
MDS806300NEV	Desk Length Removable	Elevating Legrest	16 x 16 (41 x 41)
MDS806150EV	Full Length Permanent	Swing Away	18 x 16 (46 x 41)
MDS806200EV	Full Length Permanent	Elevating Legrests	18 x 16 (46 x 41)
MDS806250EV	Desk Length Removable	Swing Away	18 x 16 (46 x 41)
MDS806300EV	Desk Length Removable	Elevating Legrest	18 x 16 (46 x 41)
MDS806400EV	Desk Length Removable	Swing Away	20 x 16 (51 x 41)
MDS806450EV	Desk Length Removable	Elevating Legrest	20 x 16 (51 x 41)

K3 BASIC

- » Elevating legrest models feature a notched, stainless steel ratchet bar to lock the legrests securely in place
- » Durable, tig-welded frame with silver vein finish
- » Comfortable nylon upholstery
- » Padded, upholstered armrests and calf pads are standard
- » Hemi-height adjustable axles raise/lower seat 2" (5 cm)
- » Weighs only 35 lbs. (16 kg)
- » 300-lb. (136 kg) weight capacity




MDS806600E

Item Number	Arm Style	Leg Style	Seat W x D: Inch (cm)
MDS806600NE	Desk Length Removable	Detachable Footrest	16 x 16 (41 x 41)
MDS806650NE	Desk Length Removable	Detachable Elevating Legrest	16 x 16 (41 x 41)
MDS806600E	Desk Length Removable	Detachable Footrest	18 x 16 (46 x 41)
MDS806650E	Desk Length Removable	Detachable Elevating Legrest	18 x 16 (46 x 41)
MDS806560E*	Swing Back Desk Length Removable	Detachable Footrest	20 x 18 (51 x 46)
MDS806565E*	Swing Back Desk Length Removable	Detachable Elevating Legrest	20 x 18 (51 x 46)

*Labeled as K4 basic

K3 LIGHTWEIGHT

- » Elevating legrest models feature a notched, stainless steel ratchet bar to lock the legrests securely in place
- » Durable, tig-welded frame with black powder coat finish
- » Comfortable nylon upholstery
- » Padded and upholstered armrests and calf pads are standard
- » Threaded seat screw inserts reduce stripping from upholstery changes
- » Hemi-height adjustable axles raise/lower seat 2" (5 cm)
- » Weighs only 33 lbs. (15 kg)
- » 300-lb. (136 kg) weight capacity



MDS806650

Item Number	Arm Style	Leg Style	Seat W x D: Inch (cm)
MDS806600N	Desk Length Removable	Detachable Footrest	16 x 16 (41 x 41)
MDS806650N	Desk Length Removable	Detachable Elevating Legrest	16 x 16 (41 x 41)
MDS806600	Desk Length Removable	Detachable Footrest	18 x 16 (46 x 41)
MDS806650	Desk Length Removable	Detachable Elevating Legrest	18 x 16 (46 x 41)
MDS806600FLA	Full Length Removable	Detachable Footrest	18 x 16 (46 x 41)
MDS806650FLA	Full Length Removable	Detachable Elevating Legrest	18 x 16 (46 x 41)

K4 BASIC

- » Adjustable-height back
- » Swing-back arms with easy flip back arm release
- » Elevating legrest models feature a notched, stainless steel ratchet bar to lock the legrests securely in place
- » Durable, tig-welded frame with gray powder coat finish
- » Comfortable nylon upholstery
- » Tool-free push button adjustable footrests
- » Hemi-height adjustable axles raise/lower seat 2" (5 cm)
- » Smooth-rolling solid flat-free latex tires
- » Weighs only 33 lbs. (15 kg)
- » 300-lb. (136 kg) weight capacity



MDS806550E

Item Number	Arm Style	Leg Style	Seat W x D: Inch (cm)
MDS806500NE	Swing Back Desk Length Removable	Detachable Footrest	16 x 16 (41 x 41)
MDS806550NE	Swing Back Desk Length Removable	Detachable Elevating Legrest	16 x 16 (41 x 41)
MDS806500E	Swing Back Desk Length Removable	Detachable Footrest	18 x 16 (46 x 41)
MDS806550E	Swing Back Desk Length Removable	Detachable Elevating Legrest	18 x 16 (46 x 41)
MDS806560E	Swing Back Desk Length Removable	Detachable Footrest	20 x 18 (51 x 46)
MDS806565E	Swing Back Desk Length Removable	Detachable Elevating Legrest	20 x 18 (51 x 46)

K4 LIGHTWEIGHT

- » Swing back arms with easy flip back arm release
- » Quick-release axles allow tool-free removal of rear wheels
- » Adjustable height back
- » Elevating legrest models feature a notched, stainless steel ratchet bar to lock the legrests securely in place
- » Padded and upholstered armrests and calf pads are standard
- » Dual axle hemi-height adjustable
- » Tough nylon upholstery
- » PLUS models come standard with anti-tippers
- » Weighs only 31 lbs. (14 kg)
- » 300- to 350-lb. (136–159 kg) weight capacity



MDS806550

Item Number	Arm Style	Leg Style	Seat W x D: Inch (cm)
MDS806500N	Swing Back Desk Length Removable	Swing-Away Detachable Footrest	16 X 16 (41 x 41)
MDS806550N	Swing Back Desk Length Removable	Swing-Away Detachable Elevating Legrest	16 X 16 (41 x 41)
MDS806500	Swing Back Desk Length Removable	Swing-Away Detachable Footrest	18 X 16 (46 x 41)
MDS806550	Swing Back Desk Length Removable	Swing-Away Detachable Elevating Legrest	18 X 16 (46 x 41)
MDS806500PLUS	Swing Back Height Adjustable Desk Length	Swing-Away Detachable Footrest	18 X 16 (46 x 41)
MDS806550PLUS	Swing Back Height Adjustable Desk Length	Swing-Away Detachable Elevating Legrest	18 X 16 (46 x 41)
MDS806560	Swing Back Desk Length Removable	Swing-Away Detachable Footrest	20 X 18 (51 x 46)
MDS806560PLUS	Swing Back Height Adjustable Desk Length	Swing-Away Detachable Footrest	20 X 18 (51 x 46)
MDS806565	Swing Back Desk Length Removable	Swing-Away Detachable Elevating Legrest	20 X 18 (51 x 46)
MDS806570	Swing Back Desk Length Removable	Swing-Away Detachable Footrest	22 X 18 (56 x 46)
MDS806575	Swing Back Desk Length Removable	Swing-Away Detachable Elevating Legrest	22 X 18 (56 x 46)

EXCEL® KIDZ PEDIATRIC



- » Adjustable telescoping handles for caregiver comfort
- » Durable vinyl upholstery
- » Treaded, flat-free tires for safe riding
- » Anti-tippers included
- » Standard dual axle allows for easy adjustment of seat height to hemi level
- » Padded and upholstered armrests and calf pads come standard
- » Flip-back arms
- » 250-lb. (113 kg) weight capacity



MDS806140PD



Item Number	Arm Style	Leg Style	Seat W x D: Inch (cm)
MDS806140PD	Swing Back Desk Length Permanent Arms	Detachable Elevating Legrests	14 x 12 (36 x 30)

HOW TO ADJUST YOUR MOBILITY AIDS



Rollators

1. Select a rollator that is appropriate to the patient's weight.
2. Measure from the ground to the crease in the back of the knee. (See A)
3. Adjust the rollator legs accordingly so the seat is at same height as A.
4. Keeping arms at your side with a natural bend at the elbow, measure from the ground to the center of the wrist. (See B)




Walkers

1. Select a walker that is appropriate to the patient's weight.
2. Keeping arms at your side with a natural bend at the elbow, measure from the ground to the center of the wrist. (See B)
3. Adjust all four legs to the same height so that the handgrips match the same height as the wrist.

Canes

1. Select a cane that is appropriate to the patient's weight.
2. For the ideal cane height, keep arms at your side with a natural bend at the elbow. Measure from the ground to the center of the wrist. (See B)
3. Loosen the collar on the cane and press the push-button lock, slide the lower tube into or out of the body of the cane until you get to the wrist-height measurement, and let the push button "click" into the closest hole to the measurement. Tighten the collar (check and re-tighten the collar regularly).



Axillary Crutches

1. Select a pair of crutches that is appropriate to the patient's weight.
2. For the overall crutch height, stand straight and measure from the ground to two to three inches below the axilla. The cushion should not be touching the axilla. (See C)
3. For hand grip height, keep arms at your side with a natural bend at the elbow, measure from the ground to the center of the wrist. (See B)

EMPOWER™ ROLLATOR

- » Seat flips up for better accessibility to counters and sinks, and for walking inside the unit to promote better gait
- » Folds completely flat for convenient storage and transport
- » Backrest flips forward and backward for a more comfortable seating area
- » All adjustment points are yellow for easy identification
- » Features a roomy memory foam seat, storage pouch under the seat, tray organizer and folding cup holder
- » 300-lb. (136 kg) weight capacity



MDS86845B

Item Number	Color	Wt. Cap.: Lbs. (kg)	Wheels: Inch (cm)	Distance Between Handles: Inch (cm)	Seat Height Min/Max: Inch (cm)	Height Min/Max: Inch (cm)	Product Weight: Lbs. (kg)
MDS86845B	Blue	300 (136)	8 (20)	17.5 (44)	20-25 (51-64)	31-37 (79-94)	17 (8)
MDS86845BLK	Black	300 (136)	8 (20)	17.5 (44)	20-25 (51-64)	31-37 (79-94)	17 (8)
MDS86845S	Silver	300 (136)	8 (20)	17.5 (44)	20-25 (51-64)	31-37 (79-94)	17 (8)

FREEDOM ULTRA LIGHT ROLLATORS

- » Weighing only 11 lbs., the Freedom rollator is the lightest rollator on the market
- » Lightweight design makes this the ideal model for travel
- » Height adjustable seat and arms allow for a truly custom fit
- » Accommodates users from 4'11"–6'4" (150–193 cm)
- » 250-lb. (113 kg) weight capacity



MDS86825SLR

Item Number	Color	Wt. Cap.: Lbs. (kg)	Wheels: Inch (cm)	Distance Between Handles: Inch (cm)	Seat Height: Min/Max: Inch (cm)	Height Min/Max: Inch (cm)	Product Weight: Lbs. (kg)
MDS86825SLR	Burgundy	250 (113)	6 (15)	17 (43)	18-23 (46-58)	29-36 (74-91)	11 (5)
MDS86825SLB	Smoky Blue	250 (113)	6 (15)	17 (43)	18-23 (46-58)	29-36 (74-91)	11 (5)
MDS86825SLBLK	Black	250 (113)	6 (15)	17 (43)	18-23 (46-58)	29-36 (74-91)	11 (5)
MDS86825SLT	Titanium	250 (113)	6 (15)	17 (43)	18-23 (46-58)	29-36 (74-91)	11 (5)
MDS86825SLWT	White	250 (113)	6 (15)	17 (43)	18-23 (46-58)	29-36 (74-91)	11 (5)

POSH ROLLATOR

- » Trendy pink frame with chic zebra-print front pouch and under-seat bag
- » Lightweight—only 11 lbs. (5 kg)
- » Features a rearview mirror
- » Durable, reinforced frame that's 40% lighter than others on the market
- » Height adjustability on handles and seat allows for a truly custom fit
- » 250-lb. (113 kg) weight capacity



MDS86835SHE

Item Number	Color	Wt. Cap.: Lbs. (kg)	Wheels: Inch (cm)	Distance Between Handles: Inch (cm)	Seat Height Min/Max: Inch (cm)	Height Min/Max: Inch (cm)	Product Weight: Lbs. (kg)
MDS86835SHE	Pink/Zebra	250 (113)	6 (15)	17 (43)	18–23 (46–58)	29–36 (74–91)	11 (5)

STEALTH ROLLATOR

- » Designed for a rugged, stylish man
- » Roomy matte black frame with 7" (18 cm) chrome hubs
- » Comes standard with a cup holder, a rear-view mirror, a utility belt and a flashlight
- » 300-lb. (136 kg) weight capacity



MDS86835HE

Item Number	Color	Wt. Cap.: Lbs. (kg)	Wheels: Inch (cm)	Distance Between Handles: Inch (cm)	Seat Height Min/Max: Inch (cm)	Height Min/Max: Inch (cm)	Product Weight: Lbs. (kg)
MDS86835HE	Matte Black	300 (136)	7 (18)	19 (48)	19–25.5 (48–65)	30.5–37.5 (77–95)	17 (8)

SIMPLICITY® ROLLATOR

- » Sleek, sporty, European design
- » Compact—so it's easy to take anywhere
- » Folds up with one hand
- » Stow-n-go clip keeps it securely closed during transport
- » Features a roomy, removable storage bag, reflective tape for safety, a bell and cup and cane holder
- » 300-lb. (136 kg) weight capacity



MDS86835EURO

Item Number	Color	Wt. Cap.: Lbs. (kg)	Wheels: Inch (cm)	Distance Between Handles: Inch (cm)	Seat Height: Inch (cm)	Height Min/Max: Inch (cm)	Product Weight: Lbs. (kg)
MDS86835EURO	Red	300 (136)	8 (20)	17.5 (44)	20 (51)	32.5-37.5 (83-95)	17 (8)
MDS86835BEURO	Teal	300 (136)	8 (20)	17.5 (44)	20 (51)	32.5-37.5 (83-95)	17 (8)

BASIC ROLLATORS

- » Provides all the necessary features at an affordable price
- » Folds easily for compact storage
- » Features a convenient storage bag
- » 250-lb. (113 kg) weight capacity



MDS86850E

Item Number	Color	Wt. Cap.: Lbs. (kg)	Wheels: Inch (cm)	Distance Between Handles: Inch (cm)	Seat Height: Inch (cm)	Height Min/Max: Inch (cm)	Product Weight: Lbs. (kg)
MDS86850E	Burgundy	250 (113)	6 (15)	17.5 (44)	20.5 (52)	31-35 (79-89)	14 (6)
MDS86850EB	Blue	250 (113)	6 (15)	17.5 (44)	20.5 (52)	31-35 (79-89)	14 (6)
MDS86850EP	Purple	250 (113)	6 (15)	17.5 (44)	20.5 (52)	31-35 (79-89)	14 (6)
MDS86850EG	Green	250 (113)	6 (15)	17.5 (44)	20.5 (52)	31-35 (79-89)	14 (6)
MDS86850EO	Orange	250 (113)	6 (15)	17.5 (44)	20.5 (52)	31-35 (79-89)	14 (6)

BASIC STEEL ROLLATORS

- » Our most affordable rollator yet
- » Constructed of a durable powder coated steel
- » Features a convenient storage bag
- » Available with 6" (15 cm) or 8" (20 cm) wheels
- » 350-lb. (159 kg) weight capacity



MDS86850EBS

Item Number	Color	Wt. Cap.: Lbs. (kg)	Wheels: Inch (cm)	Distance Between Handles: Inch (cm)	Seat Height: Inch (cm)	Height Min/Max: Inch (cm)	Product Weight: Lbs. (kg)
MDS86850ES	Red	350 (159)	6 (15)	17.5 (44)	20.5 (52)	31-35 (79-89)	19 (9)
MDS86850EBS	Blue	350 (159)	6 (15)	17.5 (44)	20.5 (52)	31-35 (79-89)	19 (9)
MDS86850EGS	Green	350 (159)	6 (15)	17.5 (44)	20.5 (52)	31-35 (79-89)	19 (9)
MDS86850ES8	Red	350 (159)	8 (20)	17.5 (44)	20.5 (52)	31-35 (79-89)	19 (9)
MDS86850EBS8	Blue	350 (159)	8 (20)	17.5 (44)	20.5 (52)	31-35 (79-89)	19 (9)
MDS86850EGS8	Green	350 (159)	8 (20)	17.5 (44)	20.5 (52)	31-35 (79-89)	19 (9)

TWO-BUTTON FOLDING WALKERS

- » Easy-to-use, push-button mechanisms let walkers fold
- » Each side folds independently to allow easy movement through narrow spaces
- » Side cross braces provide extra durability
- » Push-button, height-adjustable legs
- » Replaceable tips



MDS864104

Item Number	Description	Wheels: Inch (cm)	Wt. Cap.: Lbs. (kg)	Approximate User Height: Inch (cm)	Inside Hand Grip Width: Inch (cm)	Height Adjustment: Inch (cm)	Folded Width: Inch (cm)
MDS864104	Adult	None	300 (136)	5'5"-6'4" (165-193)	16.75 (43)	32-39 (81-99)	4.5 (11)
MDS86410W54	Adult	5 (13)	300 (136)	5'5"-6'4" (165-193)	16.75 (43)	32.5-39.5 (82-100)	4.5 (11)
MDS86410W4	Adult	3 (8)	300 (136)	5'5"-6'4" (165-193)	16.75 (43)	32-39 (81-99)	4.5 (11)
MDS86410J4	Junior	None	300 (136)	4'6"-5'5" (137-165)	16.75 (43)	25-32 (64-81)	4.5 (11)
MDS86410JW54	Junior	5 (13)	300 (136)	4'6"-5'5" (137-165)	16.75 (43)	25.5-32.5 (65-83)	4.5 (11)
MDS86410WJ4	Junior	3 (8)	300 (136)	4'6"-5'5" (137-165)	16.75 (43)	25-32 (64-81)	4.5 (11)
MDS864104B	Basic Adult	None	400 (181)	5'5"-6'4" (165-193)	16.75 (43)	32-39 (81-99)	4.5 (11)
MDS86410W54B	Basic Adult	5 (13)	400 (181)	5'5"-6'4" (165-193)	16.75 (43)	32.5-39.5 (82-100)	4.5 (11)
MDS86410W4B	Basic Adult	3 (8)	400 (181)	5'5"-6'4" (165-193)	16.75 (43)	32-39 (81-99)	4.5 (11)
MDS86410J4B	Basic Junior	None	400 (181)	4'6"-5'5" (137-165)	16.75 (43)	25-32 (64-81)	4.5 (11)
MDS86410JW54B	Basic Junior	5 (13)	400 (181)	4'6"-5'5" (137-165)	16.75 (43)	25.5-32.5 (65-83)	4.5 (11)
MDS86410WJ4B	Basic Junior	3 (8)	400 (181)	4'6"-5'5" (137-165)	16.75 (43)	25-32 (64-81)	4.5 (11)



HOME CARE PATIENT AIDS

MOBILITY

SAFETY

COMFORT

Patient Name _____

Additional Equipment Needs: _____

Date of Discharge _____

Patient Medicare/Medicaid # _____

Name of Therapist or Institution _____

Authorized by _____

Phone # _____

Optional Information _____

MOBILITY

Hybrid 2 (pg. 7)

- MDS806250NH2 Desk Length Removable Arms
Swing-Away Detachable Footrests
- MDS806250H2 Desk Length Removable Arms
Swing-Away Detachable Footrests
- MDS806300H2 Desk Length Removable Arms
Swing-Away Detachable Elevating Legrests

Freedom 2 (pg. 7)

- MDS808200F2R Red
- MDS808200F2B Blue
- MDS808200F2P Pink
- MDS808200F2S Metallic Silver

Translator (pg. 8)

- MDS808200TRR Red Swing-Away Detachable Footrests
- MDS808200TR Blue Swing-Away Detachable Footrests

Excel® K1 Basic (pg. 8)

- MDS806250NEE Desk Length Removable Arms
Detachable Footrests
- MDS806300NEE Desk Length Removable Arms
Detachable Elevating Legrests
- MDS806150EE Full Length Permanent Arms
Detachable Footrests
- MDS806200EE Full Length Permanent Arms
Detachable Elevating Footrests
- MDS806250EE Desk Length Removable Arms
Detachable Footrests
- MDS806300EE Desk Length Removable Arms
Detachable Elevating Legrests
- MDS806400EE Desk Length Removable Arms
Detachable Footrests

- MDS806450EE Desk Length Removable Arms
Detachable Elevating Legrests

Excel® K2B (pg. 9)

- MDS806250NEV Desk Length Removable Arms
Swing-Away Legrests
- MDS806300NEV Desk Length Removable Arms
Elevating Legrests
- MDS806150EV Full Length Permanent Arms
Swing-Away Legrests
- MDS806200EV Full Length Permanent Arms
Elevating Legrests
- MDS806250EV Desk Length Removable Arms
Swing-Away Legrests
- MDS806300EV Desk Length Removable Arms
Elevating Legrests
- MDS806400EV Desk Length Removable Arms
Swing-Away Legrests
- MDS806450EV Desk Length Removable Arms
Elevating Legrests

K3 Basic (pg. 9)

- MDS806600NE Desk Length Removable Arms
Detachable Footrests
- MDS806650NE Desk Length Removable Arms
Detachable Elevating Legrests
- MDS806600E Desk Length Removable Arms
Detachable Footrests
- MDS806650E Desk Length Removable Arms
Detachable Elevating Legrests
- MDS806560E Swing Back Desk Length Removable Arms
Detachable Footrests
- MDS806565E Swing Back Desk Length Removable Arms
Detachable Elevating Legrests

K3 Lightweight (pg. 10)

- MDS806600N Desk Length Removable Arms Detachable Footrests
- MDS806650N Desk Length Removable Arms Detachable Elevating Legrests
- MDS806600 Desk Length Removable Arms Detachable Footrests
- MDS806650 Desk Length Removable Arms Detachable Elevating Legrests
- MDS806600FLA Full Length Removable Arms Detachable Footrests
- MDS806650FLA Full Length Removable Arms Detachable Elevating Legrests

K4 Basic (pg. 10)

- MDS806500NE Swing Back Desk Length Removable Arms Detachable Footrests
- MDS806550NE Swing Back Desk Length Removable Arms Detachable Elevating Legrests
- MDS806500E Swing Back Desk Length Removable Arms Detachable Footrests
- MDS806550E Swing Back Desk Length Removable Arms Detachable Elevating Legrests
- MDS806560E Swing Back Desk Length Removable Arms Detachable Footrests
- MDS806565E Swing Back Desk Length Removable Arms Detachable Elevating Legrests

K4 Lightweight (pg. 11)

- MDS806500N Swing Back Desk Length Removable Arms Swing-Away Detachable Footrests
- MDS806550N Swing Back Desk Length Removable Arms Swing-Away Detachable Elevating Legrests
- MDS806500 Swing Back Desk Length Removable Arms Swing-Away Detachable Footrests
- MDS806550 Swing Back Desk Length Removable Arms Swing-Away Detachable Elevating Legrests
- MDS806500PLUS Swing Back Height Adjustable Desk Length Arms Swing-Away Detachable Footrests
- MDS806550PLUS Swing Back Height Adjustable Desk Length Arms Swing-Away Detachable Elevating Legrests
- MDS806560 Swing Back Desk Length Removable Arms Swing-Away Detachable Footrests
- MDS806560PLUS Swing Back Height Adjustable Desk Length Arms Swing-Away Detachable Footrests
- MDS806565 Swing Back Desk Length Removable Arms Swing-Away Detachable Elevating Legrests
- MDS806570 Swing Back Desk Length Removable Arms Swing-Away Detachable Footrests
- MDS806575 Swing Back Desk Length Removable Arms Swing-Away Detachable Elevating Legrests

Excel® Kidz Pediatric (pg. 11)

- MDS806140PD Swing Back Desk Length Permanent Arms Detachable Elevating Legrests

Empower™ Rollator (pg. 13)

- MDS86845B Blue
- MDS86845BLK Black
- MDS86845S Silver

Freedom Ultra Light Rollators (pg. 13)

- MDS86825SLR Burgundy
- MDS86825SLB Smoky Blue
- MDS86825SLBLK Black
- MDS86825SLT Titanium
- MDS86825SLWT White

Posh Rollator (pg. 14)

- MDS86835SHE Pink/Zebra

Stealth Rollator (pg. 14)

- MDS86835HE Matte Black

Simplicity® Rollator (pg. 15)

- MDS86835EURO Red
- MDS86835BEURO Teal

Basic Rollators (pg. 15)

- MDS86850E Burgundy
- MDS86850EB Blue
- MDS86850EP Purple
- MDS86850EG Green
- MDS86850EO Orange

Basic Steel Rollators (pg. 16)

- MDS86850ES Red
- MDS86850EBS Blue
- MDS86850EGS Green
- MDS86850ES8 Red
- MDS86850EBS8 Blue
- MDS86850EGS8 Green

Two-Button Folding Walkers (pg. 16)

- MDS864104 Adult
- MDS86410W54 Adult
- MDS86410W4 Adult
- MDS86410J4 Junior
- MDS86410JW54 Junior
- MDS86410WJ4 Junior
- MDS864104B Basic Adult
- MDS86410W54B Basic Adult
- MDS86410W4B Basic Adult
- MDS86410J4B Basic Junior
- MDS86410JW54B Basic Junior
- MDS86410WJ4B Basic Junior

Weil Knee Walker (pg. 17)

- MDS86000G2 Weil Knee Walker

Offset Handle Canes (pg. 17)

- MDS86420 Black
- MDS86420BRZW Bronze
- MDS86420CHRW Chrome
- MDS86420B Blue Ice
- MDS86420G Green Ice

Basic Aluminum Quad Canes (pg. 18)

- MDS86222W Small Base Black
- MDS86222BRZ Small Base Bronze
- MDS86222CHR Small Base Chrome
- MDS86228 Large Base Black
- MDS86228BRZ Large Base Bronze
- MDS86228CHR Large Base Chrome

Standard Crutches (pg. 18)

- MDSV80534 Adult, Tall
- MDSV80535 Adult, Med
- MDSV80536 Youth

Child Crutches (pg. 18)

- G53314-8 Child

3-in-1 Quick Fit Crutch (pg. 18)

- MDS80540 Quick Fit

Medline Elements™ Tub and Shower Chair (pg. 20)

- MDS89745ELMB Tub and Shower Chair

Medline Elements™ Tub and Shower Bench (pg. 20)

- MDS89740ELMB Tub and Shower Bench

Medline Elements™ Transfer Bench (pg. 21)

- MDS86960ELMB Transfer Bench

Medline Elements™ Bedside Commode (pg. 21)

- MDS89664ELMB Bedside Commode

Basic Bath Benches (pg. 22)

- MDS89740RW Without Back Gray
- MDS89745R With Back Gray
- MDS89740A Without Back White
- MDS89745A With Back White
- MDS89745RA With Arms White
- MDS89740RWA Without Back, With Arms White
- MDS89745PO With Commode Opening White
- MDS89740PO With Commode Opening White
- MDS89730R Bath Stool White

Non-Padded Transfer Bench (pg. 22)

- MDS86952 White

Knurled Chrome Wall Grab Bars (pg. 23)

- MDS86012CHR 12" (30 cm)
- MDS86016CHR 16" (41 cm)
- MDS86018CHR 18" (46 cm)
- MDS86024CHR 24" (61 cm)
- MDS86032CHR 32" (81 cm)

White Finish Wall Grab Bars (pg. 23)

- G87520M 12" (30 cm)
- G87530M 16" (41 cm)
- G87540M 18" (46 cm)

- G87550M 24" (61 cm)
- G87560M 32" (81 cm)
- G87570M 36" (91 cm)

Bath Tub Assist Bar (pg. 23)

- G98006 Deluxe Tub Assist Bar

Ergonomic Hand-Held Shower (pg. 24)

- G98035 Hand Held Shower

Elevated Toilet Seat With Arms (pg. 24)

- MDS80326 Extra Wide with Legs
- G30270A Locking with Padded Arms
- MDS80320S Locking with Molded Arms

Elevated Toilet Seat Without Arms (pg. 24)

- MDS80315 Universal Locking
- MDS80324S Standard Hinged Bolt-On
- MDS80324E Elongated Hinged Bolt-On
- G30260 Locking Raised
- G30250 Basic Raised Seat

Toilet Safety Rails (pg. 24)

- G30300 250 Lbs. (113 kg)
- G30300-1 250 Lbs. (113 kg)

Basic Steel Folding Commodes (pg. 25)

- G30210-4F Aluminum
- MDS89664EL Elongated Seat
- G30213-4F Steel

Basic 3-in-1 Non-Folding Commodes (pg. 25)

- G30211-4 Chrome Plated Steel
- MDS89664A Aluminum
- G30213-4 Painted Steel
- MDS89664 Basic Steel

SAFETY

Medline Basic Beds (pg. 26)

- MDR107003E Basic Full-Electric Bed
- MDR107002E Basic Semi-Electric Bed

MedLite Beds (pg. 26)

- MDR107003L MedLite Full-Electric Bed
- MDR107002L MedLite Semi-Electric Bed

COMFORT

Hygiene/Grooming (pg. 28)

- MDSB1215 Swiveling Back Scrubber with Curved Handle
- MDSL1200 Self Wipe Aid

Dressing Aids (pg. 28)

- MDSD1100 Dressing Stick
- MDSD1330 Shoehorn, Sock Remover, with Non-Slip Handle
- MDSD1310 Long Shoehorn with T-Shaped Handle
- MDSD1210 Deluxe Sock/Stocking Aid
- MDSD1110 Button Aid/Zipper Pull

Ultra Reacher (pg. 29)

- MDS86031RCE Ultra Reacher

Mattresses (pg. 29)

- MDR237827EB 6" (15 cm) Fiber Mattress
- MDR237827B 7" (18 cm) Innerspring Mattress
- MDR230982 8" (20 cm) Premium Innerspring Mattress

Bed Assist Bar (pg. 30)

- MDS6800BA Bed Assist Bar

Side Rails (pg. 30)

- MDS89694N Full Side Rail, 4 Bar Telescoping

WEIL KNEE WALKER

- » The Weil Knee Walker is a great alternative to crutches
- » Unique, 3-wheel design protects propelling foot from hitting the rear
- » Features a full-length, extra-padded legrest for maximum comfort
- » Tool free for easy height adjustment
- » Universal design allows use with either a left or right below-the-knee injury



MDS86000G2

Item Number	Description
MDS86000G2	Weil Knee Walker

OFFSET HANDLE CANES

- » Contoured handle with foam grip provides a comfortable secure grip
- » Complete with convenient wrist strap
- » Offset handle design centers the user's weight over the strongest part of the cane
- » Ergonomically designed to provide comfort and support for partial weight-bearing users
- » Push-button, height-adjustable to find the perfect height
- » Secured with a locking ring to provide quiet, rattle-free use and added safety



MDS86420BRZW

Item Number	Description	Color	Weight Cap.: Lbs. (kg)	Adjustment: Inch (cm)	Cane Weight: Lbs. (kg)
MDS86420	Foam Hand Grip	Black	300 (136)	29-38 (74-97)	0.8 (0.4)
MDS86420BRZW	Foam Hand Grip	Bronze	300 (136)	29-38 (74-97)	0.8 (0.4)
MDS86420CHRW	Foam Hand Grip	Chrome	300 (136)	29-38 (74-97)	0.8 (0.4)
MDS86420B	Foam Hand Grip	Blue Ice	300 (136)	29-38 (74-97)	0.8 (0.4)
MDS86420G	Foam Hand Grip	Green Ice	300 (136)	29-38 (74-97)	0.8 (0.4)

BASIC ALUMINUM QUAD CANES

- » Contoured handle with foam grip provides a comfortable secure grip
- » Canes with four feet allow for greater weight bearing capabilities versus single point canes
- » Small base design allows for a faster pace; large base design gives more stability at a slower pace
- » Quad canes stand up on their own




MDS86222W

Item Number	Description	Color	Weight Cap.: Lbs. (kg)	Adjustment: Inch (cm)	Footprint: Inch (cm)	Cane Weight: Lbs. (kg)
MDS86222W	Small Base	Black	300 (136)	29-38 (74-97)	6 x 8 (15 x 20)	2 (1)
MDS86222BRZ	Small Base	Bronze	300 (136)	29-38 (74-97)	6 x 8 (15 x 20)	2 (1)
MDS86222CHR	Small Base	Chrome	300 (136)	29-38 (74-97)	6 x 8 (15 x 20)	2 (1)
MDS86228	Large Base	Black	300 (136)	30-39 (76-99)	8 x 12 (20 x 30)	2.6 (1)
MDS86228BRZ	Large Base	Bronze	300 (136)	30-39 (76-99)	8 x 12 (20 x 30)	2.6 (1)
MDS86228CHR	Large Base	Chrome	300 (136)	30-39 (76-99)	8 x 12 (20 x 30)	2.6 (1)

STANDARD CRUTCHES

- » Axillary crutches are the most common type of crutches
- » Crutches are constructed of lightweight aluminum and are not made with natural rubber latex
- » I-beam construction helps reduce the wobbliness one sometimes feels when using crutches

I-beam post



Item Number	Description	Weight Capacity: Lbs. (kg)	Patient Height: Inch (cm)	Adj. Range: Inch (cm)
Standard				
MDSV80534	Adult, Tall	300 (136)	5'10"-6'6" (178-198)	52.5-60.5 (133-154)
MDSV80535	Adult, Med	300 (136)	5'1"-5'9" (155-175)	44.5-52.5 (113-133)
MDSV80536	Youth	300 (136)	4'6"-5'2" (137-157)	36.5-44.5 (93-113)
Child				
G53314-8	Child	200 (91)	3'9"-4'7" (114-137)	31.5-37.5 (80-95)

3-IN-1 QUICK FIT CRUTCH

- » One crutch fits all three sizes (youth, adult and tall adult)
- » Features patented Red-Dot easy hand grip adjustment



Item Number	Description	Weight Capacity: Lbs. (kg)	Patient Height: Inch (cm)	Adj. Range: Inch (cm)
MDS80540	Quick Fit	300 (136)	4'7"-6'7" (137-200)	38.5-59.5 (98-151)



ELEMENTS™

ADDED STABILITY FOR ANY HOME.

Everyday bathroom activities are safer and easier with Medline. For added protection, our bath safety products are infused with Microban®, which helps prevent the growth of stain- and odor-causing bacteria.*



MEDLINE ELEMENTS™ TUB AND SHOWER CHAIR

- » Microban® antimicrobial protection inhibits the growth of stain and odor-causing bacteria, mold and mildew*
- » Warm taupe color and stylish bronze finish
- » Legs, arms and back adjust for the perfect fit
- » Extra tall armrests provide exceptional stability
- » Easy, tool-free assembly
- » Includes a clip for handheld showers**
- » 350-lb. (159 kg) weight capacity

Item Number	Description
MDS89745ELMB	Tub and Shower Chair



MDS89745ELMB

MEDLINE ELEMENTS™ TUB AND SHOWER BENCH

- » Microban® antimicrobial protection inhibits the growth of stain and odor-causing bacteria, mold and mildew*
- » Warm taupe color and stylish bronze finish
- » Legs, arms and back adjust for the perfect fit
- » Extra-wide contoured seat enhances comfort
- » Easy, tool-free assembly
- » Includes a clip for handheld showers**
- » 350-lb. (159 kg) weight capacity

Item Number	Description
MDS89740ELMB	Tub and Shower Bench



MDS89740ELMB

*This information is based upon standard laboratory tests and is provided for comparative purposes to substantiate antimicrobial activity for non-public health applications. Microban technology is not designed to replace normal cleaning practices. Antimicrobial action limited to product surface.
 **Handheld shower not included.

MEDLINE ELEMENTS™ TRANSFER BENCH

- » Microban® antimicrobial protection inhibits the growth of stain and odor-causing bacteria, mold and mildew*
- » Warm taupe color and stylish bronze finish
- » Features non-slip, height-adjustable legs, side arm and suction cup feet for greater leverage and stability
- » Easy, tool-free assembly
- » Includes a clip for handheld showers**
- » 350-lb. (159 kg) weight capacity

Item Number	Description
MDS86960ELMB	Transfer Bench



MDS86960ELMB

MEDLINE ELEMENTS™ BEDSIDE COMMODE

- » Microban® antimicrobial protection inhibits the growth of stain and odor-causing bacteria, mold and mildew*
- » Warm taupe color and stylish bronze finish
- » Features a flat seat cover and skirt around the bucket to discretely conceal the commode
- » Legs are height adjustable so commode can be used over toilet to raise the seat and enhance stability—splash guard is included
- » Elongated seat opening for added comfort—convenient toilet paper holder included
- » Easy, tool-free assembly
- » 400-lb. (181 kg) weight capacity

Item Number	Description
MDS89664ELMB	Bedside Commode



MDS89664ELMB

BASIC BATH BENCHES

- » Angled legs distribute weight over a large surface area to increase patient security
- » Anodized aluminum frame is rust resistant
- » Bath benches help patients who have difficulty sitting in a standard bathtub or standing in the shower
- » Benches with a back give extra support and comfort
- » Stools without a back provide easy access and maneuverability when bathing
- » Non-slip, rubber feet resist slipping on the bathtub or shower floor
- » Legs are height-adjustable for a proper fit



MDS89745R

Item Number	Description	Color	Weight Cap.: Lbs. (kg)	Seat Height: Inch (cm)	Seat Depth: Inch (cm)	Seat Width: Inch (cm)
MDS89740RW	Without Back	Gray	250 (113)	14.5-19.5 (37-50)	12 (30)	20 (51)
MDS89745R	With Back	Gray	250 (113)	14.5-19.5 (37-50)	12 (30)	20 (51)
MDS89740A	Without Back	White	250 (113)	14.5-19.5 (37-50)	12 (30)	20 (51)
MDS89745A	With Back	White	250 (113)	14.5-19.5 (37-50)	12 (30)	20 (51)
MDS89745RA	With Arms	White	350 (158)	14.5-19.5 (37-50)	12 (30)	20 (51)
MDS89740RWA	Without Back, With Arms	White	350 (158)	14.5-19.5 (37-50)	12 (30)	20 (51)
MDS89745PO	With Commode Opening	White	300 (136)	14.5-19.5 (37-50)	14.5 (37)	16 (41)
MDS89740PO	With Commode Opening	White	300 (136)	14.5-19.5 (37-50)	14.5 (37)	16 (41)
MDS89730R	Bath Stool	White	300 (136)	14-19 (36-48)	13 (33)	13 (33)

NON-PADDED TRANSFER BENCH

- » Allows patients to enter and exit tubs safely
- » Suction cup feet grip the tub floor for added safety
- » Side arm provides extra stability and leverage
- » Height-adjustable legs accommodate most tub wall heights
- » Textured backrest and seating surface offer a non-slip surface with holes for drainage
- » Backrest attaches for right or left facing tubs easily without tools



MDS86952

Item Number	Description	Color	Weight Cap.: Lbs. (kg)	Seat Height: Inch (cm)	Seat Width: Inch (cm)	Seat Depth: Inch (cm)
MDS86952	Basic	White	300 (136)	16.5-21.5 (42-55)	17 (43)	16.5 (42)

WALL GRAB BARS

- » Locate grab bars at any angle in the bathroom to help prevent falls
- » Use grab bars in the kitchen, along stairways, or in the laundry room for additional fall prevention
- » Combine grab bars with an elevated toilet seat or transfer bench for additional stability
- » Knurled-texture steel bars are easy to grip even in wet bathing areas
- » White painted steel with white screws blend in to match most bathroom decors



MDS8601* SERIES

Standard Knurled Chrome Grab Bars

Item Number	Length: Inch (cm)	Diameter Inch: (cm)
MDS86012CHR	12 (30)	1 (3)
MDS86016CHR	16 (41)	1 (3)
MDS86018CHR	18 (46)	1 (3)
MDS86024CHR	24 (61)	1 (3)
MDS86032CHR	32 (81)	1 (3)

Basic White Finish Grab Bars

Item	Length: Inch (cm)	Diameter: Inch (cm)
G87520M	12 (30)	1 (3)
G87530M	16 (41)	1 (3)
G87540M	18 (46)	1 (3)
G87550M	24 (61)	1 (3)
G87560M	32 (81)	1 (3)
G87570M	36 (91)	1 (3)

BATH TUB ASSIST BAR

- » Give patients confidence and support with an extra grip to help them transfer into the tub
- » Step through clamp design provides unobstructed entry/exit path allowing better hand positioning
- » Tool-free installation
- » Rubber-lined pads won't mar the tub surface
- » Adjust to fit tubs 2.75"–6.25" (7–26 cm)

Item Number	Description
G98006	Deluxe Tub Assist Bar



G98006

ERGONOMIC HAND-HELD SHOWER

- » Extra long 86" (218 cm) hose allows user to bathe comfortably while seated in the shower
- » Diverter valve allows you to use your existing shower head
- » Adapts to existing shower arm, delivering both standard and pulsating spray



G98035

Item	Description
G98035	Hand Held Shower

ELEVATED TOILET SEAT

- » Makes using the toilet easier for those with trouble bending
- » Available in a variety of locking mechanisms and mounting methods



G30270A

Item Number	Description	Weight Capacity: Lbs. (kg)	Seat Elevation: Inch (cm)
With Arms			
MDS80326	Extra Wide with Legs	600 (272)	4 (10)
G30270A	Locking with Padded Arms	350 (159)	5 (12)
MDS80320S	Locking with Molded Arms	300 (136)	5 (12)
Without Arms			
MDS80315	Universal Locking	250 (113)	4.75-6.75 (12-17)
MDS80324S	Standard Hinged Bolt-On	300 (136)	4 (10)
MDS80324E	Elongated Hinged Bolt-On	300 (136)	4 (10)
G30260	Locking Raised	350 (159)	5 (12)
G30250	Basic Raised Seat	300 (113)	5 (12)

TOILET SAFETY RAILS

- » Adjustable handles rotate to accommodate a wide range of secure and comfortable positions
- » Closed cell foam armrests provide a comfortable, sure grip
- » Adjustable height ranges from 26"-31" (66-79 cm)
- » Accommodates standard or elevated toilet seats
- » 250-lb. (113 kg) weight capacity



G30300

Item Number	Weight Capacity: Lbs. (kg)	Height x Width: Inch (cm)
G30300	250 (113)	26-31 (66-79) x 18-24 (46-61)
G30300-1	250 (113)	26-31 (66-79) x 18-24 (46-61)

BASIC STEEL FOLDING COMMODES

- » Folding frame design assembles in seconds without tools
- » Commode folds down to reduce storage space by 35% and makes delivery easier
- » All models include pail, lid and splash guard
- » Clip-on seats can be removed for easy cleaning



G30213-4F

Item Number	Description	Weight Capacity: Lbs. (kg)	Seat Height Range: Inch (cm)	Seat Depth: Inch (cm)	Width Between Arms: Inch (cm)	Overall Width: Inch (cm)
G30210-4F	Aluminum	250 (113)	17.5-21.5 (44-55)	16 (41)	19 (48)	21.25 (54)
MDS89664EL	Elongated Seat	400 (182)	17-23 (43-58)	16.5 (42)	19 (48)	22 (56)
G30213-4F	Steel	350 (159)	16.5-22.5 (42-57)	16 (41)	19 (48)	21.25 (54)

BASIC 3-IN-1 NON-FOLDING COMMODES

- » Commodes can be used bedside for those who may not be mobile enough to walk to the bathroom
- » Can also be used in the bathroom as an elevated toilet seat or as toilet safety rails when combined with most standard toilets
- » All models include pail, lid and splash guard
- » Clip-on seats can be removed easily for cleaning
- » Gray powder-coated steel commodes are rust-resistant and durable
- » Aluminum commodes are lightweight and rustproof



MDS89664

Item Number	Description	Weight Capacity: Lbs. (kg)	Seat Height Range: Inch (cm)	Seat Depth: Inch (cm)	Width Between Arms: Inch (cm)	Overall Width: Inch (cm)
G30211-4	Chrome Plated Steel	350 (159)	16.5-22.5 (42-57)	16 (41)	19 (48)	21.25 (54)
MDS89664A	Aluminum	300 (136)	16-20 (41-51)	15 (38)	18.5 (47)	25.5 (65)
G30213-4	Painted Steel	350 (159)	16.5-22.5 (42-57)	16 (41)	19 (48)	21.25 (54)
MDS89664	Basic Steel	350 (159)	15-21 (38-53)	15 (38)	19 (48)	24 (61)

MEDLINE BASIC BEDS

- » Provides full automation at an economical price
- » Light weight—foot section weighs just 42 lbs. (19 kg)
- » Convenient DC motor box snaps off quickly, making delivery, set-up, and routine maintenance much easier
- » 9V battery back up allows movement during power failures (battery not included)
- » Recessed headspring supports and foot assembly reduce the chance of pinching and injuries
- » Height range: 15"–20" (38–51 cm)
- » Bed deck dimensions: 80" (203 cm) x 36" (91 cm)
- » Assembled weight: 155 lbs. (70 kg)
- » Weight capacity: 450 lbs. (204 kg)
- » Casters: four (two locking)
- » Warranty: 5 years for electrical and mechanical components; lifetime for the bed frame



MDR107002E

Item Number	Description
MDR107003E	Basic Full-Electric Bed
MDR107002E	Basic Semi-Electric Bed

MEDLITE BEDS

- » Exceptionally light and easy to use
- » Can convert to a low bed using MedLite Low Bed Conversion Kit—sold separately
- » Premium Linak® motors provide dynamic movement without disturbing a partner or surrounding residents
- » Heavy-duty hand pendant
- » Recessed headspring supports and foot assembly reduce the chance of pinching and injuries
- » Height range: 14.5"–22.5" (37–57 cm)
- » Assembled weight: 145 lbs. (66 kg)
- » Weight capacity: 450 lbs. (204 kg)
- » Casters: four (two locking)
- » Warranty: 5 years for electrical and mechanical components; lifetime for the bed frame



MDR107002L

Item Number	Description
MDR107003L	MedLite Full-Electric Bed
MDR107002L	MedLite Semi-Electric Bed

MORE COMFORT AT HOME.

Comfort is more than just a soft place to sit. It's the little things that help reduce the strain and challenges of everyday living.



HYGIENE/GROOMING

- » Aids for people needing a little help with day-to-day bathroom needs

Item Number	Description
MDSB1215	Swiveling Back Scrubber with Curved Handle
MDSL1200	Self Wipe Aid



MDSL1200



MDSB1215



MDSD1100

DRESSING AIDS

- » Assists a person putting on socks, shirts and jackets; pulling up slacks or skirts, or removing socks
- » Essential for anyone who has difficulty bending, limited range of motion or the use of only one arm or hand
- » Also ideal for those that have difficulty manipulating buttons or zippers due to arthritis or limited fine motor control



MDSD1330



MDSD1310

Item Number	Description
MDSD1100	Dressing Stick
MDSD1330	Shoehorn, Sock Remover, with Non-Slip Handle
MDSD1310	Long Shoehorn with T-Shaped Handle
MDSD1210	Deluxe Sock/Stocking Aid
MDSD1110	Button Aid/Zipper Pull



MDSD1210



MDSD1110

ULTRA REACHER

- » Specially designed with a hand grip allowing the use of all four fingers for greater grip strength
- » Non-slip rubber tabs over composite jaws latch onto items more securely
- » Weighs only 14 oz. (0.5 kg)
- » 31" (79 cm) reach




MDS86031RCE

Item Number	Description
MDS86031RCE	Ultra Reacher

MATTRESSES

- » Durable 8 oz. vinyl mattresses feature antimicrobial ticking
- » Inverted side seams for patient comfort and caretaker convenience
- » Your choice of innersprings or high performance fiber



MDR237827

Item Number	Description	Dimensions: Inch (cm)	Weight Capacity: Lbs. (kg)
MDR237827EB	6" (15 cm) Fiber Mattress	37 x 84 x 6 (94 x 213 x 15)	250 (113)
MDR237827B	7" (18 cm) Innerspring Mattress	37 x 84 x 7 (94 x 213 x 18)	250 (113)
MDR230982	8" (20 cm) Premium Innerspring Mattress	35 x 80 x 8 (89 x 203 x 20)	400 (181)

BED ASSIST BAR

- » Features a soft grip for comfort, a pocket for handy storage and full length legs and straps for added safety
- » Sturdily constructed bar offers quick and easy tool-free assembly

Item Number	Description
MDS6800BA	Bed Assist Bar



MDS6800BA

SIDE RAILS

- » Tool-free installation!
- » Attractive hammertone finish hides fingerprints and complements home décor
- » Spacing between bars is less than 4" (10 cm) for patient safety
- » MDS89694 and MDS89690 adjust in length from 48"-72" (122-183 cm)

Item Number	Description
MDS89694N	Full Side Rail, 4 Bar Telescoping




MDS89694N





Town Name Location 1

123 Main Street
Anytown, IL
00000
T: 800-123-4567
F: 800-987-6543
www.sitename.com



Town Name Location 2

123 Main Street
Anytown, IL
00000
T: 800-123-4567
F: 800-987-6543
www.sitename.com



Town Name Location 3

123 Main Street
Anytown, IL
00000
T: 800-123-4567
F: 800-987-6543
www.sitename.com



Medline Industries, Inc. One Medline Place, Mundelein, IL 60060

Medline United States
1-800-MEDLINE (633-5463)
medline.com | info@medline.com

Medline Canada
1-800-396-6996
medline.ca | canada@medline.com

Medline México
01-800-831-0898
medlinemexico.com | mexico@medline.com

FOLLOW US    