

## CHAD® Bonsai® Velocity Pneumatic Oxygen Conserver

HCPCS:



- Delivers a uniform oxygen pulse at increments: 1, 2, 3, 4, 5 and 6
- Two continuous flow (CF) settings preset @ 2 LPM and 4 LPM
- Delivers up to a 6:1 savings ratio from 14-40 breaths per minute
- Senses a breath/delivers oxygen in the first half of the inspiratory cycle
- Operating pressure range: 200-3000 PSI
- Light weight and quiet operation
- Easy-to-operate single selector knob
- Compact design fits easily in virtually any oxygen carry bag

### PRODUCTS

#### Description

OM-812	CHAD® Bonsai® Velocity Pneumatic Oxygen Conserver, EA
--------	---

### SPECIFICATIONS

#### Description

Warranty	Limited 2 Year
----------	----------------

### ACCESSORIES

#### Description

OM-900-FH	S-Handle for OM-900/800/500 Series, 1/e, EA
OM-900-FEC	E-Clip for Conserver Handles, 10/bag, BAG
SEAL100	Seal Washer,Conservers/Regulators,10/pkg, PKG