

## CHAD® Evolution Electronic Oxygen Conserver

HCPCS:



- Delivers a uniform oxygen pulse at increments: 1, 2, 3, 4, 5, 6, 7 Plus CF preset @ 2 LPM
- Minimum 2-year battery life\* utilizing 2 AA batteries
- Delivers a 5:1 savings ratio from 14-40 breaths per minute
- Senses a breath/delivers oxygen in the first half of the inspiratory cycle
- Operating pressure range: 200-3000 PSI
- Weighs only 14.9 ounces with batteries
- Single push button selects the equivalent LPM dosage settings
- Auto shut-off feature conserves battery life
- Fits easily in virtually any carry bag

### PRODUCTS

#### Description

OM-900	CHAD® Evolution Electronic Oxygen Conserver, EA
--------	-------------------------------------------------

### SPECIFICATIONS

#### Description

Warranty	Limited 2 Year
----------	----------------

### ACCESSORIES

#### Description

SEAL100	Seal Washer,Conservers/Regulators,10/pkg, PKG
OM-900-FEC	E-Clip for Conserver Handles, 10/bag, BAG
OM-900-FH	S-Handle for OM-900/800/500 Series, 1/e, EA
OM-900-BD	BatteryDoors&Screws,OM-900/OM-900M,5/pkg, EA