

D. ProBasics® Folding Canes

- Black or bronze anodized aluminum finish
- "T" Shaped handle
- Folds easily for travel or storage
- Adjustable height 32.5" to 36.5"
- Convenient carry and display pouch
- HCPCS Code: E0100

Model #	Description
PB2055B	Folding Cane - Bronze
PB2055K	Folding Cane - Black



E. ProBasics® Hemi Walker

- Anodized aluminum frame
- Lightweight compact design folds easily
- Special angled footpieces provide stability
- Adjustable height from 28" to 33"
- HCPCS Code: E0135


Model #	Description
PB6055	Hemi Walker



F. ProBasics® Aluminum Crutches

- Solid one piece lightweight anodized aluminum frame
- High impact plastic armrests and hand grips
- Fully assembled with cushions, handgrips, and rubber tips
- Push button adjustment allows for safe, accurate and quick adjustment of height
- Height adjustment scale on center tube
- HCPCS Code: E0114

Model #	Description
PB1055A	Aluminum Crutches - Adult
PB1055T	Aluminum Crutches - Tall
PB1055Y	Aluminum Crutches - Youth



Model #	Height Adjustment	Units Per Carton	Weight Capacity	Limited Warranty
PB1057	31"-40"	10	250 lb.	2 Years Frame
PB2057 PB2057K PB2057B	31"-40"	10	250 lb.	2 Years Frame
PB2058	29"-38"	20	500 lb.	2 Years Frame
PB4057 PB3057	31"-40"	6	250 lb.	2 Years Frame
PB4058 PB3058	29"-38"	2	500 lb.	2 Years Frame
PB6055	28"-33"	4	250 lb.	2 Years Frame
PB2055B PB2055K	32.5"-36.5"	20	250 lb.	2 Years Frame
PB1055A	45"-53"	8	350 lb.	2 Years Frame
PB1055T	53"-61.5"	8	350 lb.	2 Years Frame
PB1055Y	37"-46"	8	350 lb.	2 Years Frame

All references to HCPCS codes: Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare program, and should consult an attorney or other advisor to discuss specific situations in further detail.