

## Nitro Duet Rollator and Transport Chair



HCPCS:




FIGURE A




FIGURE B




FIGURE C



FIGURE D



FIGURE E




FIGURE F




FIGURE G



FIGURE H

- Combines the features of a rollator and transport chair into one unit. The user can ambulate independently or be pushed safely by a caregiver
- The tool-free two-position contoured backrest is reversible depending on whether you are using the Nitro Duet as a rollator or transport chair
- Back support height easily adjusts with tool-free buttons
- Caster fork design and large 10" front casters enhances turning radius and rolling comfort
- Cross-brace design allows for side-to-side folding and added stability
- Brake cables routed inside the frame add safety and convenience

### PRODUCTS

### Description

RTL10266DT	Nitro Duet Rollator and Transport Chair, EA
------------	---

### SPECIFICATIONS

### Description

Front Wheels	10"
Handle Height	33"   37"
Inside Hand Grip Width	18"
Number of Wheels	4
Rear Wheel	8"
Seat Depth	14"
Seat To Floor Height	20.5"
Seat Width	18"
Warranty	Lifetime Limited
Weight Capacity	300 lbs.