

Gel Grip, Aluminum Offset Canes, Height Adjustable

HCPCS: E0100




FIGURE A




FIGURE B



FIGURE C



FIGURE D

- Ergonomically designed handle with soft Gel Grip provides comfort and security
- Comes standard with premium Gel Grip and wrist strap
- Manufactured with sturdy, extruded aluminum tubing
- Easy-to-use, one-button height adjustment with locking ring Locking ring prevents rattling
- Available in Blue Crackle (main image), Red Crackle (Figure A) or Black (Figure B)
- Handle height adjusts from 30" to 39"
- Weight Capacity 300 lbs

PRODUCTS

Description

RTL10372BC	Gel Grip, Aluminum Offset Canes, Height Adjustable, CS
RTL10372BK	Gel Grip, Aluminum Offset Canes, Height Adjustable, CS
RTL10372RC	Gel Grip, Aluminum Offset Canes, Height Adjustable, CS

SPECIFICATIONS

Description

Handle Height	30" - 39"
Warranty	Limited Lifetime
Weight Capacity	300 lbs