

Therma Moist Heating Pad and Floor Display



FIGURE A



FIGURE B



FIGURE C



RETAIL BOX



INSET E

- Delivers soothing moist heat to a variety of areas
- The positioning strap with hook and loop fastener closure ensures a secure, proper fit for maximum treatment (Figure A)
- The easy-to-use hand control features four heat settings and four time settings (Figure B)
- Moist heat is more effective than dry heat because moisture transmits heat more efficiently than air
- The Therma Moist Heating Pad has an outer flannel cover that draws the air's humidity and retains it
- As the heating pad warms, the moisture is released the flannel cover and onto the treated area
- To prevent accidents, hand control will not allow unit to turn on unless the temperature is set at 0° (Figure C)
- The heating pad will automatically shut off after treatment
- Three sizes are available to provide added convenience
- Retail packaged

PRODUCTS

	Description
10890	Standard, 14" x 27", 1/ea
RTL10893	Medium, 14" x 14", 1/ea

SPECIFICATIONS

	Description
Heat Settings	115 130 150 165 Fahrenheit, 45 55 65 75 Celsius
Time Settings	15 30 45 60 Minutes
Power Consumption	90 - 130 Volts, 40 - 70 Hz
UL/cUL/TUV Listed	Yes
Warranty	Limited 1 Year