

Poly-Fly High Strength, Lightweight Wheelchair/Flyweight Transport Chair Combo

HCPCS: E1038 or K0004+E0951+E0971+E0978



FIGURE A



FIGURE B



FIGURE C



FIGURE D

- Can be used as standard self-propelled wheelchair or a transport chair all in one
- Quick-release 24" wheels can be removed to transition self-propelled chair to a transport chair (Figure A - B)
- All-aluminum frame
- Comes with two sets of wheel locks: One for use with wheelchair, one for use with transport chair
- Back folds down with deluxe back release (Figure C)
- Nylon upholstery
- Comes standard with swing-away footrests with composite foot plates and heel loops
- 8" casters
- Comes with anti-tippers
- STDELR-TF - Swing-away Elevated Leg rest, Tool-free, (Optional front riggings)
- Weight Capacity: 250 lbs

PRODUCTS

PRODUCTS	Description
TR16	Poly-Fly High Strength, Lightweight Wheelchair/Flyweight Transport Chair Combo, CS
TR18	Poly-Fly High Strength, Lightweight Wheelchair/Flyweight Transport Chair Combo, CS
TR20	Poly-Fly High Strength, Lightweight Wheelchair/Flyweight Transport Chair Combo, CS

SPECIFICATIONS

SPECIFICATIONS	Description
Back Height	16"
Carton Shipping Weight	38 lbs. 39 lbs. 41 lbs.
Overall Height w/ 24 Wheels	37"
Overall Length w/ Rigging	42"
Overall Width with 24 Wheels (TR16)	22.375" (open) 24.5" (open with mag-wheels) 12" (closed) 23.625" (open) 26.5" (open with mag-wheels) 12" (closed) 26.25" (open) 29" (open with mag-wheels) 12" (closed)
Seat Depth	16"
Seat To Floor Height	19.5"
Seat Width	16" 18" 20"
Warranty	Limited Lifetime
Weight Capacity	250 lbs
Weight w/ 24 Wheels (TR16 TR18 TR20)	33 lbs. 34 lbs 36 lbs.

PARTS

PARTS	Description
STDS2A1621	RearAxle,QuickRelease,Polyfly,1/ea, EA
STDS1012	Rear Wheel Bearing,11.5 x 32 x 9.5mm,1/e, EA
STDS2A161R	Brake,Right,Polyfly,24"wheel,1/ea, EA
STDS2A161L	Brake,Left,Polyfly,24"wheel,1/ea, EA
STDS2A163R	Brake,Right, 8"Wheel,Polyfly,1/ea, EA

PARTS	Description
STDS2A163L	Brake,Left,8"Wheel,Polyfly,1/ea, EA
STDS2A1603	8"Rear Wheel,Polyfly,1/ea, EA
STDS1005	Rear Axle Bolts,Sentra,1/pkg, EA
STDS2A1611	RearAxle Receptacle,Polyfly,1/ea, EA
STDS2A1610	SlideTube w/Bushing,Polyfly,1/ea, EA
STDS2A1686	Right Side Frame,TR16 and TR18, 1ea, EA
STDS2A1687	Left Side Frame,TR16 and TR18,1 ea, EA
STDS2A1686-20	Right Side Frame,TR20,1ea, EA
STDS2A1687-20	Left Side Frame,TR20,1ea, EA
STDS2A1612	Round Brake Pad,8",Polyfly,1/ea, EA
STDS2A161P	WheelBrakePad24"Polyfly,1 ea, EA
WCP 1287	Frame Plug,Polyfly,1/ea, EA
STDSSB	Spring Btn for Arm Rest,1/ea., EA
STDS004012	Handrim,Black,24"Wheel/1 ea.", EA
STDS2A1600	24" Wheel,BLK,HandRim TR18,1/ea, EA